



The strong system partner for the future

Product range 2019 / 2020

Quality - Partnership - Economy

Basic care in our health systems is an important pillar for the functioning of our society. We are a manufacturer and system supplier of medical consumables and consider ourselves a partner and service provider for basic care in all medical institutions. Quality, innovation and beneficial care options motivate us each day to make a small and yet significant contribution to modern and future oriented healthcare. The employees of MEDITRADER GMBH are dedicated to this core principle and fulfil these demanding tasks competently and professionally to achieve a higher health standard.

Your MEDITRADER Team





Meditrade – The strong system partner for the future

Meditrade has been successful in the market for decades as the manufacturer of the highest quality risk protection products. By purchasing the 90 year old company Karl Beese in 2010, we expanded into the sector of high quality

bandages and incontinence products. We build on our expert knowledge and constantly develop new products so that we are always able to offer the best professional solution with the maximum efficiency.



Our product range

- ➡ Gloves
- ➡ Masks, hoods and surgical clothing
- ➡ Dressings, bandages and plasters
- ➡ Individual sets and surgical requirements
- ➡ Incontinence
- ➡ Disinfection and care
- ➡ Hygiene, cleaning and cleaning products
- ➡ Ward and patient requirements

Your demands - Our motivation



We work with two logistics centres which enable us to cover the increased system partnership requirements. South (Kiefersfelden) and North

(Neuenkirchen / Hamburg). State-of-the-art logistics processes are constantly adapted to the latest customer requirements.

Good advice - more important than ever

Meditrade is a manufacturer and system partner for the complete supply of medical devices and services to hospitals and clinic-related facilities.



Our core competences

Products:

High quality and innovative quality products for medical care

Service:

On site advice from our competent team of consultants as well as a contact person who is always available in our customer service

Knowledge:

Expert knowledge in all industry relevant areas

Organisation:

Supply chain management solutions which increase profitability guarantee a comprehensive solution for our customers



Complex needs assessment - How to

Comprehensive supply coupled with a tailor-made needs assessment is still the greatest challenge for cost effective procurement.

Thanks to our high standards and proximity to the market, we have extensive knowledge that enables us to provide suitable products for your specific applications across the whole system and to optimise them further for you.



Gloves

Latex examination gloves	
Reference™	9
Gentle Skin® sensitive	9
Gentle Skin® classic (x-long)	10
Gentle Skin® grip	11
Gentle Skin® compact +	11
Gentle Tec®	12
Gentle Skin® HiRisk	13
Gentle Skin® Aloecare™	14
Gentle Skin® black	15
Nitrile examination gloves	
Nitril® 3000	16
Nitril® NextGen®	16
Nitril® Sensory®	17
Nitril® BestGen®	18
Nitril® 3000 x-long blue	19
Nitril® black	20
Nitril® magenta	21
Nitril® LIOX	22 / 23
Vinyl examination gloves	
Vinyl 2000	24
Vinyl 2000 Stretch PF	25
Cotton and hygiene gloves	
Cotton gloves	26
Froggy™	27
Élastify®	27
Finger cots	
Latex finger cots	28
Household gloves	
MyClean® easy touch	29
Sterile examination gloves	
Copolymed™	31
Gentle Skin® steril	32
Nitril® steril	33
Surgical gloves	
Reference™ OP	35
Gentle Skin® Premium OP™	36
Gentle Skin® Smooth OP	37
Gentle Skin® Micro OP®	38
Gentle Skin® Securitex® OP	39
Gentle Skin® Isopretex® OP	40
Neopretex® OP	41
Finessis® Zero®	43
Finessis® Aegis®	44 / 45
Finessis® Corium®	46
X-Guard	47

Protective clothing

Masks	
Suavel® Protec Plus	51
Suavel® Protec	51
Suavel® Comfort Plus	52
Suavel® Comfort	52
Suavel® Antifog	53
Suavel® Antifluid	53
Suavel® Sensima	54
Hoods	
Suavel® Apollo	55
Suavel® Astrid	55
Jackets / Gowns / Aprons	
BeeSana® PP protective gown	57
BeeSana® PP gown	57
BeeSana® PP/PE gown	58
BeeSana® SMS gown	58 / 59
BeeSana® SSS gown	59
PE disposable apron	60

Wound care

Lap sponges	
BeeSana® Lap sponge x-ray thread	64 / 65
Swabs & tamponades	
BeeSana® Gauze ball	66
BeeSana® Gauze ball x-ray thread	67
BeeSana® Thread swab	68
BeeSana® Preparation balls x-ray thread	69
BeeSana® Pointed swab / Pointed swab x-ray thread	70 / 71
BeeSana® Cellulose swabs	72
Alcohol swabs	72
ABE® Ribbon gauze / Ribbon gauze x-ray thread	73
Gauze and gauze cutting	
BeeSana® Gauze cutting	74
BeeSana® Medical gauze	74
Swabs	
BeeSana® Gauze swab	75 - 77
BeeSana® Gauze swab x-ray thread	78 / 79
BeeSana® Slit gauze swabs	80
BeeSana® Non-woven swab	81
BeeSana® Slit non-woven swab	82
ABE® drainal	83
ABE® tracheal	83
BeeSana® Absorbent dressing pads	84
BeeSana® Eye pad	85
Wound plaster & bandages	
ABE® Injection plaster	86
ABE® light	86
ABE® Absorbent adhesive dressing	87
ABE® derm	88
ABE® derm plus	88
Fixation plaster & non-woven material	
ABE® Cannula retention dressing I.V.	89
ABE® Fixation non-woven fabric	89
ABE® silk	90
Fixation bandages & bandages	
ABE® Tubular net bandage	91
ABE® Head bandage	91
ABE® Nose bandage	92
Suspensory	92
ABE® Gauze bandage	93
ABE® Last	94
ABE® Last crepe	95
ABE® fix color	96
ABE® Last adhesive	97
BeeSana® SafeCut	98
ABE® fix adhesion c	99
Compression & support bandages	
ABE® extra	100
ABE® extra adhesion	101
ABE® Ideal 100	102
ABE® lastic	103
ABE® universal	104
ABE® Ven	105
ABE® Padding bandage	106
ABE® compression bandage roll	107
Testicular roll	107
ABE® Tubular bandage	108
ABE® Gilchrist bandage	109
ABE® Body bandage	110
ABE® Heart catheter pressure bandage	111
ABE® Limb holder	111
Modern wound care	
PolyMem® Membrane bandage	112 - 115

Sets and surgical supplies

BeeSana® Ward catheter set	118
BeeSana® PEG bandage kit	119
BeeSana® Bandage change kit	120
BeeSana® Suture removal kit	120
BeeSana® Ward dialysis kit	121
BeeSana® Ward protection set	121

Incontinence products

Incontinence briefs	
Superslip®	124
Superslip® Air Comfort	125
Medi-Pants	125
Incontinence pads / inlays	
Liberty Light	126
Incontinence pads	
Superform®	126
Superform® Air	127
BeeSana® Absorbent pad	127 / 128
BeeSana® Rectangular pad	128
BeeSana® Male incontinence pad	129
BeeSana® Women's incontinence pad	129
Fixation pants	
Light fixation pants	130
Panty fixation pants	131
Incontinence pads	
Layered underpad	132
BeeSana® Comflock Super	133
Nursing pads	134
Baby nappies	
New Quik	135

Disinfection / Care

Disinfection	
Alcoman® +	138
Alcoman® Gel	139
Ethasept®	140
Medizid® AF	141
Medizid® Rapid Optima	142
Medizid® Rapid	143
Medizid® Rapid Disinfection wipes	144
Medizid® Rapid QF Flowpack	145
Medizid® Alpha	146
Medizid® Alpha Disinfection wipes	147
Medizid® Alpha Flowpack	148
Instrusol® AF +	149 / 150
Instrusol® Enzym	151 / 152
InstruPlus® forte	153 / 154
Care	
Gentle Med® wash lotion	155
Gentle Med® MRSA effective wash lotion	155
Gentle Med® Skin care lotion (O/W)	156
Gentle Med® Cream lotion (W/O)	156
Gentle Med® Active Gel	157
Gentle Med® Oil bath	157
Gentle Med® skin cleansing and care foam	158
Gentle Med® Moist care wipes	158
Application aids	
Dosing pumps / attachments	159
(Foam) spray heads	159
Canister tap / Canister key	159
Dosing bottle / measuring beaker	159
Disinfectant dispenser	160

Hygiene / Cleaning / Cleaning products

beopur WC fresh	164
BioSan® Fresh scent	165
Reliable hygiene in patient rooms	166- 169
Cleaning in commercial kitchens	170 - 173

Ward and patient supplies

Patient & personal hygiene	
BagBath® Haircare	176
BagBath®	176
BeeSana® Wash glove	178
Protective towel	178
Protective bibs	179
BeeSana® Non-woven disposable panty	179
Anti-slip sock universal	180
Surgical sock	180
PE shoe covers	181
Cellulose swabs	181
Layered cellulose	181
Rollicel® / Medical crepe	182
Protective sheets	182
Protective covering	183
Multi-purpose tissues	183
Bed covers	183
Universal tissues	184
Micro-fibre tissues	184
Wipes	185
Cosmetic tissues	185
V-fold paper hand towels	185
Sick Bag (BaSick Bag)	186
Kidney dish	186
Disposable comb	186
Disposable razor	187
Urine bottle	187
Urine bottle holder	187
Treatment care & accessories	
Cotton buds	188
Lemon glycerine stick	188
Wooden tongue depressor	189
Bandage clips	189
Diagnosis & accessories	
Nurse's stethoscope	190
Tourniquet	190
Dosicair®	191
Urine beaker	191
Ultrasound scan gel	191
Pain therapy & accessories	
Repose®	192
Cold-warm compress	193
Non-woven cover for compresses	193
Ward supplies & miscellaneous	
Rescue blanket	194
Medicine beaker	194
Medication dispenser	195
Medi-Box	195
Wall holder	196
Refuse bags HD	197
Kitchen towel	198
Toilet paper	198 / 199

Gloves

Protective clothing

Wound care

Sets / surgical supplies

Incontinence

Disinfection / care

Hygiene / Cleaning /
Cleaning productsWard and
patient supplies



We know a lot about gloves!

We as Meditrade see it as a priority during manufacturing of examination and surgical gloves to offer the user gloves that combine very high mechanical strength with optimal tactile sensitivity. The top priority, however, is the protection of the user and the patient. We are constantly investing into the know-how of product development and into our production sites, because we know that this represents the basic work for an optimized product range. Various models of latex gloves made of natural rubber or nitrile gloves made of nitrile butadiene rubber are complemented by very economical vinyl gloves and other alternative materials in order to always be able to offer the right products for users.

Reference™

Our powdered latex examination glove

Characteristics

- examination glove made of latex acc. DIN EN 455
- lightly powdered
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit
- suitable for contact with food



non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1121XS	100 pcs	10 boxes	06177515
S	1121S	100 pcs	10 boxes	06177521
M	1121M	100 pcs	10 boxes	06177538
L	1121L	100 pcs	10 boxes	06177544
XL	1121XL	100 pcs	10 boxes	02243899

Application

In all areas where the user is exposed to a risk of contamination. Suitable for highly infectious situations in hospitals and in private practice. Ideal for use in industrial areas. Very popular among users who tend to have moist hands.

Gentle Skin® sensitive

Our most popular latex examination glove

Characteristics

- examination glove made of latex according to DIN EN 455
- protective glove according to PPE
- skin friendly: powder-free, very low in proteins and low in accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit
- suitable for contact with food



non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1221RXST	100 pcs	10 boxes	02243563
S	1221RST	100 pcs	10 boxes	02243586
M	1221RMT	100 pcs	10 boxes	02243592
L	1221RLT	100 pcs	10 boxes	02243600
XL	1221RXLT	100 pcs	10 boxes	02243617

Application

In all areas where the user is exposed to a risk of contamination. Suitable for highly infectious situations in hospitals and in private practice. Ideal for use in the care sector, in laboratories, in the pharmaceutical industry, in cleaning and in the industrial sector.

Gentle Skin® classic / classic x-long

Our extra strong latex examination glove

Characteristics

- examination glove made of latex according to DIN EN 455
- protective glove according to PPE
- high wall thickness due to 20% more latex content than conventional latex gloves
- skin friendly: powder-free, very low in proteins and low in accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit
- suitable for contact with food
- classic x-long: with extra long shaft (at least 29 cm)



Application

In all areas where the user is exposed to a risk of contamination. Ideal for all emergency situations, rescue services and intensive care. Suitable for highly infectious situations in hospitals and for private practice. Can be used in the care sector, in laboratories, in the pharmaceutical industry, in cleaning and in the industrial sector.

Gentle Skin® classic, unsteril

Size	REF	Box contents	Shipping unit	PZN
XS	1221RXS	100 pcs	10 boxes	02243698
S	1221RS	100 pcs	10 boxes	02243706
M	1221RM	100 pcs	10 boxes	02243729
L	1221RL	100 pcs	10 boxes	02243735
XL	1221RXL	100 pcs	10 boxes	02243801

Gentle Skin® classic x-long, non-sterile

S	1223S	100 pcs	10 boxes	11349906
M	1223M	100 pcs	10 boxes	11349935
L	1223L	100 pcs	10 boxes	11349958

Gentle Skin® grip

Our easy-grip latex glove

Characteristics

- examination glove made of latex according to DIN EN 455
- protective glove according to PPE
- optimum grip safety
- skin-friendly: powder-free, very low in proteins and low in accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit
- suitable for contact with food



non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1221GRIP-XS	100 pcs	10 boxes	02470253
S	1221GRIP-S	100 pcs	10 boxes	02470282
M	1221GRIP-M	100 pcs	10 boxes	02470299
L	1221GRIP-L	100 pcs	10 boxes	02470307
XL	1221GRIP-XL	100 pcs	10 boxes	02470313

Application

Especially suitable for dental surgery and because of its excellent grip it enables the safe handling of instruments and machines.

Gentle Skin® compact +

Our latex glove for the medical and industrial sector

Characteristics

- examination glove made of latex according to DIN EN 455
- popular disposable gloves in industry
- non-slip for a secure grip
- skin-friendly: powder-free
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit
- suitable for contact with food



non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1221XSI	100 pcs	10 boxes	05522453
S	1221SI	100 pcs	10 boxes	05522447
M	1221MI	100 pcs	10 boxes	05522430
L	1221LI	100 pcs	10 boxes	05522424
XL	1221XLI	100 pcs	10 boxes	05521198

Application

Particularly popular in the industrial sector and in food processing. Ideal for household and cleaning.

Gentle Tec®

Our extra long latex glove for clean rooms and critical applications

Characteristics

- examination glove made of latex acc. to DIN EN 455
- protective glove according to clean room class ISO 4
- with extra long shaft (at least 29 cm)
- high wall thickness due to 20% more latex content than conventional latex gloves
- packed in a practical dispenser bag for additional protection against dirt
- skin friendly: powder-free, very low in proteins and low in accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit



Application

In all areas where the user is exposed to a risk of contamination. Ideal for all emergency situations, rescue services and intensive care. Suitable for highly infectious situations in hospitals and for private practice. Can be used in the care sector, in laboratories, in the pharmaceutical industry, in cleaning and in the industrial sector.

non-sterile

Size	REF	Box contents	Shipping unit	PZN
S	1229S	100 pcs	10 boxes	06966643
M	1229M	100 pcs	10 boxes	06966666
L	1229L	100 pcs	10 boxes	06966672

Gentle Skin[®] HiRisk

Our robust disposable latex glove

Characteristics

- examination glove made of latex acc. to DIN EN 455
- protective glove according to PPE
- high wall thickness for intensive protection
- skin friendly: powder-free, very low in proteins and low in accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit



Application

The high material quality provides maximum protection during heavy duty activities that require a thicker and longer glove. Especially popular and suitable for work in laboratories. Intensive protection for cytostatic preparation and handling of cleaning agents and disinfectants. Ideal for cleaning, including rinsing, washing, gardening and animal care.

non-sterile

Size	REF	Box contents	Shipping unit	PZN
S	1228S	50 pcs	10 boxes	11362924
M	1228M	50 pcs	10 boxes	11362930
L	1228L	50 pcs	10 boxes	11362947
XL	1228XL	50 pcs	10 boxes	11362953

Gentle Skin[®] Aloecare[™]

Our latex glove with aloe vera coating

Characteristics

- examination glove made of latex according to DIN EN 455
- with aloe vera containing inner coating
- skin friendly: powder-free, very low in proteins and low in accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fit / excellent fit



Application

Especially for users whose hands are additionally stressed by daily work with medical gloves over several hours and by regular use of disinfectants containing alcohol.

Aloe Vera inner coating - for the sake of your hands

The special inner coating has proven to offer the following advantages for users and especially stressed skin:

- cools
- moisturises
- noticeably relieves the typical itching and burning on rough and chapped skin
- contains more than 200 biologically active substances, including the following properties: anti-inflammatory, wound healing, anti-arthritis, analgesic

non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1231XS	100 pcs	10 boxes	05030230
S	1231S	100 pcs	10 boxes	02244516
M	1231M	100 pcs	10 boxes	02244580
L	1231L	100 pcs	10 boxes	02244597
XL	1231XL	100 pcs	10 boxes	05030247

Gentle Skin[®] black

Our black latex glove

Characteristics

- examination glove made of latex according to DIN EN 455
- skin-friendly: powder-free, very low in proteins and accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit
- suitable for contact with food



Application

The ideal latex glove, especially for hairdressers, beauticians, piercers and tattoo artists. Very popular in catering, trend chefs and the food industry as well as cleaning in general.

non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1224XS	100 pcs	10 boxes	11352676
S	1224S	100 pcs	10 boxes	11352682
M	1224M	100 pcs	10 boxes	11352736
L	1224L	100 pcs	10 boxes	11352771
XL	1224XL	100 pcs	10 boxes	11352802

Nitril® 3000

The white nitrile examination and protective glove - a resistant second skin

Characteristics

The Nitril® 3000 is a non-latex nitrile butadiene rubber (NBR) disposable glove. The material is not only very durable, but also adapts to the hand after a short time feeling like a second skin and thus guarantees fatigue-free working. The polymer inner coating makes it easy to put them on. Ambidextrous, with rolled edge, white, powder-free, micro-roughened, non-sterile, for single use.



non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1280XS	100 pcs	10 boxes	06966703
S	1280S	100 pcs	10 boxes	06967967
M	1280M	100 pcs	10 boxes	06968056
L	1280L	100 pcs	10 boxes	06968300
XL	1280XL	100 pcs	10 boxes	06968381

Application

Medical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). Suitable for cleaning agents and disinfection and for other applications in the laboratory, in the pharmaceutical industry, for professional cleaning, in catering and for handling food in general. Ideal for users with sensitive skin with diagnosed latex allergy (Type I).

Nitril® NextGen®

The extremely stretchable nitrile examination and protective glove

Characteristics

The Nitril® NextGen® is a non-latex nitrile butadiene rubber (NBR) disposable glove. Our "Innovation-Soft-Formula" gives this nitrile glove an almost latex-like, high elasticity and tear resistance. The polymer inner coating makes it easy to put them on. Ambidextrous, with rolled edge, blue, powder-free, micro-rough, non-sterile, for single use.



non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1283XS	100 pcs	10 boxes	09102239
S	1283S	100 pcs	10 boxes	09102245
M	1283M	100 pcs	10 boxes	09102251
L	1283L	100 pcs	10 boxes	09102268
XL	1283XL	100 pcs	10 boxes	09102274

Application

Medical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against chemicals and micro-organisms). Suitable for cleaning agents and disinfection and for other applications in the laboratory, in the pharmaceutical industry, for professional cleaning, in catering and for handling food in general. Ideal for users with sensitive skin with diagnosed latex allergy (Type I).

Nitril[®] Sensory[®]

Our new tactile and compact examination glove made of nitrile for your personal protection

Characteristics

The latex-free medical protective gloves made of nitrile butadiene rubber (NBR) are not only very elastic and tear resistant, but they will feel like a second skin in no time and thus guarantee fatigue-free working. The polymer inner coating allows effortless putting on. Ambidextrous, with rolled edge, white or blue, powder-free, micro-roughened, non-sterile, packed in 200 pieces, for single use.



Application

Medical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). Suitable for cleaning agents and disinfection and for other applications in the laboratory, in the pharmaceutical industry, for professional cleaning, in catering and for handling food in general. Ideal for users with sensitive skin with diagnosed latex allergy (Type I).

non-sterile

	Size	REF	Box contents	Shipping unit	PZN
blue	XS	2283XS	200 pcs	10 boxes	14310075
	S	2283S	200 pcs	10 boxes	14272305
	M	2283M	200 pcs	10 boxes	14272311
	L	2283L	200 pcs	10 boxes	14272328
	XL	2283XL	200 pcs	10 boxes	14272334
white	XS	2280XS	200 pcs	10 boxes	14306783
	S	2280S	200 pcs	10 boxes	14272268
	M	2280M	200 pcs	10 boxes	14272274
	L	2280L	200 pcs	10 boxes	14272280
	XL	2280XL	200 pcs	10 boxes	14272297

Nitril® BestGen®

The highly elastic nitrile examination and protective glove

Characteristics

The Nitril® BestGen® is a latex-free disposable nitrile rubber (NBR) glove. With our "Touch feeling formula", this type of glove surpasses most other products with excellent touch and gives it an almost latex-like high elasticity and tear resistance. The polymer inner coating makes them easy to put on. Ambidextrous with roll edge, blue, powder-free, micro-rough, non-sterile, for single use.



Application

Medical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374. (Protection against dangerous chemicals and micro-organisms). Suitable for cleaning agents and disinfection and ideal for other applications in the laboratory, in the pharmaceutical industry, for professional cleaning, in catering and for handling food in general. Ideal for users with sensitive skin with diagnosed latex allergy (Type I).

non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1286XS	100 pcs	10 boxes	10134612
S	1286S	100 pcs	10 boxes	10134641
M	1286M	100 pcs	10 boxes	10129545
L	1286L	100 pcs	10 boxes	10134658
XL	1286XL	100 pcs	10 boxes	10134670
XXL	1286XXL	100 pcs	10 boxes	11349987

Nitril[®] 3000 x-long blue

Our extra long nitrile glove

Characteristics

- examination glove made of nitrile according to DIN EN 455
- protective glove according to PPE
- with extra long shaft (about 29 cm) for optimum protection
- skin friendly: powder, protein and latex-free, low in accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit
- suitable for contact with food



Application

Due to its extra-long shaft, this glove not only provides reliable (spray) protection for the hands, but also for the forearms when handling irritating fluids. Therefore, it is ideal for handling chemicals, cytostatics, detergents and disinfectants. It is also used in the laboratory, the pharmaceutical industry and in catering. Ideal for latex allergy sufferers with diagnosed type I allergy.

non-sterile

Size	REF	Box contents	Shipping unit	PZN
S	1282S	100 pcs	10 boxes	00420825
M	1282M	100 pcs	10 boxes	00420802
L	1282L	100 pcs	10 boxes	00400857
XL	1282XL	100 pcs	10 boxes	00400834

Nitril[®] black

Our black nitrile examination and protective glove

Characteristics

The Nitril[®] black is a latex-free disposable glove made of nitrile butadiene rubber (NBR). The material is not only very durable, but also adapts to the hand after a short time feeling like a second skin and thus guarantees fatigue-free working. The polymer inner coating makes it easy to put them on. Ambidextrous, with rolled edge, black, powder-free, micro-roughened, non-sterile, for single use.



Application

Medical glove (according to DIN EN 455) and PPE (Personal Protective Equipment) tested according to EN 374 (protection against chemicals and micro-organisms). Suitable for cleaning agents and disinfection and for other applications in the laboratory, in the pharmaceutical industry, for professional cleaning, in catering and for handling food in general. Ideal for users with sensitive skin with diagnosed latex allergy (Type I).

non-sterile

Size	REF	Box contents	Shipping unit	PZN	Aid no.
S	1284S	100 pcs	10 boxes	12614730	15.25.03.1456
M	1284M	100 pcs	10 boxes	12614747	-
L	1284L	100 pcs	10 boxes	12614753	-
XL	1284XL	100 pcs	10 boxes	12614782	-

Nitril[®] Magenta

Our eye-catcher with caring silk sericin coating

Characteristics

The Nitril[®] Magenta is a latex-free examination protective glove made of nitrile butadiene rubber (NBR). The patented inner coating cares for the hands and makes them easy and effortless to put on. Ambidextrous, with rolled edge, pink, powder-free, micro-rough, non-sterile, for single use.



Application

Medical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). Ideal for cosmetic sectors. Suitable for handling cleaning agents and disinfectants, for applications in the laboratory, in the pharmaceutical industry, in catering and for handling food in general.

Silk sericin coating - for the sake of your hands

The special interior coating has proven to offer the following advantages for the user and particularly stressed skin:

- anti-aging effect
- natural humidification
- antimicrobial
- UV protection
- reduces redness, irritation and skin damage
- promotes wound healing
- collagen enrichment
- hydrating effect
- nourishes dry and sensitive skin
- firming effect

non-sterile

Size	REF	Box contents	Shipping unit	PZN
S	1287S	100 pcs	10 boxes	11584731
M	1287M	100 pcs	10 boxes	11584694
L	1287L	100 pcs	10 boxes	11584783




ANTIMICROBIAL
GLOVES™

Nitril LIOX[®]

Our first antimicrobial examination glove

Characteristics

The Nitril LIOX[®] is a non-latex, single use examination glove made from nitrile butadiene rubber. In contrast to conventional examination gloves, that only provide a passive barrier between micro-organisms and the hand, the antimicrobial effect of Nitril LIOX[®] plays an active role in reducing infection (hospital) transfer. Ambidextrous fitting, with rolled edges, blue violet, powderfree, textured, nonsterile, for single use only.



Application

Universally usable examination and protection glove for use in all areas, where the user is exposed to a contamination danger. Ideal for all emergency, paramedical and intensive care situations. Suitable for highly infectious situations in hospitals, as well as in general practice settings.

Advantages:

- Active protection against a multitude of hospital-acquired infections
- Kills up to 99,999% of selected micro-organisms
- Skin-friendliness has been confirmed via a multitude of tests

Mechanism of Action

An active pigment, imbedded on the outside of the glove, transforms naturally occurring oxygen into so called "Singlet Oxygen" as soon as the glove comes into contact with both oxygen and light. The reactive singlet oxygen molecules attack the cell membranes of gram-positive bacteria (MRSA, VRE, etc), thereby killing the affected microbes. Singlet Oxygen is already commonly used in medical settings. Amongst others, in the treatment of cancer using photodynamic therapy (PDT), due to impair tumor growth and pathological tissue alteration.



Size	REF	Box contents	Shipping Unit	PZN
XS	2284XS	200 pcs	10 boxes	15377756
S	2284S	200 pcs	10 boxes	15377762
M	2284M	200 pcs	10 boxes	15377779
L	2284L	200 pcs	10 boxes	15377785
XL	2284XL	180 pcs	10 boxes	15377816

Vinyl 2000

Our disposable vinyl examination glove

Characteristics

The Vinyl 2000 is a latex-free examination glove which thanks to its soft consistency and high tactile sensitivity is very popular among users. As a latex replacement, vinyl is ideal for users and patients with diagnosed latex (type I) and accelerator allergy (type IV). Ambidextrous, with rolled edge, non-sterile, for single use.



Application

Suitable for daily care of the elderly in nursing homes and hospitals, in private practice, in the nursing field in general and for handling non-fatty foods. Ideal for use in industrial environments. Due to its high functionality and its economic advantages it can be used where the consumption of disposable gloves is very high.

powdered, non-sterile

Size	REF	Box contents	Shipping unit	PZN
XS	1151XS	100 pcs	10 boxes	02244002
S	1151S	100 pcs	10 boxes	06177691
M	1151M	100 pcs	10 boxes	06177716
L	1151L	100 pcs	10 boxes	06177722
XL	1151XL	100 pcs	10 boxes	06177739

powder-free, non-sterile

XS	1251XS	100 pcs	10 boxes	02243907
S	1251S	100 pcs	10 boxes	02243913
M	1251M	100 pcs	10 boxes	02243936
L	1251L	100 pcs	10 boxes	02243942
XL	1251XL	100 pcs	10 boxes	02243994

Vinyl 2000 Stretch PF

Our highly elastic vinyl examination glove

Characteristics

The Vinyl 2000 Stretch PF has significantly better elasticity and density compared to traditional vinyl gloves, providing the user safer work for daily needs. As a latex substitute, vinyl is ideal for users and patients with diagnosed Latex-(type I) and accelerator allergy (type IV). Ambidextrous, with rolled edge, powder-free, non-sterile, for single use.



Application

The reasonably priced latex-free and comfortable alternative for care. Ideal for use in industrial environments. Due to its high functionality and its economic advantages it can be used where the consumption of disposable gloves is very high.

non-sterile

Size	REF	Box contents	Shipping unit	PZN
S	1253S	100 pcs	10 boxes	09544664
M	1253M	100 pcs	10 boxes	09544670
L	1253L	100 pcs	10 boxes	09544687
XL	1253XL	100 pcs	10 boxes	09544724

Cotton gloves

Our soft glove made of pure cotton

Characteristics

- glove made of pure cotton
- skin friendly protection for the hands
- permeable to air
- hard wearing due to firmly processed seams
- anatomical fit



Application

Ideal for the dermatological sector, especially when treated with very oily or slow absorbing ointments, e.g. in atopic dermatitis, very dry skin or dermatitis and to cover gauze dressings (do not use directly on open wounds!). Also suitable for the reprocessing of cleaned surgical instruments. They protect the instruments from the natural acidity and fat content of the hands.

non-sterile

Größe	REF	Box contents	Shipping unit	PZN
Gr. 6 / XS	3611	10 pairs	30 boxes	04354686
Gr. 7 / XS	3612	10 pairs	30 boxes	04354692
Gr. 8 / S	3613	10 pairs	30 boxes	04354700
Gr. 9 / S	3614	10 pairs	30 boxes	04354717
Gr. 10 / M	3619	10 pairs	30 boxes	04354568
Gr. 11 / M	3620	10 pairs	30 boxes	04354574
Gr. 12 / L	3621	10 pairs	30 boxes	04354580
Gr. 13 / L	3622	10 pairs	30 boxes	04354597
Gr. 14 / XL	3623	10 pairs	30 boxes	04354657
Gr. 15 / XL	3624	10 pairs	30 boxes	04354663

Froggy™

Our hygiene glove made of soft polyethylene

Characteristics

- transparent hygienic glove made of PE
- skin friendly: powder, protein and latex-free
- in practical glove dispenser for easy removal
- stable welded
- tear resistant due to firm hems on the sides
- ambidextrous
- in women's and men's size



non-sterile

Size	REF	Box contents	Shipping unit	PZN
Women	1192S	100 pcs	100 boxes	06177745
Men	1192L	100 pcs	100 boxes	06177751

Application

For economical use in simple activities. Perfect for protecting hands against rough external soiling and for all areas requiring sanitary standards. Ideal for chemists, supermarkets, bakeries, the food industry and petrol stations. **WARNING:** Only suitable for single use when fuelling the car with diesel fuel. Short term protection against minor, superficial diesel pollution (this PE disposable glove is not Personal protective equipment (PPE) under Regulation (EU) 2016/425).

Élastify®

Our low cost hygienic glove made of super tear resistant TPE

Characteristics

- transparent hygienic glove made of TPE (thermoplastic elastomer), powder-free
- skin friendly: free of plasticizers (phthalates), proteins and latex
- in practical glove dispenser for easy removal
- breathable and pleasant wearing comfort due to special material composition
- stable and extremely tear resistant
- ambidextrous suitable
- for single use
- suitable for contact with food



non-sterile

Size	REF	Box contents	Shipping unit	PZN
S	1193S	200 pcs	10 boxes	11362893
M	1193M	200 pcs	10 boxes	11362901
L	1193L	200 pcs	10 boxes	11362918

Application

Perfect for protecting hands against rough external soiling and for all areas requiring sanitary standards. Also for more demanding activities. Ideal for the food industry, chemists, supermarkets, bakeries, butchers and industrial sectors in general. Popular and economical alternative to vinyl gloves.

Latex finger cots

Non-sterile disposable product with rolled edge and smooth surface

Characteristics

- CE according to MPG - Directive 93/42 / EEC + 2007/47 / EC - Class 1 product
- declaration of safety for short term contact with food
- applicable quality assurance: EN ISO 13485:2012 + AC: 2012
- white finger cots made of natural latex, powder, silicone and thiuram-free
- only suitable for minimal risks
- not for activities with dangerous, sharp edged and aggressive substances



Application

To cover individual fingers for the administration of suppositories, application of creams and simple examinations. Use mostly in examination, nursing, hygiene, laboratory and industrial areas.

non-sterile

Size	REF	Box contents	Shipping unit	PZN
S	9776	100 pcs	1 pack	01333784
M	9777	100 pcs	1 pack	01333790
L	9778	100 pcs	1 pack	01333809
XL	9779	100 pcs	1 pack	01333815

MyClean® easy touch

The reusable glove
made of latex

Characteristics

- Reusable glove made of latex, powder-free
- Fully anatomically shaped
- High wearing comfort due to cotton inner flocking
- Special grip profile due to honeycomb structure on palm
- Suitable for food



Anwendung

Particularly suitable for cleaning work, rinsing, washing, gardening or animal care. Due to its suitability for use in foodstuffs, it may also be used in gastronomy, kitchens and food processing.

non-sterile

Size	REF	Box contents	Shipping unit
S	MAM74152	1 pair	100 boxes
M	MAM74153	1 pair	100 boxes
L	MAM74154	1 pair	100 boxes
XL	MAM76286	1 pair	100 boxes

Copolymed®

Our sterile copolymer examination glove

Characteristics

- sterile examination glove made of copolymer
- skin friendly: powder, protein and latex-free
- good tactile sensation
- tear resistant
- ambidextrous fitting / excellent fit
- punched on paper



Application

Ideal for short medical interventions, for sterile activities in the daily work of surgical areas and wards, as well as in industry. Performs the function of an under and over glove when additional dirt and liquid protection is required. Suitable for allergy sufferers with diagnosed type I and / or type IV allergy.

sterile, individually packed

Size	REF	Box contents	Shipping unit	PZN
S	8092S	100 pcs	12 boxes	03896416
M	8092M	100 pcs	12 boxes	03896221
L	8092L	100 pcs	12 boxes	03896215

sterile, packed in pairs

S	8091S	50 pairs	12 boxes	11362864
M	8091M	50 pairs	12 boxes	11362870
L	8091L	50 pairs	12 boxes	11362887

Gentle Skin[®] steril

Our sterile latex examination glove

Characteristics

- sterile latex examination glove according to DIN EN 455
- optimum grip
- skin-friendly: powder-free, very low in proteins and accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit
- packed individually or in pairs



Application

Suitable for all surgical areas. Wherever the user works in a sterile environment. Ideal for minor surgical procedures, dressing changes, punctures, etc...

sterile, individually packed

Size	REF	Box contents	Shipping unit	PZN
S	8022S	50 pcs	8 boxes	06402811
M	8022M	50 pcs	8 boxes	06402828
L	8022L	50 pcs	8 boxes	06402834

sterile, packed in pairs

S	8021RSW	50 pairs	4 boxes	02244479
M	8021RMW	50 pairs	4 boxes	02244485
L	8021RLW	50 pairs	4 boxes	02244491

Nitril® steril

Our sterile nitrile examination glove

Characteristics

- sterile examination glove made of nitrile according to DIN EN 455
- skin friendly: powder, protein and latex-free, low in accelerators
- best fitting
- micro-roughened for sensitive and excellent tactile sensation
- highly elastic and tear resistant
- ambidextrous fitting / excellent fit



Application

Can be used in all surgical areas. Wherever the user works in a sterile environment.

Ideal for minor surgical procedures, dressing changes, punctures, etc... Suitable for latex allergy sufferers with diagnosed type I allergy. Protective glove according to PPE.

sterile, packed in pairs

Size	REF	Box contents	Shipping unit	PZN
S	8023S	50 pairs	4 boxes	09544658
M	8023M	50 pairs	4 boxes	09544641
L	8023L	50 pairs	4 boxes	09544635
XL	8023XL	50 pairs	4 boxes	11362976

Our surgical gloves set standards in terms of protection and safety

In operations, the protection of medical personnel and patients comes first. The Meditrade surgical glove range includes a diverse selection of different glove models made of different materials, all of which are characterised by mechanical strength, tightness, grip, comfort and biocompatibility. The gloves are continuously optimised in their properties and further developed to meet the highest demands in the market.



Reference™ OP

The powdered latex surgical glove

Characteristics

The Reference™ OP is a powdered surgical glove made of a highly adaptable and very stretchable latex structure with optimum grip and micro-roughening, which gives the glove a very high comfort and, especially at the fingertips - a sensitive tactile sensation.

The optimum fit allows long term users in the operative field to work fatigue free for several hours. The light powder coating guarantees that the gloves can be put on and taken off with ease. Sterile, powdered, micro-roughened, anatomically shaped, with rolled edge, for single use.

Application

Surgical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). Suitable in all surgical areas with sterile requirements and wherever the user moves in a sterile environment and detects a risk of contamination. Ideal for highly infectious situations in hospitals, in private practice, in laboratories or in the pharmaceutical industry. Also suitable for double-gloving.



sterile, packed in pairs

Size	REF	Box contents	Shipping unit	PZN
6	93216	50 pairs	6 boxes	02243824
6,5	932165	50 pairs	6 boxes	02243830
7	93217	50 pairs	6 boxes	02243847
7,5	932175	50 pairs	6 boxes	02243853
8	93218	50 pairs	6 boxes	02243876
8,5	932185	50 pairs	6 boxes	02243882
9	93219	50 pairs	6 boxes	05030313



Gentle Skin[®] Premium OP[™]

The premium quality latex surgical glove

Characteristics

The Gentle Skin[®] Premium OP[™] is a surgical glove made of a highly adaptable and very stretchable latex structure with optimum grip and micro-roughening, which gives the glove a very high comfort and, especially at the fingertips - a sensitive tactile sensation. The special polymer coating forms an additional barrier between skin and latex and noticeably improves putting on. The easy handling and excellent product quality to increase the feeling of safety considerably, make the Gentle Skin[®] Premium OP[™] rightly one of the most popular powder-free surgical gloves on the market. Sterile, powder-free, micro-roughened, anatomically shaped, with rolled edge, for single use.

Application

Surgical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). Particularly suitable for vascular surgery, orthopaedics, abdominal surgery and in all surgical areas with sterile requirements. Wherever the user works in a sterile environment and detects a risk of contamination. Ideal for all highly infectious situations in hospitals, in private practice, in laboratories or in the pharmaceutical industry. Also suitable for double-gloving.



sterile, packed in pairs

Size	REF	Box contents	Shipping unit	PZN
5,5	902155W	50 pairs	4 boxes	05030299
6	90216W	50 pairs	4 boxes	02244114
6,5	902165W	50 pairs	4 boxes	02244278
7	90217W	50 pairs	4 boxes	02244290
7,5	902175W	50 pairs	4 boxes	02244315
8	90218W	50 pairs	4 boxes	02244338
8,5	902185W	50 pairs	4 boxes	02244427
9	90219W	50 pairs	4 boxes	02244433



Gentle Skin[®] Smooth OP

The latex surgical glove with a smooth texture

Characteristics

The Gentle Skin[®] Smooth OP is a surgical glove with a very supple and stretchy latex structure, which gives the glove a very high wearing comfort and sensitive tactile sensation. The mentioned properties coupled with the anatomical fit guarantee the user fatigue-free and safe work in the surgical daily routine. The special polymer coating forms an additional barrier between skin and latex and noticeably improves their putting on. Very good grip for safe instrument operation. Sterile, powder-free, anatomically shaped, with rolled edge, for single use.

Application

Surgical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). Particularly suitable for vascular surgery, orthopaedics, abdominal surgery and in all surgical areas with sterile requirements. Wherever the user works in a sterile environment and detects a risk of contamination. Ideal for all highly infectious situations in hospitals, in private practice, in laboratories or in the pharmaceutical industry. Also suitable for double-gloving.



sterile, packed in pairs

Size	REF	Box contents	Shipping unit	PZN
6	90266	50 pairs	4 boxes	11154775
6,5	902665	50 pairs	4 boxes	11154781
7	90267	50 pairs	4 boxes	11154798
7,5	902675	50 pairs	4 boxes	11154806
8	90268	50 pairs	4 boxes	11154812
8,5	902685	50 pairs	4 boxes	11154829
9	90269	50 pairs	4 boxes	11154835



Gentle Skin[®] Micro OP[®]

The touch sensitive latex surgical glove

Characteristics

The Gentle Skin[®] Micro OP[®] was specially developed for fine surgical applications. It has optimum touch sensitivity and is slightly thinner when compared to other surgical gloves. Despite the thinner wall thickness, the highest level of safety is guaranteed by the fulfilment of all standards and norms.

The colour and the exceptionally good fit of the glove allow fatigue-free working, especially in the field of ophthalmology. Sterile, powder-free, micro-roughened, anatomically shaped, with rolled edge, for single use.

Application

Surgical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). Surgical glove for the highest demands in areas where particularly sensitive tactile sensation is required.

Therefore particularly suitable for use in micro-surgery, neurology and ophthalmology. Generally in all surgical areas with sterile requirements and wherever the user moves in a sterile environment and recognizes a risk of contamination. Ideal for all highly infectious situations in hospitals, in



private practice, in laboratories or in the pharmaceutical industry. Also suitable for double-gloving.

sterile, packed in pairs

Size	REF	Box contents	Shipping unit	PZN
6	90316	50 pairs	4 boxes	06972460
6,5	903165	50 pairs	4 boxes	06972589
7	90317	50 pairs	4 boxes	06972595
7,5	903175	50 pairs	4 boxes	06972626
8	90318	50 pairs	4 boxes	06972655
8,5	903185	50 pairs	4 boxes	06972661



Gentle Skin® Securitex® OP

The surgical latex glove set with indicator function

Characteristics

The Gentle Skin®-Securitex® set consists of 2 pairs of surgical gloves with indicator function and is therefore ideal for surgical procedures. The special colour indicator function and the reliable fluid indication of our glove set offer double safety for the user. The inner glove has a striking green signal colour meaning the surgeon is immediately warned if the outer glove is damaged. If the damage is not visible to the naked eye, the special liquid indication offers additional security. As: If the glove is damaged and gets in contact with liquid of any kind, the perforated area will immediately turn deep green. The outer glove is a bit thicker, to minimise the risk of such damage in advance. It is brownish to bring out the indicator function as effectively as possible by means of high contrast. Sterile, powder-free, micro-roughened, anatomically shaped, with rolled edge, for single use.


Application

Surgical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). Ideally suited for surgical and highly infectious interventions in which a particularly high level of protection is required, in hospitals, in private practice, in laboratories or in the pharmaceutical industry. Ideal for all sterile requiring surgical areas and wherever the user moves in a sterile environment and detects a risk of contamination.



sterile set, packed as 2 pairs

Size	REF	Box contents	Shipping unit	PZN
6	90516	25 sets	4 boxes	11349970
6,5	905165	25 sets	4 boxes	03851224
7	90517	25 sets	4 boxes	03851230
7,5	905175	25 sets	4 boxes	03851247
8	90518	25 sets	4 boxes	03851253
8,5	905185	25 sets	4 boxes	03851276



Gentle Skin® Isopretex® OP

The latex-free surgical glove made of fine polyisoprene

Characteristics

The Gentle Skin® Isopretex® OP is a latex-free surgical glove made of synthetic rubber. This glove is similar in quality to a powder-free latex surgical glove and offers an excellent alternative for latex allergy sufferers. A special formula gives the user an incredibly high degree of flexibility, softness and comfort.

The properties mentioned coupled with the anatomical fit guarantee the user fatigue-free and safe work in the day to day of the operating theatre. Its outstanding elasticity, suppleness and perfectly balanced grip make the Gentle Skin® Isopretex® OP the perfect glove for lengthy procedures.

A unique inner coating guarantees perfect putting on, both in dry and wet conditions. Sterile, powder-free, micro-roughened, anatomically shaped, with rolled edge, for single use.

Application

Surgical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). Due to its latex-free composition especially suitable for emergency and rescue services and especially for the health care of children. Generally in all surgical areas with sterile requirements and wherever the user moves in a sterile environment and detects a risk of contamination. Ideal for all highly infectious situations in hospitals, in private practice, in laboratories or in



the pharmaceutical industry. Also suitable for double-gloving.

sterile, packed in pairs

Size	REF	Box contents	Shipping unit	PZN
5,5	952255	50 pairs	4 boxes	11190854
6	95226	50 pairs	4 boxes	11190860
6,5	952265	50 pairs	4 boxes	11190877
7	95227	50 pairs	4 boxes	11190883
7,5	952275	50 pairs	4 boxes	11190908
8	95228	50 pairs	4 boxes	11190914
8,5	952285	50 pairs	4 boxes	11190920
9	95229	50 pairs	4 boxes	11190937



Neopretex® OP

The latex-free surgical glove made of flexible polychloroprene

Characteristics

The Neopretex® OP is a latex-free surgical glove made of polychloroprene, also known as Neoprene. Due to its properties as an "artificial substitute" for natural rubber latex, its highly elastic formula makes it a suitable alternative for latex allergy sufferers. The material is not only very durable, but also adapts to the hand after a short time feeling like a second skin and thus guarantees fatigue-free working. Sterile, powder-free, micro-roughened, anatomically shaped, with rolled edge, for single use.

Application

Surgical glove (according to DIN EN 455) and PPE (Personal protective equipment) tested according to EN 374 (protection against hazardous chemicals and micro-organisms). As they are latex-free, they are especially suitable for emergency and rescue services and above all for the health care of children. Generally in all surgical areas with sterile requirements and wherever the user works in a sterile environment and detects a risk of contamination. Ideal for all highly infectious situations in hospitals, in private practice, in laboratories or in the pharmaceutical industry.



Also suitable for double-gloving.

sterile, packed in pairs

Size	REF	Box contents	Shipping unit	PZN
6	95216	50 pairs	4 boxes	02426890
6,5	952165	50 pairs	4 boxes	02426909
7	95217	50 pairs	4 boxes	02426915
7,5	952175	50 pairs	4 boxes	02426921
8	95218	50 pairs	4 boxes	02426938
8,5	952185	50 pairs	4 boxes	02426944
9	95219	50 pairs	4 boxes	02426950

Finessis®

An innovative surgical glove program that reinvents

CLEANNES, COMFORT and REACTIVE SAFETY.




Made of **FLEXYLON™** - The technology of the future


Flexylon™ is a compact poured film without any accelerators or crosslinkers.

This technique eliminates almost all puncture holes and other micro-defects and thus achieves the highest quality standards (AQL = 0.1)!

Crosslinking technology



Flexylon™
Molecular
Layer technology



Finessis® Zero®

Performance & safety
General surgery

FLEXYLON™ I Without accelerators I Synthetic surgical gloves

Finessis® Zero® is the next generation of surgical gloves, made without any accelerators and thus achieving new levels of PURITY. Using FLEXYLON™ layer technology eliminates the risk of Type I skin allergy and restricts Type IV skin allergy to an unprecedented level.

Special features:

- WITHOUT accelerators
- Best AQL in the sector (AQL = 0.10)
- Made from FLEXYLON™ (layer technology)
- High electrical resistance
- Low stress thanks to optimum anatomical fit

Size	REF	Box contents	Shipping unit
6	90276	50 pairs	4 boxes
6,5	902765	50 pairs	4 boxes
7	90277	50 pairs	4 boxes
7,5	902775	50 pairs	4 boxes
8	90278	50 pairs	4 boxes
8,5	902785	50 pairs	4 boxes
9	90279	50 pairs	4 boxes

COLOUR

White



Finessis®

Finessis® Aegis®

Reactive protection

**Trauma & high risk
Orthopaedics**



FLEXYLON™ | Without accelerators | Synthetic surgical gloves

Aegis® is a revolutionary protective system. It is a surgical glove equipped with a reaction initiated sterilisation system. The Finessis® Aegis® was designed inside the glove with an anti-microbial layer made from high performance elastomers between the sleeve. The disinfection system will only come into effect IF the glove has been PERFORATED. In this case, the activated gel is released in the APPROPRIATE QUANTITY and at the LOCATION where the perforation originated.

Special features:

- Reaction-initiated sterilisation system
- WITHOUT accelerators
- Antimicrobial interlayer that releases activated gel in case of perforation
- Best AQL in the sector (AQL = 0.10)
- Made from FLEXYLON™ (layer technology)
- High electrical resistance
- Low stress thanks to optimum anatomical fit

COLOUR


Light green / white



Size	REF	Box contents	Shipping unit
6	90286	40 pairs	4 boxes
6,5	902865	40 pairs	4 boxes
7	90287	40 pairs	4 boxes
7,5	902875	40 pairs	4 boxes
8	90288	40 pairs	4 boxes
8,5	902885	40 pairs	4 boxes
9	90289	40 pairs	4 boxes

Finessis®

Reactive protection



Reactive anti-microbial intermediate layer

Outer layer made of FLEXYLON™

Inner layer of FLEXYLON™




Figure 1: A sharp object presses on the outer layer without breaking it.



Figure 2: Under increasing pressure, the thin membranes tear, that separate the drops of the gel. A larger contiguous amount of anti-microbial, gel-like fluid is formed.






Figure 3: Once the sharp object has penetrated, it is enveloped by the solution and viruses or bacteria are effectively sterilised.



Finessis® Corium®

Micro-surgery & underglove



FLEXYLON™ I Without accelerators I Synthetic surgical gloves


Finessis® Corium® is the glove that guarantees the highest tactility for double-gloving and micro-surgery. It is designed to provide improved tactile sensitivity. Corium® is perfect as an underglove for maximum comfort without losing any of its grip. Combined with the Finessis® Zero® for a perfect fit, the Finessis® Double gloving system provides an impermeable barrier to any surgical application.

Special features:

- WITHOUT accelerators
- Best AQL in the sector (AQL = 0.10)
- Made with FLEXYLON™ (layer technology)
- Smooth surface
- High wearing comfort with improved tactility
- High electrical resistance
- Low load thanks to optimum anatomical fit

COLOUR

● Green



Size	REF	Box contents	Shipping unit
6	90246	50 pairs	4 boxes
6,5	902465	50 pairs	4 boxes
7	90247	50 pairs	4 boxes
7,5	902475	50 pairs	4 boxes
8	90248	50 pairs	4 boxes
8,5	902485	50 pairs	4 boxes
9	90249	50 pairs	4 boxes

Finessis®

X-Guard

Our radiation protection glove with shielding effect

Characteristics

- sterile, lead-enriched latex glove according to DIN EN 455
- protective glove according to PPE
- with shielding effect against X-ray scattering radiation
- high wall thickness for intensive protection
- pleasant to wear
- excellent tactile sensation
- elastic and tear resistant
- anatomically shaped



Application

For all medical interventions with X-ray support as a protection against scattered radiation and for all situations in which the user is exposed to ionizing radiation. For interventional cardiology, radiology, urology, gastroenterology, neuroradiology, orthopaedics, nuclear medicine, pain therapy, catheterisation, endoscopy, intensive care and electrophysiology.

sterile, packed in pairs, wall thickness 0.22 ± 0.02 mm

Size	REF	Box contents	Shipping unit	PZN
6	RR1-119.00.01	5 pairs	10 boxes	11607515
6,5	RR1-119.00.02	5 pairs	10 boxes	11607521
7	RR1-119.00.03	5 pairs	10 boxes	11607538
7,5	RR1-119.00.04	5 pairs	10 boxes	11607544
8	RR1-119.00.05	5 pairs	10 boxes	11607550
8,5	RR1-119.00.06	5 pairs	10 boxes	11607567
9	RR1-119.00.07	5 pairs	10 boxes	11607573

sterile, packed in pairs, wall thickness 0.30 ± 0.02 mm

6	RR2-119.00.01	5 pairs	10 boxes	11607596
6,5	RR2-119.00.02	5 pairs	10 boxes	11607604
7	RR2-119.00.03	5 pairs	10 boxes	11607610
7,5	RR2-119.00.04	5 pairs	10 boxes	11607627
8	RR2-119.00.05	5 pairs	10 boxes	11607633
8	RR2-119.00.06	5 pairs	10 boxes	11607662
9	RR2-119.00.07	5 pairs	10 boxes	11607679



More safety in the operating theatre and on the ward



With our products, we make an important contribution to preventing infection paths and germ transfer. Personnel, patients and visitors should feel safe and comfortable in dealing with each other. Wearing comfort, breathability and individual adjustment are important aspects of contemporary and suitable clothing.


Especially within the field of PPE (Personal protective equipment) it is not enough to rely on industry standardised products, but rather high standards of medical quality must be maintained.

For example, when making Suavel face masks, we use the method of Ultrasonic sealing to avoid any separation of the thread.

We always use materials that define the current state of the art, and therefore are consistently able to offer the best product.

More safety in the operating theatre and on the ward

Characteristics of all Suavel® Surgical masks



Suavel[®] Protec Plus

Our high filtration surgical mask with elastic ear loops

Characteristics

- surgical mask according to EN 14683, type II
- bacterial filter efficiency $\geq 99\%$
- 3 ply
- elastic ear loops
- green
- for single use
- real surgical fleece for a long wearing period
- increased surface for better protection
- very skin friendly
- high wearing comfort



Application

For standard and long term applications and all applications in the operating theatre and on the ward for the protection of users and patients. Guaranteed high wearing comfort even with prolonged use.

non-sterile

REF	Colour	Box contents	Shipping unit	PZN
80-955MP	green	50 pcs	18 boxes	02244605

Suavel[®] Protec

Our proven surgical mask with elastic ear loops

Characteristics

- surgical mask according to EN 14683, type II
- bacterial filter efficiency $\geq 98\%$
- 3 ply
- elastic ear loops
- green
- for single use



Application

For standard and long term applications and all applications in the operating theatre and on the ward for the protection of users and patients. Guaranteed high wearing comfort even with prolonged use.

non-sterile

REF	Colour	Box contents	Shipping unit	PZN
80-900	green	50 pcs	18 boxes	00068771

Suavel[®] Comfort Plus Suavel[®] Comfort

Our high filtration surgical mask for tying

Characteristics

- surgical mask according to EN, 14683 type II
- bacterial filter efficiency $\geq 99\%$
- 3 ply
- for tying
- extra long straps
- blue
- for single use
- real surgical fleece for a long wearing period
- increased surface for better protection
- very skin friendly
- high wearing comfort



Application

For standard and long term applications and all applications in the operating theatre and on the ward for the protection of users and patients. Guaranteed high wearing comfort even with prolonged use.

non-sterile

REF	Colour	Box contents	Shipping unit	PZN
80-410M	blue	50 pcs	12 boxes	02244611

Our proven surgical mask for tying

Characteristics

- surgical mask according to EN 14683, type II
- bacterial filter efficiency $\geq 98\%$
- 3 ply
- for tying
- extra long straps
- blue
- for single use



Application

For standard and long term applications and all applications in the operating theatre and on the ward for the protection of users and patients. Guaranteed high wearing comfort even with prolonged use.

non-sterile

REF	Colour	Box contents	Shipping unit	PZN
80-400	blue	50 pcs	12 boxes	00068788

Suavel[®] Antifog

Our surgical mask for wearers of glasses

Characteristics

- with effective glasses fitting insulation
- Surgical mask according to EN 14683, type II
- Bacterial filter efficiency $\geq 98\%$
- 3 ply
- for tying
- extra long straps
- blue
- for single use
- real surgical fleece for a long wearing period
- increased surface for better protection
- very skin friendly
- high wearing comfort



Effective glasses fog insulation



Application

Especially for wearers of glasses. For standard and long-term applications and all applications in the operating theatre and on the ward for the protection of users and patients. Guaranteed high wearing comfort even with prolonged use.

non-sterile

REF	Colour	Box contents	Shipping unit	PZN
80-470M	blue	50 pcs	12 boxes	00559457

Suavel[®] Antifluid

Our surgical mask with protective visor for highly infectious situations

Characteristics

- the highest protection against infection through integrated protective visor for mouth, eyes and nose
- liquid refracting layer
- Surgical mask according to EN 14683, Type II R
- Bacterial filter efficiency $\geq 98\%$
- 4 ply
- for tying
- extra long straps
- blue
- for single use
- real surgical fleece for a long wearing period
- increased surface for better protection
- very skin friendly
- high wearing comfort



The highest protection against infection thanks to integrated protective visor for mouth, eyes and nose, liquid refracting layer



Application

For highly infectious situations, protects users against blood and/or wound secretion splashes as well as during rinsing. Ideal for standard and long term applications and all applications in the operating theatre and on the ward. Guaranteed high wearing comfort even with prolonged use.

non-sterile

REF	Colour	Box contents	Shipping unit	PZN
80-455M	blue	50 pcs	12 boxes	00559517

Suavel[®] Sensima

Our surgical mask for allergy sufferers

Characteristics

- hypo-allergenic
- surgical mask according to EN 14683, type II for sensitive users
- bacterial filter efficiency $\geq 98\%$
- 3 ply
- for tying
- extra long straps
- white
- for single use
- real surgical fleece for a long wearing period
- increased surface for better protection
- high wearing comfort



Hypoallergenes
Lagenmaterial

Application

Ideal for allergy sufferers and sensitive users. For standard and long term applications and all applications in the operating theatre and on the ward for the protection of users and patients. Guaranteed high wearing comfort even with prolonged use.

non-sterile

REF	Colour	Box contents	Shipping unit	PZN
80-440M	white	50 pcs	12 boxes	11362700

Suavel® Apollo

Our astronaut shape surgical hood

Characteristics

- high quality surgical astronaut shape hood
- made of tear resistant polypropylene non-woven fabric
- with extra long straps and neck protection
- highest wearing comfort with good fit
- breathable for optimum heat regulation
- with integrated sweatband
- hypo-allergenic

Application

For longer operations to cover the hair and reduce perspiration. Especially suitable for beards.



non-sterile

REF	Colour	Box contents	Shipping unit	PZN
70-131G	green	100 pcs	10 boxes	11362597

Suavel® Astrid

Our universal bouffant cap

Characteristics

- high quality bouffant cap
- made of soft and lightweight polypropylene non-woven material
- with double stitched rubber spring
- highest wearing comfort with good fit
- breathable for the best possible heat regulation
- hypo-allergenic



Application

Ideal for use in hospitals and the care area for nurses and patients.

Wherever hair coverage is required, including sectors such as food and catering.

Suavel® Astrid, non-sterile

REF	Size	Colour	Box contents	Shipping unit	PZN
70-121B	53 cm	blue	100 pcs	10 boxes	11362611
70-121G	53 cm	green	100 pcs	10 boxes	11362628
70-121W	53 cm	white	100 pcs	10 boxes	11362605

Suavel® Astrid XL, non-sterile

70-124G	60 cm	green	100 pcs	10 boxes	11362634
70-124W	60 cm	white	100 pcs	10 boxes	11362640



BeeSana[®]

PP protective gown 35 g

The high performance

protective gown for personnel and visitors

Characteristics

Protective gown made of polypropylene non-woven material in 35 g/m² thickness. The material is easy to wear, breathable and thus ensures the best temperature compensation. Tie straps at the neck and at the waist and cuffs as sleeve ends for the highest wearing comfort. Fluid repellent, extremely durable, non-sterile, for single use.

Application

Ward and visitor gowns to protect against germs in the hospital and in the care of the elderly for personnel, patients and visitors.



non-sterile

REF	Size	Colour	Box contents	Shipping unit	PZN
3686	140 x 140 cm	white	10 pcs	10 boxes	09747456
3700	140 x 140 cm	green	10 pcs	10 boxes	09747462

BeeSana[®]

PP gown 23 g

The classic

protective gown for personnel and visitors

Characteristics

Protective gown made of polypropylene non-woven material in 23 g/m² thickness. Closure bands at the neck and waist. Sleeve ending with sewn-in, elastic rubber straps or cuffs. Lightweight, breathable material, very tear resistant and at the same time moisture repellent, non-sterile, for single use.

Application

Ward and visitor gowns to protect against germs in the hospital and in the care of the elderly for personnel, patients and visitors.



non-sterile, with elastic rubber straps

REF	Size	Colour	Box contents	Shipping unit	PZN
3684	110 x 140 cm	white	10 pcs	10 boxes	11509391
3685	110 x 140 cm	green	10 pcs	10 boxes	11509416
3699	120 x 140 cm	yellow	10 pcs	10 boxes	09747479
3729	130 x 150 cm	white	10 pcs	10 boxes	11527609

non-sterile, with cuffs

4685	110 x 140 cm	green	10 pcs	10 boxes	11363036
4699	120 x 140 cm	yellow	10 pcs	10 boxes	11363042

BeeSana® PP/PE gown 26 g

The waterproof

partially laminated protective gown for personnel and visitors

Characteristics

Protective gown made of PP with PE coating in the front and arm area in 26 g/m² thickness. Closure bands at the neck and waist. Cuffs as sleeve ends. Breathable material, very tear resistant and at the same time waterproof in the coated area, non-sterile, for single use.

The BeeSana® PP/PE gown was successfully tested according to EN 14126 as protective clothing against infectious agents.

Application

Ward and visitor gowns for applications with high safety requirements to protect against germs in the hospital and in the care of the elderly for personnel, patients and visitors.



non-sterile

REF	Size	Colour	Box contents	Shipping unit	PZN
4698	136 x 136 cm	yellow	10 pcs	10 boxes	11642842

BeeSana® SMS gown 40 g

Maximum protection

Protective gown for extra protection against liquids and particles

Characteristics

Protective gown made of 3 ply, soft SMS composite material. The combined layers guarantee a high barrier against liquids, germs and light chemical splashes. The gown is also opaque, lint-free, tear resistant and breathable at the same time for optimum comfort. Neck and waist fastening bands as well as jersey cuffs as sleeve ends, non-sterile, for single use.

The BeeSana® SMS gown 40 g was successfully tested for its resistance to germ penetration in moist conditions (based on DIN EN ISO 22610) and in dry conditions (based on DIN EN ISO 22612) and on its virus impermeability (according to DIN EN 14126) and is thus perfectly suitable for handling MRSA or noroviruses.

Application

Ideal for highly infectious situations. To protect against germs in emergency and isolation wards as well as in the care of the elderly and in the hospital for personnel, patients and visitors.



non-sterile

REF	Size	Colour	Box contents	Shipping unit	PZN
3597	115 x 150 cm	blue	10 pcs	5 boxes	11652154

BeeSana[®]

SMS gown 15 g

The economical effective protective gown for personnel and visitors

Characteristics

Protective gown made of 3 ply, soft SMS composite material. The combined layers guarantee a high barrier against liquids, germs and light chemical splashes. The gown is also opaque, lint-free, tear resistant and breathable at the same time for optimum comfort. Fastening bands at the neck and waist and cuffs as sleeve ends, non-sterile, for single use.

Application

Ideal for highly infectious situations. To protect against germs in emergency and isolation wards as well as in the care of the elderly and in the hospital for personnel, patients and visitors.



BeeSana[®]

SSS gown 9 g

The ultra-light economical protective gown for personnel and visitors

Characteristics

Protective gown made of 3 ply, soft SSS composite material in thickness 9 g/m². By using the most innovative materials, it has been possible to create a highly economical product with similar properties of protective gowns from significantly higher weight classes. The gown is lint-free and breathable for optimum comfort. Sleeve ending with sewn-in, elastic rubber bands. Fastening bands at neck and waist, non-sterile, for single use.

Application

Ward and visitor gowns to protect against germs in the hospital and in the care of the elderly for personnel, patients and visitors.



non-sterile

REF	Size	Colour	Box contents	Shipping unit	PZN
3598	140 x 140 cm	blue	10 pcs	10 boxes	11518680

non-sterile

REF	Size	Colour	Box contents	Shipping unit	PZN
3681	110 x 140 cm	green	10 pcs	10 boxes	11518674

PE disposable apron

Our PE apron in various models

Characteristics

- disposable apron made of LDPE
- tested according to EN 14126
- with fastening bands at the waist for optimum fit
- pleasant wearing comfort due to lightweight material
- hard wearing and moisture impermeable
- in practical bag or blocked for hanging



Application

Optimum protection of clothing during nursing activities and in all areas with high hygiene requirements, including in commercial kitchens and hospitals.

non-sterile, in a bag, for single removal

REF	Size	Thickness	Colour	Box contents	Shipping unit	PZN
3591	140 x 75 cm	30 µ	blue	50 pcs	10 boxes	11857643
3690	125 x 81 cm	16 µ	white	100 pcs	10 boxes	04353379
3691	140 x 81 cm	16 µ	white	100 pcs	10 boxes	04353422
3692	160 x 81 cm	16 µ	white	100 pcs	10 boxes	04353474

non-sterile, blocked, for tearing off individually

3687	125 x 81 cm	16 µ	white	50 pcs	10 boxes	04353385
3688	140 x 81 cm	16 µ	white	50 pcs	10 boxes	04353468
3694	160 x 81 cm	16 µ	white	50 pcs	10 boxes	04353480



Wound treatment and dressings - the supporting pillars of medical care

Our claim is always to offer the optimum product for the respective wound or the current wound state, since a wound dressing alone can not fulfil all therapeutic requirements.

These include all non-sterile and sterile gauze and non-woven swabs as well as suction and ointment swabs. In the field of therapeutic dressings, we have all application relevant products to be able to offer the right product for compressing, supporting and relieving.

To provide even greater patient safety in the surgical area all our sterile products are equipped with the highest safety standards. Whether it be woven X-ray contrast strips, counting control cards or the highest quality control with proof of tracking from the procurement of raw materials to the finished product in double peel packaging.

BeeSana® Lap sponge x-ray thread

Our lap sponge in many variants

Characteristics

- pre-washed abdominal gauze made of medical gauze acc. EN 14079
- with X-ray contrast chip
- 20 thread
- soft and extremely absorbent
- specified formats: Measurements BEFORE washing



Application

For absorbing and stopping blood and body fluids during surgical procedures.

For a variety of purposes, the abdominal dressings are available in different sizes, layers and colours.

non-sterile, 5 pcs each, bundled with sleeve

		REF	Size	Loops	Box contents	Shipping unit	PZN	
4 ply	white	1613	20 x 30 cm	Yes	4 x 5 pcs	36 boxes	13506907	
		1547	40 x 40 cm	Yes	4 x 5 pcs	18 boxes	13506913	
		1616	45 x 45 cm	Yes	4 x 5 pcs	12 boxes	13506936	
		1648	45 x 45 cm	No	4 x 5 pcs	12 boxes	13506942	
		1670	50 x 60 cm	Yes	4 x 5 pcs	8 boxes	13506959	
	unwashed	1602	50 x 60 cm	No	4 x 5 pcs	10 boxes	13506965	
		1465	8 x 90 cm	Yes	4 x 5 pcs	30 boxes	13506994	
	green	1614	20 x 30 cm	Yes	4 x 5 pcs	36 boxes	13507002	
		1617	45 x 45 cm	Yes	4 x 5 pcs	12 boxes	13507019	
		1607	45 x 45 cm	No	4 x 5 pcs	12 boxes	13507025	
1696		8 x 90 cm	Yes	4 x 5 pcs	30 boxes	13507031		
	blue	4593	45 x 45 cm	Yes	4 x 5 pcs	12 boxes	13507048	
6 ply	white	1612	20 x 30 cm	Yes	4 x 5 pcs	30 boxes	13507054	
		1502	40 x 40 cm	Yes	4 x 5 pcs	15 boxes	13507083	
		1618	45 x 45 cm	Yes	4 x 5 pcs	10 boxes	13507108	
	unwashed	1520	45 x 45 cm	Yes	4 x 5 pcs	10 boxes	13507120	
		1659	45 x 45 cm	No	4 x 5 pcs	10 boxes	13507137	
		1671	50 x 60 cm	Yes	4 x 5 pcs	6 boxes	13507143	
		1449	8 x 90 cm	Yes	4 x 5 pcs	27 boxes	13507166	
			green	1615	20 x 30 cm	Yes	4 x 5 pcs	30 boxes
	24 ply	1496	20 x 30 cm	Yes	4 x 5 pcs	10 boxes	13507189	
		1608	45 x 45 cm	Yes	4 x 5 pcs	10 boxes	13507203	
		1651	45 x 45 cm	No	4 x 5 pcs	10 boxes	13507226	
	24 ply	1583	50 x 60 cm	Yes	2 x 5 pcs	10 boxes	13507232	
		1603	50 x 60 cm	No	4 x 5 pcs	6 boxes	13507249	
		4058	8 x 90 cm	Yes	4 x 5 pcs	27 boxes	13507255	
	blue	4592	45 x 45 cm	Yes	4 x 5 pcs	10 boxes	13507261	
4594		45 x 45 cm	No	4 x 5 pcs	10 boxes	13507278		
8 ply	white	1506	45 x 45 cm	Yes	2 x 5 pcs	10 boxes	13507284	



double sterile, with self-adhesive counting cards

				REF	Size	Loops	Box contents	Shipping unit	PZN	
4 ply	in pack of 2	white	double peel	1664	45 x 45 cm	Yes	15 x 2 pcs	4 boxes	13166417	
			double peel	4515	50 x 60 cm	Yes	24 x 2 pcs	2 boxes	13166423	
			double peel	1682	8 x 90 cm	Yes	30 x 2 pcs	4 boxes	13166446	
		green	double peel	1665	45 x 45 cm	Yes	15 x 2 pcs	4 boxes	13166452	
			double peel	4504	50 x 60 cm	Yes	24 x 2 pcs	2 boxes	13166469	
			double peel	1684	8 x 90 cm	Yes	30 x 2 pcs	4 boxes	13166475	
		in pack of 5	white	double peel	1571	20 x 30 cm	Yes	15 x 5 pcs	4 boxes	13614399
				double peel	1584	40 x 40 cm	Yes	15 x 5 pcs	2 boxes	13166481
				double peel	1593	45 x 45 cm	Yes	15 x 5 pcs	2 boxes	13166498
unwashed, 24 ply	green	double peel	1597	30 x 30 cm	No	18 x 5 pcs	2 boxes	13166297		
		double peel	1572	20 x 30 cm	Yes	15 x 5 pcs	4 boxes	13166280		
		double peel	1594	45 x 45 cm	Yes	15 x 5 pcs	2 boxes	13166334		
	in pack of 40	white	wrapping cloth	1624	40 x 40 cm	Yes	3 x 40 pcs	2 boxes	13166340	
6 ply	in pack of 2	white	double peel	1666	45 x 45 cm	Yes	15 x 2 pcs	4 boxes	13166357	
			double peel	1677	45 x 45 cm	No	15 x 2 pcs	4 boxes	13166363	
			double peel	1685	8 x 90 cm	Yes	30 x 2 pcs	4 boxes	14218665	
	in pack of 2	green	double peel	1667	45 x 45 cm	Yes	15 x 2 pcs	4 boxes	13166386	
			double peel	4501	50 x 60 cm	Yes	20 x 2 pcs	2 boxes	13166185	
			double peel	1680	8 x 90 cm	Yes	30 x 2 pcs	4 boxes	13166191	
	in pack of 5	white	double peel	1586	40 x 40 cm	Yes	12 x 5 pcs	2 boxes	13166216	
			double peel	1595	45 x 45 cm	Yes	12 x 5 pcs	2 boxes	13166222	
			double peel	1591	45 x 45 cm	No	12 x 5 pcs	2 boxes	13166239	
			double peel	1669	8 x 90 cm	Yes	24 x 5 pcs	2 boxes	13166245	
	in pack of 5	green	double peel	1574	20 x 30 cm	Yes	9 x 5 pcs	4 boxes	13166251	
			double peel	1587	40 x 40 cm	Yes	12 x 5 pcs	2 boxes	13166268	
			double peel	1596	45 x 45 cm	Yes	12 x 5 pcs	2 boxes	13166400	
			wrapping cloth	1695	45 x 45 cm	No	12 x 5 pcs	2 boxes	13166392	

BeeSana® Gauze ball

Our absorbent gauze ball

Characteristics

- gauze ball made from medical gauze to EN 14079
- innovative wrapping system from one piece of gauze cutting
- cutting edges on the inside for maximum safety
- 20 thread count
- very absorbent



Application

For wound cleansing and haemostasis as well as for extensive skin disinfection.

Due to the different swab sizes, large and small wounds can be cared for efficiently and quickly.

non-sterile, in a polybag

	REF	Size	Shape	Box contents	Shipping unit	PZN
hazelnut size	1102	10 x 10 cm	round	500 pcs	8 boxes	01364017
walnut size	4460	15 x 16 cm	round	500 pcs	2 boxes	00822216
plum size	4461	20 x 20 cm	round	250 pcs	16 boxes	00822239
plum size	1106	20 x 20 cm	Little hat	250 pcs	16 boxes	00822222
egg size	4462	25 x 26 cm	round	250 pcs	10 boxes	00822251
egg size	1109	25 x 26 cm	Little hat	250 pcs	10 boxes	00822245
X large	1348	30 x 30 cm	round	100 pcs	20 boxes	11862578
XX large	1113	30 x 40 cm	round	100 pcs	15 boxes	01363940
fist size	1116	40 x 50 cm	round	100 pcs	10 boxes	00822268



sterile, in peel pouch

in pack of 2	plum size	4356	20 x 20 cm	round	100 x 2 pcs	8 boxes	13506801
in pack of 3	plum size	4357	20 x 20 cm	round	80 x 3 pcs	8 boxes	00137696
	egg size	4308	25 x 26 cm	round	60 x 3 pcs	8 boxes	13506818
	fist size	4085	45 x 45 cm	round	35 x 3 pcs	4 boxes	13506824
in pack of 2 + 2 packed in soft peel pouch with 2 cells (perforation to tear of)	walnut size	1074	15 x 16 cm	round	45 x (2 +3) pcs	8 boxes	13506830
	plum size	1084	20 x 20 cm	round	50 x (2+2) pcs	8 boxes	01364129
	plum size	1164	20 x 20 cm	Little hat	50 x (2+2) pcs	8 boxes	01364106
in pack of 2+3 packed in soft peel pouch with 2 cells (perforation to tear of)	plum size	1131	20 x 20 cm	round	40 x (2+3) pcs	8 boxes	01364135
	plum size	1189	20 x 20 cm	Little hat	40 x (2+3) pcs	8 boxes	01364075
	egg size	1141	25 x 26 cm	round	35 x (2+3) pcs	8 boxes	-
	XX large	4216	30 x 40 cm	round	35 x (2+3) pcs	4 boxes	13506847
in pack of 5	plum size	1184	20 x 20 cm	round	50 x 5 pcs	8 boxes	13506853
	egg size	4209	25 x 26 cm	round	40 x 5 pcs	8 boxes	13506876
	egg size	1165	25 x 26 cm	Little hat	40 x 5 pcs	8 boxes	01364069
	XX large	4211	30 x 40 cm	round	40 x 5 pcs	4 boxes	13506882
in pack of 10	plum size	4215	20 x 20 cm	round	30 x 10 pcs	10 boxes	07393706
in pack of 10	egg size	4219	25 x 26 cm	round	40 x 10 pcs	4 boxes	-

BeeSana® Gauze ball x-ray thread

Characteristics

- gauze ball made from medical gauze to EN 14079
- woven X-ray contrast thread (BaSO₄)
- innovative wrapping system from one piece of gauze cutting
- cutting edges on the inside for maximum safety
- 20 thread count
- very absorbent
- for secure documentation in the surgical area



non-sterile, in a polybag

	REF	Size	Shape	Box contents	Shipping unit	PZN
plum size	1154	20 x 20 cm	round	250 pcs	16 boxes	00822297
egg size	1155	25 x 26 cm	round	250 pcs	10 boxes	00822311
egg size	1110	25 x 26 cm	Little hat	250 pcs	10 boxes	00822305
X large	1349	30 x 30 cm	round	100 pcs	20 boxes	13508243
XX large	1114	30 x 40 cm	round	100 pcs	15 boxes	01363868
XXX large	1115	40 x 40 cm	round	100 pcs	10 boxes	13508266
fist size	1117	40 x 50 cm	round	100 pcs	10 boxes	00822328
dabber swab	1119	60 x 80 cm	round	50 pcs	8 boxes	13508295



double sterile, with self-adhesive counting cards

in pack of 5	egg size	4324	25 x 26 cm	round	25 x 5 pcs	8 boxes	13165145
	XX large	4326	30 x 40 cm	round	30 x 5 pcs	4 boxes	13165151
	fist size	4325	40 x 50 cm	round	20 x 5 pcs	4 boxes	13165168
in pack of 10	walnut size	1180	15 x 16 cm	round	30 x 10 pcs	8 boxes	13165174
	plum size	1181	20 x 20 cm	round	18 x 10 pcs	8 boxes	13165180
	egg size	1075	25 x 26 cm	round	30 x 10 pcs	4 boxes	11703503
	egg size	1186	25 x 26 cm	Little hat	30 x 10 pcs	4 boxes	13165197
	X large	1182	30 x 30 cm	round	20 x 10 pcs	4 boxes	13165211
	XX large	4319	30 x 40 cm	round	18 x 10 pcs	4 boxes	13165292
	fist size	4333	40 x 50 cm	round	12 x 10 pcs	4 boxes	13165228
in pack of 20	plum size	4329	20 x 20 cm	round	25 x 20 pcs	4 boxes	13165234
	egg size	4337	25 x 26 cm	round	18 x 20 pcs	4 boxes	13165240
	egg size	4012	25 x 26 cm	Little hat	18 x 20 pcs	4 boxes	13165257
	X large	4362	30 x 30 cm	round	15 x 20 pcs	4 boxes	-
	XX large	4320	30 x 40 cm	round	20 x 20 pcs	2 boxes	13165263
in pack of 30	egg size	1159	25 x 26 cm	Little hat	12 x 30 pcs	4 boxes	13231480

BeeSana[®] Thread swab

Characteristics

- gauze swabs according to EN 14079
- with green return thread
- innovative wrapping system from one piece of gauze cutting
- cutting edges on the inside for maximum safety
- 20 thread count
- very absorbent

Application

For blood and secretion absorption.

Safe retrieval of the swab by means of the thread.



non-sterile, without x-ray thread

	REF	Size	Shape	Thread	Box contents	Shipping unit	PZN
egg size	4015	25 x 26 cm	round	single	100 pcs	6 boxes	13578278

non-sterile, with x-ray thread

	REF	Size	Shape	Thread	Box contents	Shipping unit	PZN
plum size	4043	20 x 20 cm	round	double	200 pcs	5 boxes	12390911

BeeSana® Preparation balls x-ray thread

Characteristics

- gauze ball made of 100% cotton according to EN 14079
- woven X-ray contrast thread (BaSO₄)
- 24 ply
- innovative wrapping system from one piece of gauze cutting
- high absorption capacity



Application

For free preparation of tissue layers and anatomical structures.

The X-ray contrast thread makes these swabs X-Ray detectable and meets the required safety standard.

non-sterile

	REF	Size	Box contents	Shipping unit	PZN
in a polybag	1169	4 x 4 cm	100 pcs	50 boxes	13583440
	1170	6 x 6 cm	100 pcs	40 boxes	01363839
	1171	8 x 8 cm	100 pcs	30 boxes	01363822
	1173	10 x 10 cm	100 pcs	20 boxes	12386275
	1172	12 x 12 cm	100 pcs	15 boxes	01363845
pre-counted to 10 pcs in paper bag for set packer	1170P10	6 x 6 cm	100 x 10 pcs	10 boxes	13583457
	1171P10	8 x 8 cm	100 x 10 pcs	8 boxes	13583463
	1172P10	12 x 12 cm	100 x 10 pcs	6 boxes	13583486
pre-counted to 10 pcs in counting box for set packer	1169B10	4 x 4 cm	146 x 10 pcs	1 carton	-
	1170B10	6 x 6 cm	146 x 10 pcs	1 carton	-
	1171B10	8 x 8 cm	146 x 10 pcs	1 carton	-
	1173B10	10 x 10 cm	148 x 10 pcs	1 carton	-
	1172B10	12 x 12 cm	148 x 10 pcs	1 carton	-



double sterile, with self-adhesive counting cards

in pack of 10 pcs in paper bag	4089	6 x 6 cm	60 x 10 pcs	2 boxes	13165731
	4311	8 x 8 cm	60 x 10 pcs	2 boxes	13165754
	4339	12 x 12 cm	50 x 10 pcs	2 boxes	13165777
in pack of 20 pcs in paper bag	4331	6 x 6 cm	60 x 20 pcs	2 boxes	13165748
	4335	8 x 8 cm	60 x 20 pcs	2 boxes	13165760
in pack of 10 pcs in counting box	4343	4 x 4 cm	20 x 10 pcs	9 boxes	13585373
	4340	6 x 6 cm	20 x 10 pcs	9 boxes	12731105
	4341	8 x 8 cm	20 x 10 pcs	9 boxes	12731111
	4342	12 x 12 cm	10 x 10 pcs	4 boxes	12731128

BeeSana[®] Pointed swab

Characteristics

- gauze ball made of 100% cotton according to EN 14079
- 20 thread count
- conically shaped from just one piece of gauze
- high absorption capacity



Application

For blood and secretion, especially in the ENT area. For haemostasis or as a tamponade.

without x-ray thread, double sterile with counting cards

	REF	Size (unrolled)	Length	Box contents	Shipping unit	PZN
in pack of 10	4170	large, 15 x 16 cm	approx. 13.5 cm	40 x 10 pcs	8 boxes	-

BeeSana[®] Pointed swab x-ray thread

Characteristics

- gauze ball made of 100% cotton according to EN 14079
- 20 thread count
- conically shaped from just one piece of gauze
- high absorption capacity
- woven X-ray contrast thread (BaSO₄)



Application

For blood and secretion, especially in the ENT area. For haemostasis or as a tamponade.

non-sterile

REF	Size (unrolled)	Length	Box contents	Shipping unit	PZN
1238	large, 15 x 16 cm	approx. 13.5 cm	500 pcs	20 boxes	13509403



double sterile, with self-adhesive counting cards

	REF	Size (unrolled)	Length	Box contents	Shipping unit	PZN
in pack of 10	4004	small, 6 x 8 cm	approx. 8 cm	60 x 10 pcs	8 boxes	13165783
	4161	large, 15 x 16 cm	approx. 13.5 cm	40 x 10 pcs	8 boxes	13165808

BeeSana[®] Cellulose swabs

Our cellulose swab off the roll

Characteristics

- swabs of high bleached, absorbent and non-fibrous cellulose according to DAB
- ready to use
- pre-punched



Application

For local disinfection before injections and infusions.

non-sterile

	REF	Size	Box contents	Shipping unit	PZN
	1190	4 x 5 cm	2 rolls of 500 pcs	16 boxes	07210048
	1278	8 x 10 cm	2 rolls of 125 pcs	10 boxes	01364253
reduced number of germs	1177	4 x 5 cm	1 roll of 500 pcs	20 boxes	01331070
Cellulose swab box	9871		1 pc	1 pack	01363756

Alcohol swabs

Our ready-to-use swab for skin cleansing

Characteristics

- non-woven alcohol swabs
- soaked in 70% isopropyl alcohol
- ready to use

Application

To cleanse the skin before subcutaneous injections and capillary blood collection.



non-sterile

REF	Size	Box contents	Shipping unit	PZN
4465	approx. 3 x 3 cm (folded) approx. 6 x 3 cm (unfolded)	100 pcs	20 boxes	09391893

ABE[®] Ribbon gauze

Our pure cotton ribbon gauze

Characteristics

- medical gauze according to EN 14079
- 100 % cotton
- woven edges
- white
- rolled
- woven X-ray contrast thread(BaSO₄)



Application

As a tamponade in natural body openings.

sterile, pack of 1

	REF	Size	Box contents	Shipping unit	PZN
24 ply	4300	1 cm x 5 m	100 x 1 pc	8 boxes	13574814
	4301	2 cm x 5 m	7 x 1 pc	8 boxes	13574820

ABE[®] Ribbon gauze x-ray thread

Our pure cotton ribbon gauze

Characteristics

- medical gauze according to EN 14079
- 100 % cotton
- woven edges
- 20 thread count
- 4 x ply
- white
- rolled
- woven X-ray contrast thread (BaSO₄)



Application

Tamponing of wound cavities or fistula courses. Absorbing body fluids in hard to reach wound and body cavities.

non-sterile

	REF	Size	Box contents	Shipping unit	PZN
20 thread count	1446	6 cm x 2 m	40 pcs	10 boxes	12386298

double sterile, 1 pc. with counting cards

	REF	Size	Box contents	Shipping unit	PZN
20 thread count	4312	6 cm x 2 m	25 x 1 pc	4 boxes	02214231

BeeSana[®] Gauze cutting

Our medical gauze for individual wound dressings

Characteristics

- gauze cutting according to EN 14079
- 20 thread count
- very absorbent
- breathable and tear-resistant at the same time
- wound friendly and pleasant on the skin



Application

For the individually prepared wound dressing. For external wound care of secreting and / or contaminated wounds. For blood absorption and secretion.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1314	6 x 6 cm	8000 pcs	1 carton	13719804
1330	20 x 20 cm	500 pcs	10 boxes	-
1335	30 x 40 cm	500 pcs	6 boxes	-

BeeSana[®] Medical gauze

Characteristics

- medical gauze according to EN 14079
- 20 ply

Application

For making individual dressings.

non-sterile

80 cm wide, 1 ply, 4 ply laid

REF	Size	Box contents	Shipping unit	PZN
1301	80 cm x 40 m	5 pcs	2 boxes	04852674

double sterile, with self-adhesive counting cards

80 cm wide, 5 ply, placed on 16 x 20 cm, packed as 3 pcs sterile

REF	Length	Box contents	Shipping unit	PZN
4179	100 cm	8 x 3 pcs	2 boxes	13581004



BeeSana[®] Gauze swab

Our classic gauze with selvedge cutting edges

Characteristics

- gauze swab made from medical gauze according to EN 14079
- selvedge cutting edges, 100% cotton
- folded, 17 thread count
- very absorbent
- not fibrous
- breathable and tear resistant at the same time
- wound friendly and pleasant on the skin



Application

For general wound care of secreting wounds and absorption of body fluids. As a carrier for ointments. For minor surgery in the ambulance, on the ward and in the first aid area.

non-sterile

	REF	Size	Box contents	Shipping unit	PZN	
8 ply	1000	5 x 5 cm	100 pcs	10 boxes	04852705	
	1006	7.5 x 7.5 cm	100 pcs	10 boxes	04852711	
	1012	10 x 10 cm	100 pcs	10 boxes	04852728	
	1053	10 x 20 cm	100 pcs	5 boxes	01330633	
pre-counted to 10 pcs for set packer	2712	10 x 10 cm	200 x 10 pcs	1 carton	-	
12 ply	1002	5 x 5 cm	100 pcs	10 boxes	04852734	
	1008	7.5 x 7.5 cm	100 pcs	10 boxes	04852740	
	1014	10 x 10 cm	100 pcs	10 boxes	04852757	
	1038	10 x 13 cm	100 pcs	10 boxes	04852763	
	1023	10 x 20 cm	100 pcs	10 boxes	04852786	
	pre-counted to 10 pcs for set packer	2708	7.5 x 7.5 cm	200 x 10 pcs	1 carton	-
	2714	10 x 10 cm	200 x 10 pcs	1 carton	-	
2723	10 x 20 cm	100 x 10 pcs	1 carton	-		
16 ply	1004	5 x 5 cm	100 pcs	10 boxes	04852792	
	1010	7.5 x 7.5 cm	100 pcs	10 boxes	04852800	
	1016	10 x 10 cm	100 pcs	10 boxes	04852817	
	1052	10 x 13 cm	100 pcs	10 boxes	04852823	
	1025	10 x 20 cm	100 pcs	5 boxes	11849589	
pre-counted to 10 pcs for set packer	2716	10 x 10 cm	200 x 10 pcs	1 carton	-	
32 ply	1018	10 x 10 cm	50 pcs	10 boxes	-	
	1091	12.5 x 20 cm	50 pcs	4 boxes	-	
64 ply	20 thread count	1235	10 x 13 cm	25 pcs	8 boxes	11849626

BeeSana[®] Gauze swab

Our classic gauze with selvedge cutting edges



sterile

		REF	Size	Box contents	Shipping unit	PZN	
8 ply	in pack of 2	1043	5 x 5 cm	25 x 2 pcs	27 boxes	04852846	
		1044	7.5 x 7.5 cm	25 x 2 pcs	16 boxes	04852852	
		1045	10 x 10 cm	25 x 2 pcs	16 boxes	04852869	
		1046	10 x 20 cm	25 x 2 pcs	8 boxes	04852875	
	in pack of 5	4600	5 x 5 cm	80 x 5 pcs	8 boxes	13458367	
		4601	7.5 x 7.5 cm	80 x 5 pcs	8 boxes	01330662	
		1062	10 x 10 cm	35 x 5 pcs	6 boxes	01330679	
	in pack of 10	4606	10 x 10 cm	50 x 10 pcs	4 boxes	12348949	
	12 ply	in pack of 2	4130	5 x 5 cm	25 x 2 pcs	24 boxes	07216329
			4131	7.5 x 7.5 cm	25 x 2 pcs	18 boxes	07216335
4132			10 x 10 cm	25 x 2 pcs	12 boxes	07216341	
4133			10 x 20 cm	25 x 2 pcs	6 boxes	07216358	
in pack of 5		1064	5 x 5 cm	30 x 5 pcs	12 boxes	01330722	
		1065	7.5 x 7.5 cm	30 x 5 pcs	6 boxes	01330739	
		1063	10 x 10 cm	30 x 5 pcs	6 boxes	01330745	
		1066	10 x 20 cm	15 x 5 pcs	6 boxes	01330751	
in pack of 10		4616	5 x 5 cm	50 x 10 pcs	8 boxes	12348955	
		4617	7.5 x 7.5 cm	30 x 10 pcs	8 boxes	12349044	
		1026	10 x 10 cm	40 x 10 pcs	4 boxes	12348694	
		4619	10 x 20 cm	20 x 10 pcs	4 boxes	12349015	
in pack of 20		4622	10 x 10 cm	20 x 20 pcs	4 boxes	12349021	
in pack of 30		4640	7.5 x 7.5 cm	25 x 30 pcs	4 boxes	12349185	
16 ply		in pack of 4	1092	7.5 x 7.5 cm	64 x 4 pcs	8 boxes	12348783
		in pack of 10	4173	10 x 10 cm	30 x 10 pcs	4 boxes	13582564
	in pack of 40	4618	10 x 10 cm	18 x 40 pcs	2 boxes	12349009	
32 ply	in pack of 2	4628	10 x 10 cm	30 x 2 pcs	8 boxes	12349038	



sterile, rolled gauze swab

		REF	Size	Box contents	Shipping unit	PZN
16 ply	in pack of 2	1146	20 x 20 cm	20 x 2 pcs	4 boxes	13582587



double sterile, with self-adhesive counting cards

16 ply	in pack of 10	1058	10 x 10 cm	25 x 10 pcs	4 boxes	13582601
	in pack of 60	4644	10 x 10 cm	8 x 60 pcs	2 boxes	12349216
	in pack of 100	4645	5 x 5 cm	15 x 100 pcs	4 boxes	12349222
32 ply	in pack of 5	4636	12.5 x 20 cm	12 x 5 pcs	4 boxes	12349239

BeeSana® Gauze swab x-ray thread

Our classic with X-ray contrast thread

Characteristics

- gauze swabs made from medical gauze according to EN 14079
- woven X-ray contrast thread (BaSO₄)
- selvedge cutting edges
- folded, 17 thread count
- very absorbent
- breathable and tear resistant at the same time
- wound friendly and pleasant on the skin
- non-sterile in sterilisable paper packaging
- double sterile packed with self-adhesive counting cards
- for safe documentation in the surgical area, wound friendly and pleasant on the skin



Application

For surgical invasive use in the surgical area.
For haemostasis and wound care.

non-sterile

	REF	Size	Box contents	Shipping unit	PZN
8 ply	1013	10 x 10 cm	100 pcs	10 boxes	12348599
	1037	10 x 20 cm	100 pcs	5 boxes	12348719
12 ply	1003	5 x 5 cm	100 pcs	10 boxes	12348582
	1009	7.5 x 7.5 cm	100 pcs	10 boxes	06613703
	1015	10 x 10 cm	100 pcs	10 boxes	06613726
	1039	10 x 13 cm	100 pcs	10 boxes	12348731
	1024	10 x 20 cm	100 pcs	5 boxes	06613732
	pre-counted to 10 pcs for set packer	2715	10 x 10 cm	200 x 10 pcs	1 carton
	2724	10 x 20 cm	100 x 10 pcs	1 carton	-
16 ply	1005	5 x 5 cm	100 pcs	10 boxes	06613749
	1011	7.5 x 7.5 cm	100 pcs	10 boxes	06613755
	1017	10 x 10 cm	100 pcs	10 boxes	06613761
	1054	10 x 13 cm	100 pcs	10 boxes	12348665
	1034	10 x 20 cm	100 pcs	5 boxes	12348702
	pre-counted to 10 pcs for set packer	2717	10 x 10 cm	200 x 10 pcs	1 carton
	2727	20 x 20 cm	60 x 10 pcs	1 carton	-
24 ply	1090	10 x 20 cm	50 pcs	5 boxes	12348642
32 ply	1019	10 x 10 cm	50 pcs	10 boxes	12348607
	20 thread count	1021	10 x 10 cm	50 pcs	10 boxes
pre-counted to 10 pcs for set packer	2719	10 x 10 cm	100 x 10 pcs	1 carton	-



double sterile, with self-adhesive counting cards

		REF	Size	Box contents	Shipping unit	PZN
12 ply	in pack of 5	1149	10 x 10 cm	25 x 5 pcs	8 boxes	12348754
		1150	10 x 20 cm	30 x 5 pcs	4 boxes	12348760
	in pack of 10	4392	5 x 5 cm	30 x 10 pcs	8 boxes	12348932
		1174	7.5 x 7.5 cm	40 x 10 pcs	4 boxes	12348843
		1049	10 x 10 cm	30 x 10 pcs	4 boxes	11703472
		1168	10 x 20 cm	15 x 10 pcs	4 boxes	12348777
	in pack of 20	1176	7.5 x 7.5 cm	20 x 20 pcs	4 boxes	12348866
		4314	10 x 10 cm	15 x 20 pcs	4 boxes	12348820
		1041	10 x 20 cm	20 x 20 pcs	2 boxes	12348659
16 ply	in pack of 10	4368	5 x 5 cm	30 x 10 pcs	8 boxes	12348990
		1071	7.5 x 7.5 cm	30 x 10 pcs	4 boxes	12348613
		1072	10 x 10 cm	25 x 10 pcs	4 boxes	12348636
		1142	10 x 20 cm	15 x 10 pcs	4 boxes	12348808
	in pack of 20	4351	10 x 10 cm	12 x 20 pcs	4 boxes	12348984
	in pack of 40	4204	10 x 10 cm	15 x 40 pcs	2 boxes	12348814
32 ply	in pack of 10	4318	10 x 10 cm	15 x 10 pcs	4 boxes	12348837

BeeSana[®] Slit gauze swab

Our classic with Y-slit

Characteristics

- gauze swab made from medical gauze according to EN 14079
- Y-slit
- selvedge cutting edges
- folded, 17 thread count
- breathable and tear resistant at the same time
- wound friendly and pleasant on the skin



Application

For wound care in drainages, cannulas, tracheotomies and extensions.

non-sterile

	REF	Size	Box contents	Shipping unit	PZN
8 ply	4291	7.5 x 7.5 cm	100 pcs	10 boxes	12491084
12 ply	4192	5 x 5 cm	100 pcs	10 boxes	12491090
	4122	7.5 x 7.5 cm	100 pcs	10 boxes	12491109
	4160	10 x 10 cm	100 pcs	10 boxes	12491115
16 ply	4159	5 x 5 cm	100 pcs	10 boxes	12491121



sterile, pack of 2

8 ply	4158	10 x 10 cm	100 x 2 pcs	8 boxes	03870954
12 ply	4115	5 x 5 cm	25 x 2 pcs	16 boxes	01916669
	4116	7.5 x 7.5 cm	25 x 2 pcs	16 boxes	01916652
	4117	10 x 10 cm	25 x 2 pcs	16 boxes	01916646

BeeSana® Non-woven swab

Our non-woven swabs for external wound care, sterile and non-sterile

Characteristics

- non-woven viscose and polyester swab
- very absorbent
- breathable and tear resistant at the same time
- wound friendly and especially soft on the skin
- not fibrous



Application

For general wound care of secreting wounds and absorption of body fluids.

As a carrier for ointments for smaller, surgical interventions in the ambulance, on the ward and in the area of first aid.

non-sterile

		REF	Size	Box contents	Shipping unit	PZN	
4 ply	30 g	1789	5 x 5 cm	100 pcs	10 boxes	01450254	
		1790	7.5 x 7.5 cm	100 pcs	10 boxes	01450225	
		pre-counted to 10 pcs for set packer	2790	7.5 x 7.5 cm	200 x 10 pcs	1 carton	-
			1791	10 x 10 cm	100 pcs	10 boxes	01450231
		pre-counted to 10 pcs for set packer	2791	10 x 10 cm	200 x 10 pcs	1 carton	-
			1792	10 x 20 cm	100 pcs	5 boxes	01445773
	40 g		1124	5 x 5 cm	100 pcs	10 boxes	04852881
			1125	7.5 x 7.5 cm	100 pcs	10 boxes	04852898
			1126	10 x 10 cm	100 pcs	10 boxes	04852906
			large carton	11261	10 x 10 cm	100 pcs	76 boxes
		1127	10 x 20 cm	100 pcs	5 boxes	04852912	
6 ply	30 g	1796	5 x 5 cm	100 pcs	10 boxes	01330886	
		1797	7.5 x 7.5 cm	100 pcs	10 boxes	01330900	
		1798	10 x 10 cm	100 pcs	10 boxes	01330917	
		1799	10 x 20 cm	100 pcs	5 boxes	01330923	
			pre-counted to 10 pcs for set packer	2798	10 x 10 cm	200 x 10 pcs	1 carton
			2799	10 x 20 cm	100 x 10 pcs	1 carton	-

sterile

4 ply	30 g	in pack of 2	1137	5 x 5 cm	25 x 2 pcs	27 boxes	03036843	
			1138	7.5 x 7.5 cm	25 x 2 pcs	16 boxes	03036895	
			1139	10 x 10 cm	25 x 2 pcs	16 boxes	03036903	
			1140	10 x 20 cm	25 x 2 pcs	8 boxes	03036949	
		in pack of 5	4376	7.5 x 7.5 cm	100 x 5 pcs	8 boxes	12350768	
	4377		10 x 10 cm	60 x 5 pcs	8 boxes	12350774		
		in pack of 10	1219	10 x 10 cm	25 x 10 pcs	8 boxes	12350691	
	40 g	in pack of 100	4229	5 x 5 cm	30 x 100 pcs	4 boxes	12350716	
	6 ply	30 g	in pack of 2	1800	5 x 5 cm	25 x 2 pcs	27 boxes	01330834
				1801	7.5 x 7.5 cm	25 x 2 pcs	16 boxes	01330840
				1802	10 x 10 cm	25 x 2 pcs	16 boxes	01330857
				1803	10 x 20 cm	25 x 2 pcs	8 boxes	03036955
		in pack of 5	4225	5 x 5 cm	120 x 5 pcs	8 boxes	12350662	
4226			7.5 x 7.5 cm	70 x 5 pcs	8 boxes	12350679		
4227			10 x 10 cm	40 x 5 pcs	8 boxes	12350685		

BeeSana[®] Slit non-woven swab

Our slit non-woven swab with Y-slit

Characteristics

- non-woven viscose and polyester swab
- Y-slit
- very absorbent
- breathable and tear resistant at the same time
- wound friendly and especially soft on the skin
- not fibrous



Application

For wound care in drainages, cannulas, tracheotomies and extensions.

non-sterile

	REF	Size	Box contents	Shipping unit	PZN
4 ply, 30 g	1793	5 x 5 cm	100 pcs	10 boxes	01330768
	1794	7.5 x 7.5 cm	100 pcs	10 boxes	01330774
	1795	10 x 10 cm	100 pcs	10 boxes	01330780



sterile, pack of 2

4 ply, 30 g	1232	5 x 5 cm	25 x 2 pcs	27 boxes	01364684
	1234	7.5 x 7.5 cm	25 x 2 pcs	16 boxes	01364678
	1233	10 x 10 cm	25 x 2 pcs	16 boxes	01364661
6 ply, 30 g	4105	10 x 10 cm	25 x 2 pcs	28 boxes	13507290

ABE[®] drainal

Our slit hole compress with metal liner coating

Characteristics

- non-woven slit swab with hole
- soft and absorbent non-woven layer
- gentle on skin and wound friendly
- metallic coating prevents sticking to the wound, allowing a painless dressing change
- hole diameter: approx. 5.5 - 11 mm



Application

For external wound care in small calibre drains and infusion catheters. Apply metal coating to the skin.

sterile, pack of 1

REF	Size	Box contents	Shipping unit	PZN
4118	6 x 7 cm	50 x 1 pc	6 boxes	01333867

ABE[®] tracheal

Slit hole compress with metal coating

Characteristics

- non-woven material swabs, slit with hole
- soft and absorbent non-woven layer
- gentle on skin and wound friendly
- metallic coating prevents sticking to the wound, allowing a painless dressing change
- latex-free, DEHP and PVC
- sterile packed
- hole diameter: approx. 13 - 30 mm



Application

Slit hole compress for the treatment of catheters, stomata and drainage. The hole in the swab is also slit, thus ABE[®] tracheal makes a perfect fit to the drainage or stoma and the exudate can be quickly absorbed into the absorbent non-woven back. Apply metal coating to the skin.

sterile, pack of 1

REF	Size	Box contents	Shipping unit	PZN
4119	8 x 10 cm	50 x 1 pc	6 boxes	–

BeeSana[®] Absorbent dressing pad

Our absorbent dressing pad for high fluid absorption

Characteristics

- wound friendly absorbent dressing pad
- absorbent core made of 100% cellulose
- completely covered with non-adhesive non-woven material
- blue top to mark the back and as a hydrophobic barrier against wetting and as a soiling protection of laundry and clothing



Application

For wound treatment of severely secretive wounds in all medical areas. Absorbent dressing pads are also ideal for padding wounds.



non-sterile, in a polybag

REF	Size	Box contents	Shipping unit	PZN
1227	10 x 10 cm	25 pcs	16 boxes	01363816
1228	10 x 20 cm	25 pcs	8 boxes	07607202
1224	10 x 30 cm	25 pcs	15 boxes	12146619
1225	15 x 20 cm	25 pcs	12 boxes	12146625
1231	15 x 25 cm	25 pcs	7 boxes	00137727
1229	20 x 20 cm	25 pcs	12 boxes	08696337
1244	20 x 30 cm	25 pcs	5 boxes	12146631
1230	20 x 40 cm	25 pcs	6 boxes	08442499

sterile, pack of 1 sealed

1354	10 x 10 cm	60 x 1 pcs	10 boxes	08416059
1355	10 x 20 cm	30 x 1 pcs	12 boxes	06193514
1388	15 x 25 cm	30 x 1 pcs	8 boxes	12146602
1356	20 x 20 cm	30 x 1 pcs	6 boxes	08416042
1357	20 x 40 cm	25 x 1 pcs	4 boxes	07737904

BeeSana[®] Eye pad

Our eye pad with cushioning effect

Characteristics

- oval, soft and fully sealed dressing pad
- absorbent inner filling
- liquid is passed through the polyester non-woven material directly into the interior of the pad



Application

For absorption of eye fluid. For padding and eye bandages. Warning: In the case of injuries of the external eye as well as in applications after operations, a sterile eye pad / wound covering is to be used. Due to the oval shape, this wound dressing fits the eye socket particularly well.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
4195	5.6 x 7.2 cm	50 pcs	10 boxes	02214099
4155	5.6 x 7.2 cm	500 pcs	1 pack	12348872
4194	7 x 8.5 cm	50 pcs	10 boxes	02435297

sterile, pack of 1 sealed

4196	5 x 7 cm	25 x 1 pc	9 boxes	02214107
4156	7 x 8.5 cm	25 x 1 pc	6 boxes	02214136

ABE[®] Injection plaster

Our fabric injection plaster

Characteristics

- ready-to-use injection plaster
- with wound pad
- polyester mesh prevents the dressing pad sticking to the wound
- coated with skin friendly polyacrylate adhesive
- flexible and very robust due to woven fabric structure



Application

For minor injuries such as abrasions and superficial cuts as well as for care after blood extraction.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
9525	2 x 4 cm	500 pcs	20 Boxes	01664860

ABE[®] light

Our elastic wound dressing with wound bandage

Characteristics

- white, slightly elastic wound bandage
- centrally placed, absorbent wound pad
- Non-adhesive wound pad thanks to polyester mesh cover
- breathable, water vapour permeable carrier material, coated with skin friendly polyacrylate adhesive
- ideal fit thanks to high flexibility of the non-woven carrier



Application

For initial and follow-up care of small injuries, especially on particularly sensitive parts of the body. Ideal fit for joint wounds.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
9547	4 cm x 5 m	1 pc	14 boxes	06790233
9548	6 cm x 5 m	1 pc	10 boxes	06790256
9549	8 cm x 5 m	1 pc	8 boxes	06790262

ABE[®] Absorbent adhesive dressing

Our absorbent and skin friendly absorbent adhesive dressing

Characteristics

- self-adhesive absorbent dressing made of white non-woven material
- with skin-friendly polyacrylate adhesive
- with central, absorbent wound pad
- no sticking to the wound
- resistant to secreting
- permeable to water vapour and air



Application

For post-operative wound care or as a sterile wound dressing in the ward area for blood and secretion as well as for the protection of wounds against (cross) infections. Well suited for patients with sensitive skin.

sterile, pack of 1 sealed

REF	Size	Box contents	Shipping unit	PZN
4750	5 x 7 cm	100 x 1 pc	12 boxes	01452997
4752	6 x 10 cm	50 x 1 pc	6 boxes	01126022
4754	8 x 10 cm	50 x 1 pc	6 boxes	01446011
4756	10 x 10 cm	50 x 1 pc	6 boxes	03081973
4760	15 x 10 cm	50 x 1 pc	6 boxes	01445980
4758	12 x 10 cm	50 x 1 pc	6 boxes	01446040
4762	20 x 10 cm	50 x 1 pc	6 boxes	01445968
4764	25 x 10 cm	25 x 1 pc	6 boxes	01445951
4766	30 x 10 cm	25 x 1 pc	6 boxes	01445939
4768	35 x 10 cm	25 x 1 pc	6 boxes	01446057

ABE[®] derm

Transparent film wound bandage without wound dressing

Characteristics

- elastic, semi-permeable polyurethane film
- germ proof and waterproof, but permeable to oxygen and water vapour
- coated with hypo-allergenic polyacrylate adhesive
- wound environment visible at all times through the hypo-allergenic, transparent film
- can be removed with no pain and without leaving residue

Application

Well suited for post-operative wound care and to control the wound environment, wound care and catheter fixation. Also use directly on primary wounds with low exudation for wound monitoring.



sterile, pack of 1 sealed

REF	Size	Box contents	Shipping unit	PZN
9570	4 x 5 cm	100 x 1 pc	40 boxes	06790196
9571	6 x 7 cm	100 x 1 pc	40 boxes	06790204
9572	10 x 12 cm	50 x 1 pc	12 boxes	06790210
9573	15 x 20 cm	50 x 1 pc	15 boxes	06790227

ABE[®] derm plus

Transparent film wound bandage with wound dressing

Characteristics

- elastic, semi-permeable polyurethane film
- germ proof and waterproof, however oxygen and water vapour permeable
- coated with hypo-allergenic polyacrylate
- integrated, non-adhesive wound pad
- wound environment visible at all times by the hypo-allergenic, transparent film
- can be removed without pain and without leaving residue

Application

Well suited for post-operative wound care and to control the wound environment, wound care and catheter fixation. Also to be used directly on primary wounds with low to moderate exudation for wound monitoring.



sterile, pack of 1 sealed

REF	Size	Box contents	Shipping unit	PZN
9580	5 x 7 cm	50 x 1 pc	24 boxes	04589917
9581	6 x 10 cm	50 x 1 pc	64 boxes	04589923
9583	12 x 10 cm	50 x 1 pc	48 boxes	04592351
9584	15 x 10 cm	50 x 1 pc	12 boxes	04656720
9585	20 x 10 cm	50 x 1 pc	24 boxes	04656737

ABE[®] Cannula retention dressing I.V.

Our skin friendly cannula retention plaster

Characteristics

- self-adhesive retention plaster
- with central slit
- made of soft non-woven carrier
- integrated, non-adhesive wound pad
- hypo-allergenic polyacrylate adhesive with good adhesion



Application

For sterile wound care of the puncture site of intravenous cannulas and at the same time for fast and secure fixation. Well suited for patients with sensitive skin.

sterile, pack of 1 sealed

REF	Size	Box contents	Shipping unit	PZN
4457	8 x 6 cm	50 x 1 pc	20 boxes	12589854

ABE[®] Fixation non-woven fabric

Our skin friendly fixation non-woven fabric for a secure fixation

Characteristics

- self-adhesive fixation non-woven fabric made of white non-woven fabric
- stretchable in the transverse direction (at least 100%)
- particularly snug and soft
- good adhesion
- permeable to water vapour and air
- can be removed with no pain and without leaving residues
- hypo-allergenic polyacrylate adhesive



Application

For full area dressing coverage and fixation in all medical areas. For attachment of cannulas, drainage, probes, catheters, measuring instruments, etc. Due to the elasticity and individual cutting ideal for fixation on joints, round and conical body parts. Well suited for patients with sensitive skin.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
4470	5 cm x 10 m	1 pc	50 boxes	06168663
4471	10 cm x 10 m	1 pc	24 boxes	06168686
4472	15 cm x 10 m	1 pc	72 boxes	06168692
4473	20 cm x 10 m	1 pc	52 boxes	06168700

ABE[®] silk

Our synthetic silk roll plaster

Characteristics

- fixation plaster made of white cellulose acetate
- with hypo-allergenic polyacrylate adhesive
- serrated edge and therefore easy to tear longitudinally and transversely
- breathable
- permeable to water vapour
- water and dirt repellent
- easy to write on with a ballpoint pen



Application

For the fixation of wound dressings, bandages, drains, catheters and cannulas.

non-sterile, with protection ring

REF	Size	Box contents	Shipping unit	PZN
9531	2.5 cm x 5 m	24 pcs	10 boxes	11732887

non-sterile, without protective ring

9535	1.25 cm x 5 m	36 pcs	10 boxes	11734544
9536	2.5 cm x 5 m	24 pcs	10 boxes	11745536
9537	5 cm x 5 m	12 pcs	30 boxes	-
9540	1.25 cm x 9.14 m	36 pcs	10 boxes	03916834
9541	2.5 cm x 9.14 m	24 pcs	10 boxes	03916840

ABE® Tubular net bandage

Our highly elastic tubular net bandage

Characteristics

- highly elastic tubular bandage
- easy and quick handling
- can be cut at any point
- no tearing due to the binding
- seamless to avoid pressure points
- good anatomical adjustment for a wrinkle-free fit
- secure and permanent fixation, even with extreme movements



Application

For the fixation of all kinds of dressings and padding material in pressure-prone areas. Using scissors, a hole can be cut into the mesh tube so that uninjured limbs (e.g. the thumb in a hand bandage) are outside the bandage. Particularly suitable for use on joints, conical and round body parts.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
5164	Size 3 Forearm (2.2 cm x 25 m)	1 pc	10 boxes	03065804
5166	Size 4 Upper arm (2.5 cm x 25 m)	1 pc	8 boxes	03065810
5170	Size 5 Leg (4 cm x 25 m)	1 pc	5 boxes	03065827
5172	Size 6 Head (5.5 cm x 25 m)	1 pc	4 boxes	03065833
5176	Size 7 Torso (6.2 cm x 25 m)	1 pc	4 boxes	03065856
5180	Size 8 Torso XXL (8 cm x 25 m)	1 pc	3 boxes	03065862

ABE® Head bandage

Our ready-to-use net head bandage

Characteristics

- elastic, air-permeable head quick bandage
- quick and easy handling
- time saving compact solution
- to be put on as a cap
- secure fixation
- high patient comfort
- elasticity approx. 80 %
- material composition:
Polyester, elastane, cotton



Application

Fixation of pads / dressings. For the quick treatment of accident victims with head injuries, for wound care after operations and for use in the emergency department.

non-sterile

REF	Size (unstretched)	Box contents	Shipping unit	PZN
5193	M (approx. 18 x 10 cm)	25 pcs	5 boxes	12350828
5194	L (approx. 23 x 10 cm)	25 pcs	5 boxes	12554169

ABE[®] Nose bandage

Our flexible nose slings

Characteristics

- bi-elastic nose bandage
- with elastic ear loops or elastic band
- easy to put on
- centrally integrated, highly absorbent dressing pad
- dressing pad holds absorbed fluid inside and clings to the upper lip in an optimum way



Application

In the ENT area after nasal surgery with secretion.

Easy to put on thanks to the ear loops or through a knot at the back of the head.

non-sterile, with elastic band

REF	Box contents	Shipping unit	PZN
4039	20 pcs	10 boxes	01333749

non-sterile, with ear loops

REF	Box contents	Shipping unit	PZN
9814	20 pcs	10 boxes	11885562
9818	500 pcs	1 pack	13579533

Suspensory

Special bandage

Characteristics

- with detachable, individually adjustable pouch
- belt: 100 % cotton
- very skin friendly
- folding buckle contains iron. May also contain nickel and brass.
- fully washable

Application

To protect and relieve the testicles after inflammation or surgery.



non-sterile

REF	Size	Box contents	Shipping unit	PZN
1954SUS	Size 5	1 pc	1 pack	-
1955SUS	Size 6	1 pc	1 pack	-
1956SUS	Size 7	1 pc	1 pack	-
1957SUS	Size 8	1 pc	1 pack	-
1958SUS	Size 9	1 pc	1 pack	-
1959SUS	Size 10	1 pc	1 pack	-
1960SUS	Size 11	1 pc	1 pack	-
1961SUS	Size 12	1 pc	1 pack	-

ABE[®] Gauze bandage

Our traditional fixation bandage

Characteristics

- rigid fixation bandage according to DIN 61631
- smooth structure
- woven edges
- made of skin-compatible raw materials
- 65% viscose / 35% cotton



Application

For the fixation of bandages of all kinds. Wide range of application in all medical fields such as emergency bandage for first aid or as a carrier fabric for zinc paste bandages.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1400	4 cm x 4 m	20 pcs	16 boxes	01331029
1401	6 cm x 4 m	20 pcs	16 boxes	01331035
1402	8 cm x 4 m	20 pcs	12 boxes	01331041
1403	10 cm x 4 m	20 pcs	20 boxes	01331058

ABE[®] Last

Our elastic conforming bandage with a smooth structure

Characteristics

- highly elastic conforming bandage according to DIN 61634
- fine, smooth structure
- wrinkle-free and non-slip fit
- elasticity min. 135%
- stretched length: 4 m
- 54 % polyester / 46 % Viscose



Application

For the fixation of bandages of all kinds. For the fixation of wound dressings, for the attachment of support splints and as light support and compression bandage. Particularly suitable for use on joints, conical and round body parts. To immobilise the injured body parts and limbs.

non-sterile, pack of 20

REF	Size	Box contents	Shipping unit	PZN
1409	4 cm x 4 m	20 pcs	28 boxes	04853248
1410	6 cm x 4 m	20 pcs	24 boxes	01364282
1411	8 cm x 4 m	20 pcs	20 boxes	04853260
1412	10 cm x 4 m	20 pcs	15 boxes	04853277
1413	12 cm x 4 m	20 pcs	6 boxes	04853283

ABE[®] Last crepe

Our elastic conforming bandage with creped structure

Characteristic

- elastic conforming bandage according to DIN 61634
- creped structure
- effective adhesion of the individual layers
- made of skin friendly raw materials
- open, air permeable structure
- not susceptible to ointments and perspiration
- stretchability approx. 100 %
- stretched length: 4 m
- 70 % viscose / 30 % polyamide



Application

For the fixation of bandages of all kinds. For the fixation of wound dressings, for the attachment of support splints and as light support and compression bandage. Particularly suitable for use on joints, conical and round body parts. Well suited for patients with sensitive skin.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1749	4 cm x 4 m	40 pcs	12 boxes	07198318
1750	6 cm x 4 m	20 pcs	16 boxes	07198324
1751	8 cm x 4 m	20 pcs	12 boxes	07198330
1752	10 cm x 4 m	20 pcs	10 boxes	07198347
1753	12 cm x 4 m	20 pcs	8 boxes	07198353

double sterile, pack of 1 with 2 counting cards

4741	6 cm x 4 m	30 x 1 pcs	4 boxes	13655731
4742	8 cm x 4 m	30 x 1 pcs	4 boxes	13655642
4587	10 cm x 4 m	30 x 1 pcs	4 boxes	13655688
4588	12 cm x 4 m	30 x 1 pcs	4 boxes	13655725

ABE[®] fix color

Our elastic and cohesive smooth fixation bandage in colour

Characteristics

- cohesive and elastic conforming bandage
- smooth structure
- non-slip
- breathable and skin friendly
- elasticity according to DIN 61632: approx. 60%
- 61% polyamide / 39% cotton



Application

For the fixation of wound dressings. Ideal for dressings on many moving body parts. Especially suitable for immobilising the wrist and ankle.

non-sterile

REF	Size	Colour	Box contents	Shipping unit	PZN
19461	6 cm x 20 m	blue	1 pc	36 boxes	08520615
19463	6 cm x 20 m	green	1 pc	36 boxes	08520638
19465	6 cm x 20 m	red	1 pc	36 boxes	08520650

ABE[®] Last adhesive

Our elastic and cohesive bandage

Characteristics

- cohesive and elastic conforming bandage
- creped structure
- no additional fixation necessary
- good adhesion by means of micro dot impregnation
- large pore for good skin climate and high wearing comfort
- stretchability approx. 60 %
- stretched lengths: 4 m or 20 m
- 70 % polyester / 30 % polyamide



Application

For the fixation of bandages of all kinds. For the fixation of wound dressings, for the attachment of support splints and as light support and compression bandage. Particularly suitable for use on joints, conical and round body parts. To immobilise the injured body parts and limbs.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1289	4 cm x 4 m	1 pc	120 boxes	04019746
1391	6 cm x 4 m	1 pc	36 boxes	07198436
1392	8 cm x 4 m	1 pc	36 boxes	07198442
1393	10 cm x 4 m	1 pc	36 boxes	07198459
1394	12 cm x 4 m	1 pc	36 boxes	07198465
1288	4 cm x 20 m	1 pc	22 boxes	02435498
1383	6 cm x 20 m	1 pc	16 boxes	04853308
1384	8 cm x 20 m	1 pc	12 boxes	04853314
1385	10 cm x 20 m	1 pc	10 boxes	04853320
1386	12 cm x 20 m	1 pc	10 boxes	04853337

double sterile, pack of 1 with counting card

4793	12 cm x 4 m	15 x 1 pc	4 boxes	-
------	-------------	-----------	---------	---

BeeSana[®] SafeCut

Our cohesive finger bandage

Characteristics

- adhesive and elastic conforming bandage
- non-slip due to fine latex application
- actively breathable
- water resistant and tear resistant
- stretchability approx. 100 %
- bandage is easy to tear off



Application

For the fixation of bandages of all kinds. For individual and quick first aid for finger and joint injuries as well as torn fingernails. For the fixation of wound dressings as well as a light support and compression bandage. To immobilise the injured body parts and limbs.

Product advantage: easy to tear off (no scissors necessary).



non-sterile

REF	Size	Colour	Box contents	Shipping unit	PZN
1844	2.5 cm x 4.5 m	blue	1 pc	252 boxes	10544882

Fold the support pad (fixation possible with and without pad)



1. Take the end of the bandage and fold it over a few times until you have enough protective pad.
2. This support cushion is fixed in the intended place with the rest of the bandage, so that stability is achieved.
3. Apply the dressing and simply tear it off. Ready to go.

ABE[®] fix adhesion c

Our elastic and adhesive bandage, creped, individually cellophane wrapped

Characteristics

- adhesive and elastic conforming bandage
- non-slip
- breathable and skin-friendly
- elasticity: approx. 100%
- 52% polypropylene / 41% latex / 7% elastane
- Bandage is easy to tear off



Application

For the fixation of bandages of all kinds, especially joints. Excellent for all support and relief bandages.

Product advantage: easy to tear off (no scissors necessary).



non-sterile, individually wrapped in cellophane

REF	Size	Colour	Box contents	Shipping unit	PZN
MEP1396	5 cm x 4.5 m	blue	240 x 1 pc	1 pack	-
MEP1397	7.5 cm x 4.5 m	blue	192 x 1 pc	1 pack	-
MEP1398	10 cm x 4.5 m	blue	144 x 1 pc	1 pack	-

ABE[®] extra

Our short stretch bandage for extra strong support and compression bandages

Characteristics

- particularly strong, textile elastic compression bandage
- with short stretch / elasticity approx. 80%
- extraordinarily durable bandage with permanent compression effect
- open, air-permeable structure
- extended length: 5 m
- 100 % cotton
- with bandage clips



Application

For very strong compression with high working pressure and low resting pressure. Suitable for the opposite dressing technique according to Pütter. In support of the venous valve function, after varicosclerosis, in all acute and chronic venous stasis oedema, leg ulcers and thrombophlebitis. Well suited for patients with sensitive skin.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1421	8 cm x 5 m	1 pc	72 boxes	01364402
1422	10 cm x 5 m	1 pc	36 boxes	01364394
1735	6 cm x 5 m	10 pcs	10 boxes	08490759
1736	8 cm x 5 m	10 pcs	10 boxes	08490765
1737	10 cm x 5 m	10 pcs	5 boxes	08490771
1738	12 cm x 5 m	10 pcs	5 boxes	08490788

ABE[®] extra adhesion

Our adhesive, permanently elastic compression bandage with short stretch

Characteristics

- particularly strong, permanently elastic compression bandage
- with short stretch / elasticity approx. 80%
- adhesive
- non-slip
- selvedge
- actively breathable
- sterilisable with gas and steam at 121° C, X-Ray and gamma
- stretched length: 5 m
- 88 % cotton / 11 % polyamide / 1 % polyurethane



Application

For fixation, support and relief bandages. For chronic venous insufficiency, in primary and secondary varicosity, leg ulcers, oedema, sclerotherapy, thrombophlebitis, phlebothrombosis, postoperative care after venous surgery. Also prophylactic in the initial and sequential care of sports injuries of the musculoskeletal system.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1776	6 cm x 5 m	1 pc	64 boxes	01340927
1777	8 cm x 5 m	1 pc	64 boxes	01340933
1778	10 cm x 5 m	1 pc	32 boxes	01340956
1779	12 cm x 5 m	1 pc	32 boxes	01340979

ABE® Ideal 100

Our textile elastic ideal bandage made of 100% cotton with a short stretch

Characteristics

- textile elastic bandage according to DIN 61632
- elasticity approx. 90%
- raw white
- 100 % cotton
- stretched length approx. 5 m
- very skin friendly
- selvedge
- washable



Application

For relief bandages and after treatment of dislocations and fractures. Well suited for patients with sensitive skin.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1589	4 cm x 5 m	10 pcs	10 boxes	01364299
1423	6 cm x 5 m	10 pcs	10 boxes	01364477
1424	8 cm x 5 m	10 pcs	10 boxes	01364448
1425	10 cm x 5 m	10 pcs	5 boxes	04853484
1426	12 cm x 5 m	10 pcs	5 boxes	04853490
1576	15 cm x 5 m	10 pcs	5 boxes	01364460
1577	20 cm x 5 m	10 pcs	5 boxes	01364454
1429	30 cm x 5 m	12 pcs	1 pack	01331727

double sterile, pack of 1 with counting card

4771	6 cm x 5 m	25 x 1 pc	100 pcs	-
4772	8 cm x 5 m	25 x 1 pc	100 pcs	-
4773	10 cm x 5 m	20 x 1 pc	80 pcs	-
4774	12 cm x 5 m	15 x 1 pc	60 pcs	-
4775	15 cm x 5 m	20 x 1 pc	40 pcs	-
4776	20 cm x 5 m	18 x 1 pc	36 pcs	-

Surgical - Set ABE-Ideal® 100

4414	3 bandages*	15 x 3 pcs	2 boxes	13715918
------	-------------	------------	---------	----------

*Lengths: 10 cm, 12 cm and 15 cm (pack of 3 in a pouch)

ABE[®] lastic

Our permanently elastic short stretch bandage

Characteristics

- permanently elastic bandage for medium compression
- elasticity approx. 90%
- open, air-permeable structure
- stretched length: approx. 5 m
- latex-free and skin friendly
- 60% cotton / 38% polyester / 2% elastane



Application

For all dressings where light compression is desired, especially for compressive wound support and relief bandages. To immobilise the injured body parts and limbs. For phlebological and lymphological indications. For the fixation of wound dressings. Well suited for patients with sensitive skin.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1780	4 cm x 5 m	10 pcs	10 boxes	02423265
1781	6 cm x 5 m	10 pcs	10 boxes	02423259
1782	8 cm x 5 m	10 pcs	10 boxes	02423242
1783	10 cm x 5 m	10 pcs	10 boxes	02423236
1784	12 cm x 5 m	10 pcs	5 boxes	02423213
1785	15 cm x 5 m	10 pcs	5 boxes	02423207
1786	20 cm x 5 m	10 pcs	5 boxes	02423182

double sterile, pack of 1 with counting card

4779	6 cm x 5 m	30 x 1 pc	4 boxes	13655671
4780	8 cm x 5 m	25 x 1 pc	4 boxes	13655659
4781	10 cm x 5 m	20 x 1 pc	4 boxes	13655702
4782	12 cm x 5 m	15 x 1 pc	4 boxes	13655665

ABE[®] universal

Our permanently elastic all purpose bandage with medium strength stretch

Characteristics

- permanently elastic bandage for medium compression
- elasticity approx. 110%
- open, air-permeable structure
- stretched length: approx. 5 m
- latex-free and skin friendly
- 67% cotton / 30% polyester / 3% elastane



Application

For all dressings where moderate compression is desired, especially for compressive wound dressings.

Suitable for traumatological indications as well as for the treatment of venous disorders. For the fixation of wound dressings.

To immobilise the injured body parts and limbs and to support and relieve joints.

Well suited for patients with sensitive skin.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1414	6 cm x 5 m	10 pcs	16 boxes	00761561
1415	8 cm x 5 m	10 pcs	12 boxes	00761590
1416	10 cm x 5 m	10 pcs	10 boxes	00761526
1432	12 cm x 5 m	10 pcs	8 boxes	00761532
1433	15 cm x 5 m	10 pcs	6 boxes	00761549
1434	20 cm x 5 m	10 pcs	5 boxes	00761555

double sterile, pack of 1 with counting card

4722	6 cm x 5 m	25 x 1 pc	4 boxes	13655636
4723	8 cm x 5 m	25 x 1 pc	4 boxes	13655694
4572	10 cm x 5 m	20 x 1 pc	4 boxes	13655748
4560	12 cm x 5 m	15 x 1 pc	4 boxes	13655719

ABE[®] Ven

Our permanently elastic compression bandage with a long stretch

Characteristics

- permanently elastic
- elasticity approx. 175%
- insensitive to perspiration and ointments
- with latex-free end fixation
- skin colour
- 93% cotton / 7% polyurethane
- stretched length approx. 7 m
- very skin friendly
- washable
- with bandage clips



Application

For adjustable compression, absorption of haematomas, pregnancy varicosities, dislocations, distorsions and contusions. Well suited for patients with sensitive skin.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1563	8 cm x 7 m	1 pc	64 boxes	04854147
1564	10 cm x 7 m	1 pc	32 boxes	04854153
1565	12 cm x 7 m	1 pc	32 boxes	04854176

ABE[®] Padding bandage

Our synthetic wadding for padding



Characteristics

- dabbing elastic synthetic wadding from non-absorbent, crimped polyester fibres
- high elasticity is maintained even under continuous mechanical stress and under strong influence of moisture
- wrinkle-free fit due to non-slip self-adhesion
- easy to apply
- can be torn off by hand
- sterilisable

Application

As padding material under plaster, support and compression bandages.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1701	6 cm x 3 m	1 pc	108 boxes	06627958
1806	8 cm x 3 m	1 pc	72 boxes	04351788
1702	10 cm x 3 m	1 pc	72 boxes	06627964
1703	15 cm x 3 m	1 pc	36 boxes	06627970

sterile, pack of 1

1705	10 cm x 3 m	15 x 1 pc	4 boxes	13715901
------	-------------	-----------	---------	----------

ABE® compression bandage roll

Special bandage

Characteristics

Compression bandage made of polyurethane, in a knitted tube. Diameter 7 cm. Free of latex, PVC and DEHP, for single use.

Application

To be used as a pressure pad, e.g. after cardiac catheterisations.



non-sterile

REF	Size	Box contents	Shipping unit	PZN
4037	approx. 7 cm Ø x 10 cm	12 pcs	9 boxes	06620850

Testicular roll

Characteristics

The testicular roll is made of a foam with a knitted tube cover.

Application

It serves to elevate the testicles after surgery or oedema.

An addition, the legs are spread apart due to the special form of the testicle roll.



non-sterile

REF	Size	Box contents	Shipping unit	PZN
4035	Ø 8 cm x 10 cm	6 pcs	8 boxes	-

ABE[®] Tubular bandage

Our fine mesh knitted fabric tube for universal use

Characteristics

- tubular bandage 100% cotton
- longitudinal and transversally elastic
- easy handling
- can be cut at any location
- ladder proof
- seamless to avoid pressure points
- good anatomical fit for wrinkle-free fit
- reliable protection for wounds and sensitive skin
- breathable and absorbent at the same time



Application

As a pressure-free dressing as well as underlay and coating with zinc paste or plaster casts for skin protection. For sedentary bandages and extension bandages. For pulling over splints. For the fixation of bandages of all kinds.

non-sterile

REF	Size	Application	Box contents	Shipping unit	PZN
5181	Size 1 1.5 cm x 20 m	finger / toe	1 pc	24 boxes	01331354
5182	Size 2 2.5 cm x 20 m	several fingers / toes	1 pc	16 boxes	01331360
5183	Size 3 4 cm x 20 m	hand / forearm / child's arm	1 pc	12 boxes	01331377
5184	Size 4 6 cm x 20 m	upper arm / lower leg	1 pc	28 boxes	01331383
5185	Size 5 8 cm x 20 m	thigh	1 pc	22 boxes	01331420
5186	Size 6 10 cm x 20 m	head / child's torso	1 pc	18 boxes	01331437
5187	Size 7 12 cm x 20 m	torso	1 pc	16 boxes	01331443
5189	Size 9 21 cm x 20 m	torso XXL	1 pc	16 boxes	12350840

ABE[®] Gilchrist bandage

Our universal shoulder immobiliser

Characteristics

- universal shoulder immobiliser
- comfortable wearing properties
- easy and secure to fix



Application

For injuries of the shoulder and upper arm area for immobilisation or the fixation of the shoulder joint, e.g. after shoulder dislocation (shoulder luxation), for injuries of the clavicular joint (AC joint), in the upper arm, clavicular and lateral collarbone fractures. Occasionally even after shoulder operations (e.g. shoulder joint reflection) for immobilisation.

non-sterile

REF	Size	Box contents	Shipping unit	PZN
1914	universal	1 pc	12 boxes	01363437

ABE[®] Body bandage

Our post-operative support bandage

Characteristics

- support and relief bandage
- extra wide padding
- quick and secure putting on
- wide, adjustable Velcro fastener
- adaptation of compression to therapeutic requirements



Application

For post-operative support and treatment of rib fractures.

non-sterile

REF	Size	(Scope)	Box contents	Shipping unit	PZN
1890	Size 1	65 - 75 cm	1 pc	15 boxes	01363727
1891	Size 2	75 - 85 cm	1 pc	15 boxes	01363710
1892	Size 3	85 - 95 cm	1 pc	15 boxes	01363704
1893	Size 4	95 - 105 cm	1 pc	15 boxes	01363696
1894	Size 5	105 - 115 cm	1 pc	15 boxes	01363609
1895	Size 6	115 - 125 cm	1 pc	15 boxes	01363555
1904	Size 7	125 - 135 cm	1 pc	15 boxes	02435357

ABE[®] Heart catheter pressure bandage

Our pressure bandage for cardiac catheter

Characteristics

- quick and safe putting on
- adjustable velcro closure
- adaptation of compression to therapeutic requirements
- colour coded



Application

For compression of punctured femoral arteries / veins after invasive diagnosis and therapy. Pressure bandage for radiology.

non-sterile

REF	Size (universal)	Box contents	Shipping unit	PZN
2907	left (colour code: red)	1 pc	12 boxes	04351742
2908	right (colour code: blue)	1 pc	12 boxes	04351759

ABE[®] Limb holder

Our skin friendly cuff for adults

Characteristics

- cuff made of soft special material
- very skin friendly due to cotton pad
- velcro fastener for easy fixation
- no additional fastening tape necessary, since already firmly connected with cuff



Application

For short term immobilisation.



non-sterile

REF	Size	Box contents	Shipping unit	PZN
1879	for adults	10 pcs	12 boxes	12553046



PolyMem®

The multifunctional membrane dressing for all indications

- cleanse
- moisturizes
- absorbs
- fills

Cleansing

The continuous release of an innovative wound cleanser (F68) gently dissolves deposits and necrosis and supports the self-cleansing forces of the wound.

Moistening

PolyMem® keeps the wound bed moist. This creates an optimal wound environment for faster wound healing and reduces wound pain. PolyMem® does not adhere to the wound and can be easily removed with virtually no pain.

Absorbance

PolyMem® wicks away up to ten times its weight in exudate. It remains stable and can be removed without leaving any residue.

Fills

PolyMem® gently expands to fill and conform to the wound.

The portfolio for all requirements

Depending on the requirements, PolyMem® can be used both as a primary dressing and as a secondary dressing. It combines the requirements of a wide variety of wound dressings such as alginates, hydrocolloids etc. in a single dressing.

Clinical studies show that the PolyMem® membrane dressing provides faster wound healing compared to other dressings and has been successfully used in wounds where previous treatment was unsuccessful.

Benefits for you and your patients

Simplification of wound treatment

- PolyMem® can be used throughout the wound healing process.
- PolyMem® is suitable for all types of wounds.
- Due to its special properties, PolyMem® can be left on the wound for several days.

Creation of an optimal wound environment

- PolyMem® gently dissolves deposits and necrosis.
- PolyMem® absorbs excess exudate.
- PolyMem® keeps the wound bed constantly moist and reduces wound pain.



PolyMem® Product overview

EAN	Dimensions	Dimensions of the dressing pad	Packaging	PZN
-----	------------	--------------------------------	-----------	-----

PolyMem® Wound pad with self-adhesive PU film

PM203	5 x 5 cm	2,5 x 2,5 cm	5 x 20 St./VE	02919842
PM405	10 x 13 cm	5 x 8 cm	4 x 15 St./VE	00045184
PM606	15 x 15 cm	9 x 9 cm	4 x 15 St./VE	00045190
PM3412	32 x 10 cm	25 x 5 cm	1 x 12 St./VE	03480609

PolyMem® Wound pad with self-adhesive fixation fleece

PM7203	5 x 5 cm	2,5 x 2,5 cm	5 x 20 St./VE	02919919
PM7405	10 x 13 cm	5 x 8 cm	4 x 15 St./VE	00045238
silver PM1766	15 x 15 cm	9 x 9 cm	2 x 15 St./VE	02811122
PM7606	15 x 15 cm	9 x 9 cm	4 x 15 St./VE	00045244
PM7031	8 x 2,5 cm	2,5 x 2,5 cm	5 x 20 St./VE	02919954
PM7042	10 x 5 cm	4 x 5 cm	5 x 20 St./VE	00239445

PolyMem® Non-adhesive wound pad


silver PM1022	5 x 5 cm	5 x 5 cm	5 x 20 St./VE	05389959
PM5022	5 x 5 cm	5 x 5 cm	5 x 20 St./VE	07590660
PM5033	8 x 8 cm	8 x 8 cm	4 x 15 St./VE	00045356
silver PM1044	10 x 10 cm	10 x 10 cm	2 x 15 St./VE	02811300
PM5044	10 x 10 cm	10 x 10 cm	4 x 15 St./VE	00045362
PM5055	13 x 13 cm	13 x 13 cm	2 x 15 St./VE	00045379
silver PM1077	17 x 19 cm	17 x 19 cm	1 x 15 St./VE	02671655
PM5077	17 x 19 cm	17 x 19 cm	2 x 15 St./VE	00045385
silver PM1124	10 x 32 cm	10 x 32 cm	1 x 12 St./VE	02671661
PM5124	10 x 32 cm	10 x 32 cm	1 x 12 St./VE	00960177
PM5244	10 x 61 cm	10 x 61 cm	2 x 4 St./VE	00045416
PM5824	20 x 60 cm	20 x 60 cm	1 x 2 St./VE	02258754

PolyMem® Non-adhesive pad dressing

PM5035	8 x 8 cm	8 x 8 cm	4 x 10 St./VE	10837058
silver PM1045	10 x 10 cm	10 x 10 cm	2 x 8 St./VE	04923411
PM5045	11 x 11 cm	11 x 11 cm	2 x 10 St./VE	03367876
silver PM1088	20 x 20 cm	20 x 20 cm	2 x 5 St./VE	00151093
PM5088	20 x 20 cm	20 x 20 cm	2 x 5 St./VE	00151101

PolyMem® Wound Filler

silver PM1331	2,5 x 8 cm	2,5 x 8 cm	5 x 14 St./VE	05389965
silver PM1333	8 x 8 cm	8 x 8 cm	2 x 10 St./VE	02811292
silver PM1814	1 x 35 cm	1 x 35 cm	2 x 6 St./VE	00156943
PM5731	2,5 x 8 cm	2,5 x 8 cm	5 x 14 St./VE	09266484
PM5733	8 x 8 cm	8 x 8 cm	4 x 10 St./VE	00045439
PM5712	8 x 30 cm	8 x 30 cm	1 x 12 St./VE	03508839
silver PM1788	20 x 20 cm	20 x 20 cm	2 x 5 St./VE	08468228
PM5788	20 x 20 cm	20 x 20 cm	2 x 5 St./VE	07590648



Shapes of PolyMem®

Product overview

EAN	Designation	Dimensions	Dimensions of the dressing pad	Packaging	PZN
PM3709	Sacral	18 x 20 cm	11 x 12 cm	2 x 10 pcs./Unit	02671891
PM1709	Sacral Silver	18 x 20 cm	11 x 12 cm	2 x 10 pcs./Unit	02671885
PM8086	Oval	17 x 21 cm	10 x 14,5 cm	2 x 10 pcs./Unit	02671744
PM1886	Oval Silver	17 x 21 cm	10 x 14,5 cm	2 x 10 pcs./Unit	02671879
PM8053	Oval	9 x 13 cm	5 x 7,5 cm	2 x 15 pcs./Unit	02671721
PM1853	Oval Silver	9 x 13 cm	5 x 7,5 cm	2 x 15 pcs./Unit	02671738
PM8023	Oval	5 x 7,5 cm	2,5 x 5 cm	5 x 20 pcs./Unit	02671690
PM1823	Oval Silver	5 x 7,5 cm	2,5 x 5 cm	5 x 20 pcs./Unit	02671709
PM5335	Tube	9 x 9 cm	9 x 9 cm	2 x 15 pcs./Unit	00045422
PM1401	Finger / Toe silver	Size 1	XS	5 x 6 pcs./Unit	09221990
PM1402	Finger / Toe silver	Size 2	S	5 x 6 pcs./Unit	09222009
PM1403	Finger / Toe silver	Size 3	M	5 x 6 pcs./Unit	09222015
PM1404	Finger / Toe silver	Size 4	L	5 x 6 pcs./Unit	10143930
PM1405	Finger / Toe silver	Size 5	XL	5 x 6 pcs./Unit	10143947
PM4401	Finger / Toe	Size 1	XS	5 x 6 pcs./Unit	09063752
PM4402	Finger / Toe	Size 2	S	5 x 6 pcs./Unit	09063775
PM4403	Finger / Toe	Size 3	M	5 x 6 pcs./Unit	09063781
PM4404	Finger / Toe	Size 4	L	5 x 6 pcs./Unit	10143887
PM4405	Finger / Toe	Size 5	XL	5 x 6 pcs./Unit	10143918





BeeSana® Sets - high efficiency for maximum individualisation

Our set solutions enable efficient and safe use on the wards and in the operating theatre. The various combinations of individual components in high quality packaging allow our customers the time saving and safe use of our products. The objective of our range is the economic application of the different individual components under the constant premise of high quality patient care.

Of course, there are set solutions both as ready made standard sets as well as custom made individual sets.

Individual set solutions

The contents of each set are compiled individually according to the requirements of our customers.

Together with our customers, our trained and competent sales staff develop individual solutions for the highest demands and special interventions. We also offer our customers solutions for minimum requirements. Our entire product portfolio of dressings, covers, protective clothing and other materials are the basis for this. Of course, all our sets for the surgical area are equipped with double, self-adhesive counting cards for safe documentation in the patient record. By using our sets, material usage, preparation time, waste volume and logistics costs are reduced. In addition, you ensure that all components are available for use in the optimum number. Another advantage is the exact assignment to the patient and intervention of the materials used. The use of set solutions leads to a significant optimisation of the process costs and a high level of customer profitability.

Standardised set solutions

Based on our many years of experience in the development of customised set solutions for our customers, we have been able to develop various standard sets for various procedures both in the operating theatre and on the ward. Sets especially adapted to the needs of the various operations facilitate the preparation and implementation of standard interventions and processes for our customers. The use of set solutions leads to a significant optimisation of the process costs and a high level of customer profitability.

BeeSana® Ward catheter set

Characteristics

The ward catheter set is used for hygienic intermittent catheterisation of bladders and is sterile for immediate use on the ward.



Application

For catheterisation of the bladder.

REF	Contents	Box contents	Shipping unit	PZN
4265	1 x rectangular dish 650 ml 1 x swab dish 125 ml 2 x gauze swabs 7.5 x 7.5 cm, 8 ply 6 x gauze ball 20 x 20 cm, plum size 1 x plastic tweezers 1 x slit cloth 60 x 60 cm 1 x cover cloth 60 x 60 cm	21 sets	2 boxes	04136073
4266	1 x rectangular dish 650 ml 1 x swab dish 125 ml 2 x gauze swabs 7.5 x 7.5 cm, 8 ply 6 x gauze ball 20 x 20 cm, plum size 2 x examination gloves size M, latex 1 x plastic tweezers 1 x slit cloth 60 x 60 cm 1 x cover cloth 60 x 60 cm	21 sets	2 boxes	04136050
4267	1 x rectangular dish 650 ml 1 x swab dish 125 ml 2 x gauze swabs 7.5 x 7.5 cm, 8 ply 6 x gauze ball 20 x 20 cm, plum size 2 x examination gloves size L, latex 1 x plastic tweezers 1 x slit cloth 60 x 60 cm 1 x cover cloth 60 x 60 cm	21 sets	2 boxes	04136067
4535	1 x kidney dish 1 x disposable tweezers, blue 1 x slit cloth 60 x 60 cm, blue 3 x latex gloves size M, powder-free 2 x gauze swabs 7.5 x 7.5 cm, 8 ply 1 x swab dish 4 x gauze ball 20 x 20 cm, plum size, round enveloped in: 1 x cover cloth 55 x 60 cm	14 sets	2 boxes	-

BeeSana® PEG bandage kit

Characteristics

The PEG bandage kit is used to apply a bandage hygienically in a percutaneous endoscopic gastrostomy and is sterile for immediate use on the ward.



Application

For external wound care.

REF	Contents	Box contents	Shipping unit	PZN
4546	2 x gauze swabs 7.5 x 7.5 cm, 12 ply 2 x slit non-woven swabs 7.5 x 7.5 cm, 7.5 cm, 4 ply 4 x gauze ball 20 x 20 cm, plum size 1 x ABE® Fixation non-woven fabric 12 x 15 cm	30 sets	8 boxes	-
4653	2 x non-woven slit swabs 7.5 x 7.5 cm, 4 ply, 30g 2 x non-woven swabs 7.5 x 7.5 cm, 6 ply, 30g 1 x ABE® Fixation non-woven fabric 12 x 15 cm	10 sets	30 boxes	14399023
4577	2 x slit non-woven swabs 7.5 x 7.5 cm, 7.5 cm, 4 ply 2 x non-woven swabs 7.5 x 7.5 cm, 7.5 cm, 6 ply 1 x ABE® Fixation non-woven fabric 12 x 15 cm	40 sets	8 boxes	01333962
4589	2 x slit non-woven swabs 7.5 x 7.5 cm, 7.5 cm, 4 ply 2 x non-woven swabs 7.5 x 7.5 cm, 7.5 cm, 6 ply 1 x ABE® Fixation non-woven fabric 12 x 15 cm	30 sets	8 boxes	04152190
4836	2 x gauze swabs 7.5 x 7.5 cm, 8 ply 1 x drainage pad (Al/non-woven material) 6 x 7 cm, with slit hole 1 x ABE® Fixation non-woven fabric 12 x 15 cm 1 x fixation holder size 2, small, with velcro fastener	40 sets	8 boxes	15995170

BeeSana[®] Bandage change kit

Characteristics

The bandage change kit is used for hygienic dressing changing and is sterile for immediate use on the ward.



Application

For external wound care.

REF	Contents	Box contents	Shipping unit	PZN
4820	2 x gauze swabs 10 x 10 cm, 8 ply 4 x gauze ball 20 x 20 cm, plum size	25 sets	4 boxes	00549660
5413	2 x gauze swabs 5 x 5 cm, 16 ply 4 x gauze swabs 10 x 10 cm, 12 ply	20 sets	4 boxes	-

BeeSana[®] Suture removal kit

Characteristics

The suture removal kit is applied to remove threads hygienically and is sterile for immediate use on the ward.



Application

For external wound care.

REF	Contents	Box contents	Box contents	PZN
5450	4 x gauze swabs 10 x 10 cm, 12 ply 4 x gauze ball 25 x 26 cm, egg size 1 x suture removal scissors 11.5 cm 2 x disposable tweezers 130 mm, blue 1 x cover 45 x 75 cm, 2 ply, green	15 sets	4 boxes	-

BeeSana® Ward dialysis kit

Characteristics

The ward dialysis kit is used for external hygienic wound care after dialysis and is sterile for immediate use on the ward.



Application

For external wound care.

sterile

REF	Contents	Box contents	Shipping unit	PZN
4545	3 x non-woven swabs 5 x 5 cm, 30 g, 4ply 2 x gauze ball 20 x 20 cm, plum size peel pouch with 2 cells	40 sets	8 boxes	-

BeeSana® Ward protection set

Characteristics

The ward protection set is intended to prevent patient-to-patient contamination or infection of medical and nursing staff.



Application

For protection against MRSA.

non-sterile

REF	Contents	Box contents	Shipping unit	PZN
4074	1 x PP gown 1 x nurse hood 1 x OP mask 2 x CPE shoe covers 2 x examination gloves latex size M	80 sets	1 pack	-



Our incontinence products are oriented towards quality of life of those affected

Today there are a whole range of conservative and surgical treatments for urinary and faecal incontinence that can cure or at least alleviate the condition. Modern incontinence products provide a pleasant solution for people suffering from incontinence who can no longer be treated. From an application point of view, it is by no means sufficient to provide nurses and patients with just any incontinence product. This is why we go to great lengths when developing and producing Meditrade incontinence products to come up with detailed solutions that also provide discretion for the user. Our broad range enables cost-effective use with a wide range of sizes, absorbency levels and qualities, always ensuring specific and acceptable patient care.

Superslip®

Anatomically shaped incontinence slip with flexible cuffs

Characteristics

Modern incontinence briefs with excellent comfort and safety thanks to elastic back fastener and double adhesive velcro closure. This allows repeated opening and closing without damaging the outer layer. The odour binding and pH skin neutral double cushion technology with super-absorber achieves high fluid absorption with constant absorption speed. Hydrophobic side cuffs as an outlet barrier as well as a wetness indicator to indicate the optimum changeover time all offer maximum protection. The soft and skin friendly inner non-woven material prevents skin irritation.

Application

For severe urinary and faecal incontinence and persons with limited mobility who have a high need of care.



Superslip® Classic

REF	Size	Colour code	Hips	Absorbency (according to ISO 11948-1)	Box contents	Shipping unit	PZN	Aid no.
2431	Medium	white / blue	80 - 115 cm	approx. 1900 ml	20 pcs	4 boxes	07757522	15.25.03.1118
2432	Large	white / green	115 - 150 cm	approx. 2100 ml	20 pcs	4 boxes	07757551	15.25.03.2112

Superslip® Comfort

2406	Small	white	50 -80 cm	approx. 1800 ml	15 pcs	6 boxes	07757485	15.25.03.0066
2407	Medium	white / blue	80 - 115 cm	approx. 2450 ml	15 pcs	8 boxes	07757580	15.25.03.1119
2408	Large	white / green	115 - 150 cm	approx. 2600 ml	15 pcs	6 boxes	07757605	15.25.03.2113
2414	X-Large	white / violet	> 150 cm	approx. 2700 ml	15 pcs	4 boxes	07757634	15.25.03.2312

Superslip® Extra

2435	Medium	white / blue	80 - 115 cm	approx. 2850 ml	15 pcs	4 boxes	10063817	15.25.03.1277
2436	Large	white / green	115 - 150 cm	approx. 3000 ml	15 pcs	4 boxes	09747019	15.25.03.2314

Superslip® Super

2403	Medium	white / blue	80 - 115 cm	approx. 3000 ml	15 pcs	4 boxes	09754143	15.25.03.1278
2404	Large	white / green	115 - 150 cm	approx. 3200 ml	15 pcs	4 boxes	00351917	15.25.03.2313

Superslip® Air Comfort

Breathable, anatomically shaped incontinence slip with flexible cuffs

Characteristics

Modern incontinence briefs with excellent comfort and safety thanks to elastic back fastener and double adhesive velcro closure. This allows repeated opening and closing without damaging the outer layer. The odour binding and pH skin neutral double cushion technology with super-absorber achieves a high fluid absorbency while maintaining the absorption speed. Hydrophobic side cuffs as an outlet barrier as well as a wetness indicator to indicate the optimum changeover time all offer maximum protection. The soft and skin friendly inner non-woven material prevents skin irritation.

Application

For severe urinary and faecal incontinence and persons in need of a high degree of care with limited mobility.



REF	Size	Colour code	Hips	Absorbency (according to ISO 11948-1)	Box contents	Shipping unit	PZN	Aid no.
2437	Medium	white / blue	80 - 115 cm	approx. 1750 ml	15 pcs	8 boxes	04917043	15.25.03.1276
2438	Large	white / green	100 - 150 cm	approx. 2000 ml	15 pcs	6 boxes	04918189	15.25.03.2311

Medi-Pants

Discrete disposable briefs with the highest wearing comfort

Characteristics

Incontinence briefs for protection against bladder weakness, which can be put on and taken off like normal underwear. Enables unrestricted freedom of movement and customary independence. The odour binding and fast absorbing core reliably and fluid tightly encloses liquid. Breathable, skin friendly, for women and men, for single use.

Application

To protect against bladder weakness. For mobile or people with dementia without additional restrictions in everyday life. Especially suitable for use in toilet training. Wear just like normal underwear.



REF	Size	Hips	Box contents	Shipping unit	PZN
2411	Medium	approx. 70 - 110 cm	14 pcs	6 boxes	11141726
2412	Large	approx. 100 - 140 cm	14 pcs	6 boxes	11141732
2413	X-Large	approx. 120 - 165 cm	14 pcs	6 boxes	11141749

Liberty Light

Anatomically shaped incontinence pad

Characteristics

Insert with odour binding and pH-neutral skin. Anatomical fit with soft, hydrophobic cuffs in the crotch area. With the adhesive tape, the hypo-allergenic insert is easy to fix in the underwear.

Application

Suitable for mild to moderate bladder weakness, as a protection after vaginal surgery or as a comfortable leak-proof inlay after childbirth.



REF	Size	Colour code	Dimensions (approx.)	Absorbency (according to ISO 11948-1)	Box contents	Shipping unit	PZN	Aid no.
2052	Normal	pink	28 x 10.5 cm	min. 350 ml	22 pcs	10 boxes	07305197	15.25.01.5175
2053	Extra	blue	31.5 x 12 cm	min. 400 ml	18 pcs	10 boxes	07305180	15.25.01.5177
2054	Super	green	43 x 22 cm	min. 850 ml	30 pcs	5 boxes	07305174	15.25.01.5176

Superform[®]

Anatomically shaped incontinence pad with double cushions

Characteristics

Anatomically shaped incontinence pad with excellent comfort and security for the person affected. The odour binding and pH skin neutral double cushion technology with super-absorber achieves high fluid absorption with constant absorption speed. Hydrophobic lateral cuffs as an outlet barrier as well as a moisture indicator to indicate the optimum change time provide maximum protection. The soft and skin friendly inner non-woven material prevents skin irritation.

Application

Suitable for every degree of incontinence in bedridden and mobile persons. **Combination note:** For secure fixation we recommend using it with our washable PANTY or LIGHT fixation pants.



REF	Size	Colour code	Dimensions	Absorbency (according to ISO 11948-1)	Box contents	Shipping unit	PZN	Aid no.
2047	Classic	blue	63 x 30 cm	approx. 1200 ml	30 pcs	5 boxes	07520263	15.25.01.0076
2048	Plus	blue	63 x 30 cm	approx. 1350 ml	30 pcs	5 boxes	07414302	15.25.01.1096
2049	Extra	green	70 x 37 cm	approx. 1450 ml	30 pcs	4 boxes	07414288	15.25.01.1105
2050	Super	green	70 x 37 cm	approx. 1850 ml	20 pcs	4 boxes	07414271	15.25.01.1065
2051	Maxi	white	70 x 37 cm	approx. 2450 ml	20 pcs	4 boxes	07414070	15.25.01.2101

Superform[®] Air

Breathable, anatomically shaped incontinence pad with double cushion

Characteristics

Anatomically shaped, breathable incontinence pad with excellent comfort and high security. The odour binding and pH skin neutral double cushioning technology with super-absorber achieves a high fluid absorption with constant absorption speed. Hydrophobic side cuffs as an outlet barrier as well as a wetness indicator to indicate the optimum changeover time all offer maximum protection. The soft and skin friendly inner non-woven material prevents skin irritation.

Application

Suitable for every degree of incontinence in bedridden and mobile persons. **Combination note:** For secure fixation we recommend using it with our washable PANTY or LIGHT fixation pants.



REF	Size	Colour code	Absorbency (according to ISO 11948-1)	Box contents	Shipping unit	PZN	Aid no.
20481	Plus	white / blue	approx. 1250 ml	30 pcs	5 boxes	08801840	15.25.01.0098
20491	Extra	white / yellow	approx. 1400 ml	30 pcs	4 boxes	08801857	15.25.01.1125
20501	Super	white / green	approx. 1800 ml	20 pcs	4 boxes	08801863	15.25.01.2141
20511	Maxi	white / violet	approx. 2200 ml	20 pcs	4 boxes	08834851	15.25.01.2142

BeeSana[®] Absorbent pad

Highly absorbent pad for efficient incontinence and nursing

characteristics

- absorbent pad made of cellulose flock
- with soft non-woven cover and impermeable, non-slip outer foil
- with super-absorber
- absorbs moisture quickly and completely
- protects laundry and clothes reliably from soiling
- easy to put on in combination with Light or Panty fixation pants
- latex-free



Application

For light incontinence / light stool for catheter wearers, bedridden and mobile persons. Ideal as a universal absorbent pad in nursing. Well suited for patients with sensitive skin.

	REF	Size	Packaging	Box contents	PZN
with adhesive tape	2022	12 x 38 cm	50 pcs	8 boxes	02214142
without adhesive tape	2024	18 x 40 cm	60 pcs	3 boxes	02214165

BeeSana® Absorbent pad

Our high capacity and extra large pad for incontinence care

Characteristics

- absorbent pad made of cellulose flock
- with soft non-woven cover
- particularly absorbent
- easy to put on in combination with Light or Panty fixation pants
- secure and complete coverage of the genital area
- latex-free



Application

For light incontinence for bedridden persons. Well suited for patients with sensitive skin.

REF	Size	Box contents	Shipping unit	PZN
2119	16 x 55 cm	20 pcs	5 boxes	02214219

BeeSana® Rectangular pad

Our absorbent pad for light incontinence

Characteristics

- absorbent pads made of cellulose flock
- with soft non-woven cover
- particularly absorbent
- easy to put on in combination with Light or Panty fixation pants
- latex-free



Application

For light incontinence for bedridden persons. As a maternity pad for the postpartum period. Suitable as protection after vaginal interventions.

Well suited for patients with sensitive skin.

REF	Size	Box contents	Shipping unit	PZN
2101	35 x 11 cm	25 pcs	8 boxes	12147808

BeeSana[®] Male incontinence pad

Discrete and anatomically shaped pad for men

Characteristics

Male incontinence pad with Super-absorber absorbent core. Free of latex, PVC and DEHP. Ergonomic cut for the male anatomy for secure protection with maximum comfort. Inconspicuous to wear. With adhesive strips easy to fix in the underwear.



Application

For men with mild urinary incontinence, e.g. as a result of operations in the genital area, as a secure protection for predominantly active persons.

REF	Size	Box contents	Shipping unit	PZN
2020	8.5 x 40 cm	21 pcs	7 boxes	05880164

BeeSana[®] Women's incontinence pad

Our hygienic women's incontinence pad

Characteristics

- absorption pad made of soft, highly absorbent cellulose flock
- skin-friendly, silky soft non-woven cover
- latex-free



Application

As a protection after vaginal interventions or as an inlay / maternity pad after delivery.

REF	Size	Box contents	Shipping unit	PZN
2118	approx. 7 x 23 cm	20 pcs	18 boxes	01091275

Light Fixation pants

For the fixation of incontinence pads

Characteristics

High quality, latex-free fixation pants with soft, linen-like textile structure. Permanent size marking by colour stripes. Washable at 60 degrees and suitable for tumble drying. Even after frequent washing, remains form and dimensionally stable.



Application

For safe and easy-to-wear fixation of all incontinence pads and inlays for mobile persons.

REF	Size	Colour	Hips	Box contents	Shipping unit	PZN
2221	Medium	green	approx. 80 - 110 cm	50 pcs	6 boxes	11513151
2222	Large	white	approx. 90 - 130 cm	50 pcs	6 boxes	11513211
2223	X-Large	blue	approx. 110 - 170 cm	50 pcs	6 boxes	11513228
2224	XX-Large	orange	approx. 120 - 180 cm	50 pcs	6 boxes	11513234

Panty Fixation pants

With leg attachment for the fixation of incontinence pads

Characteristics

High quality, latex-free fixation pants with high wearing comfort due to soft, cloth-like textile structure. Permanent size marking by colour stripes. Due to the leg and the fine stitch pattern, a feeling similar to a fine tights is achieved. Washable at 60 degrees and suitable for tumble drying. Even after frequent washing, remains form and dimensionally stable.



Application

For safe and easy-to-wear fixation of all incontinence pads and inlays for mobile persons.

REF	Size	Colour	Hips	Box contents	Shipping unit	PZN
0711500301	Medium	blue	approx. 70 - 100 cm	100 pcs	6 boxes	09081879
0711500501	Large	brown	approx. 85 - 120 cm	100 pcs	6 boxes	09081885
0711500701	X-Large	green	approx. 105 - 130 cm	100 pcs	6 boxes	09081891
0711500901	XX-Large	green-blue	approx. 115 - 140 cm	100 pcs	6 boxes	09081916

Layered underpad

Proven underpad from environmentally friendly layers of cellulose for needs-based use

Characteristics

Underpad with absorbent pads made of environmentally friendly, 100% unbleached recycled cellulose for rapid and full moisture absorption. The body friendly dry non-woven with non-slip outer film ensures a needs-based and efficient application.

Application

Underpads for an even greater feeling of security when working with patients wearing incontinence products on the body. They help in the home care sector, to increase the comfort of care, to facilitate daily work and not least they serve as bed protection to reduce the amount of laundry. We supply them as layer and flock underpads in different qualities.



REF	Size	Number of layers	Shipping unit	PZN
104606M	40 x 60 cm	6 ply	200 pcs	02426826
104608M	40 x 60 cm	8 ply	200 pcs	00559612
104610M	40 x 60 cm	10 ply	200 pcs	00559629
104612M	40 x 60 cm	12 ply	100 pcs	06178087
104620M	40 x 60 cm	20 ply	100 pcs	00559693
106606M	60 x 60 cm	6 ply	150 pcs	00560213
106610M	60 x 60 cm	10 ply	100 pcs	00560236
106906M	60 x 90 cm	6 ply	100 pcs	00560259
106908M	60 x 90 cm	8 ply	100 pcs	00560331
106910M	60 x 90 cm	10 ply	100 pcs	00560609
106912M50	60 x 90 cm	12 ply	50 pcs	00561046
106920M	60 x 90 cm	20 ply	50 pcs	00561052

BeeSana[®] Comflock Super

Underpad made of soft, absorbent cellulose flock

Characteristics

Underpad with absorbent body made of cellulose flock for rapid and full-surface moisture absorption. The body friendly dry non-woven with non-slip outer film ensures a needs-based and efficient application.

Application

Underpad for an even greater feeling of security when working with patients wearing incontinence products on the body. They help in the home care sector, to increase the comfort of care, to facilitate daily work and not least they serve as bed protection to reduce the amount of laundry. We supply them as layer and flock underpads in different qualities.



REF	Size	Liquid absorbency	Box contents	Shipping unit	PZN
2064	40 x 60 cm	750 ml ± 10 %	25 pcs	4 boxes	04136185
2065	60 x 60 cm	900 ml ± 10 %	30 pcs	6 boxes	04136104
2066	60 x 90 cm	1590 ml ± 10 %	25 pcs	6 boxes	04136156

Nursing pads

Our extra soft inlay for secure laundry protection during breastfeeding

Characteristics

- anatomically shaped
- particularly soft and skin friendly
- extra absorbent for a long lasting dry feeling
- cellulose core absorbs escaping milk quickly and reliably
- moisture impermeable and breathable outer non-woven material
- inconspicuous, secure fit

Application

For safe fluid intake and protection during breastfeeding after birth.



REF	Size	Box contents	Shipping unit	PZN
2068	10.5 cm	30 pcs	12 boxes	02097937

New Quik

Baby nappy in premium quality

Characteristics

Baby nappies with double side leakage protection and resealable Velcro straps. An elastic back cuff provides a secure, comfortable fit and a high level of protection against moisture leakage. The super absorbent inner cushion keeps your baby's skin healthy and absorbs moisture that could otherwise cause a rash. Very skin friendly, because the ultra-soft and breathable textile outer skin of the nappy protects the delicate skin of your baby and helps to prevent skin irritation. The highly elastic side wings allow an optimum fit and very good freedom of movement. For single use.



Application

Baby nappies in all sizes for newborns to pre-school children.

REF	Size	Baby weight (in kg)	Box contents	Shipping unit	PZN
2077	Mini	from 3 - 6	22 pcs	6 boxes	-
2070	Midi	from 4 - 9	20 pcs	6 boxes	-
2071	Maxi	from 7 - 18	18 pcs	6 boxes	-
2078	Junior	from 12 - 25	16 pcs	6 boxes	-
2079	X-Large	from 15 - 30	14 pcs	6 boxes	-



Disinfection and hygiene products with a broad spectrum of action and modern ingredients

The core elements of our broad range of products reduce the risk of infection, standardize prevention methods and provide lasting skin care.

Our disinfection and cleaning products are subject to the strictest quality standards and are constantly and extensively tested. We offer thorough optimisation and cost savings in the form of hygiene designs in line with our products at our production site so that you can keep up-to-date with current hygiene standards.

Note: The recommendations for our products are based on scientific tests and are given to the best of our knowledge. Further recommendations, e.g. in terms of material compatibility, are only possible separately. Our recommendations do not exempt our customers from testing our products independently to ensure they are suitable for the intended purposes and procedures. In that regard, we can not assume any liability.

Alcoman[®] +

The proven hand disinfectant as a biocide

Characteristics

Ready-to-use, highly effective alcohol based antiseptic for hygienic and surgical hand disinfection. Very good skin compatibility. Especially suitable for frequent applications. Dermatologically tested. Acts quickly and thoroughly.

Application

For hygienic hand disinfection: Rub portions of 3 ml on hands to keep them moist for 30 seconds. For surgical hand disinfection: thoroughly moisten dry hands and forearms with Alcoman[®] +. Repeat if necessary to keep hands and forearms moist for at least 5 minutes. It should be noted that all finger gaps, fingernails and skin folds are also completely wetted over the entire exposure time.

Spectrum of action

Comprehensively effective against bacteria and yeasts. Bactericidal, fungicidal, tuberculocidal, virucidal and therefore completely effective against all enveloped viruses (BVDV / HCV / vaccinia incl. HBV / HIV). Effective against rota virus and noro virus and MRSA. Effective against Adeno in 2 minutes.

Shelf life

Alcoman[®] + can be used for up to 12 months after opening the original container for the first time. The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 and 25°C.

Regulatory requirements

Alcoman[®] + is not a medicinal product and is covered by the Biocidal Products Directive. Biocide safe use. Always read the label and product information before use. VAH listed. Biocide of Product type 1 ("Biocidal Products for Human Hygiene")

Contents

Composition (included in 100 g): 70 ml of 2-propanol



REF	Box contents	Shipping unit	PZN
01007B	150 ml	24 boxes	13814980
00978B	500 ml	18 boxes	13814997
00983B	1 L	10 boxes	13815040
00984B	5 L	2 canisters	13815057

Alcoman® Gel

Our hand disinfection in handy gel version

Product type

Alcoman® Gel is a ready-to-use gel based alcoholic rub-in preparation for hygienic hand disinfection. The application is carried out without the addition of water.

Ideal for travelling.

Preparation properties

- acts quickly and thoroughly in 30 sec.
- thoroughly effective against bacteria and yeast fungi
- limited virucidal against all enveloped viruses: BVDV / vaccinia (incl. HCV / HBV / HIV) and noro viruses (MNV)
- evaluated according to the standard methods of DGHM and DIN EN 1500 as well as the RKI recommendations on effectiveness against viruses
- without colourings
- very good skin compatibility even with frequent use
- VAH listed

Applications

For all areas with increased demands on hand hygiene: Clinics, medical practices, ambulances, care facilities, industrial, commercial and laboratory areas, leisure and bathing facilities. Ideal for mobile activities.

Type and duration of application

Ready to use. For hygienic hand disinfection according to DIN EN 1500 min. 3 ml for 30 seconds thoroughly rub in undiluted and without adding water.

Avoid weak spots.

Regulatory requirements

Biocide of Product type 1 ("Biocidal products for human hygiene").

Use disinfectant with care. Always read the label and product information before use.

Shelf life

Alcoman® Gel can be used up to 12 months after the first opening of the original container. The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) between 5 - 25 ° C.



Composition:

The contents of 100 ml are:

70 ml of 2-propanol;

Other ingredients:

Gel former, moisturiser, perfume.

REF	Box contents	Shipping unit	PZN
00965	150 ml	24 bottles	11169713
00966	500 ml	18 bottles	11169736

Accessories (see page 160)

Dispenser pump for 500 ml bottle

Application: Conc. / Exposure time

Efficacies	15 sec.	30 sec.	1 min. / 60 sec.	1,5 min. / 90 sec.	2 min. / 120 sec.
Hygienic hand disinfection according to EN 1500 / VAH		+			
MRSA	+				
Limited virucidal according to EN 14476					
BVDV / vaccinia (including HBV, HCV, HIV)			+		
Noro (MNV)			+		
Vaccinia Ankara (MVA) / all enveloped viruses	+				

Ethasept®

Our alcoholic hand sanitizer for personal hygiene

Product properties:

Ethasept® is a very quick and comprehensive, ready to use, skin-friendly hand disinfectant **without perfumes and colourings**, which is perfect for hygienic hand disinfection. **Regenerating agents** soften the skin. Skin compatibility is confirmed by dermatological tests.

Areas of application:

For hygienic and surgical hand disinfection for all areas of human hygiene with high demands on hand hygiene: Clinics, medical practices, ambulances, care facilities, industrial and laboratory areas.

Type and duration of application:

Ready to use. Keep hands / forearms moist during application. Observe standard methods for hygienic (EN 1500) and surgical (EN 12791) hand disinfection.

Preparation properties:

- acts very quickly and thoroughly as a bactericide, yeasticide, tuberculocide
- very quick and effective against all enveloped viruses BVDV / vaccinia (incl. HCV / HBV / HIV)
- thoroughly effective against the non-enveloped viruses rota, adeno, polio and noro virus (virucidal effective according to EN 14476)
- effective against MRSA
- perfume and colouring-free
- very good skin compatibility - Dermatologically tested by: "very good"
- also examined for skin disinfection
- evaluated according to the standard methods of DGHM, EN 1500 and EN 12791 as well as the RKI recommendations for efficacy against viruses
- VAH listed
- IHO listed

Hygienic hand disinfection:

Reduction of transient germs in 15 seconds Rub in undiluted hands with at least 3 ml.

Surgical hand disinfection:

Reduction of transient and resident germs: 3 ml portions on hands for 90 seconds rub to keep hands moist during this exposure time.

REF	Box contents	Shipping unit	PZN
01003	150 ml	24 bottles	09102222
01004	500 ml	18 bottles	09102216
01005	1 l	10 bottles	09102191

Shelf life

Ethasept® can be used for up to 12 months after opening the original container for the first time. The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 and 25° C.

Composition:

100g contains:
Active ingredient: 73.5 g Ethanol; Other constituent materials: Purified water, moisturiser / surfactants.



Regulatory requirements

Ethasept® is not a medicinal product and is covered by the Biocides Directive. Use biocides safely. Always read the label and product information before use. VAH listed. Biocide of Product type 1 ("Biocidal Products for Human Hygiene") Biocide Reg. No. N-44236.

Application: Conc. / Exposure time

Efficacies	15 sec.	30 sec.	1 min. / 60 sec.	1,5 min. / 90 sec.	2 min. / 120 sec.
Hygienic hand disinfection: Bactericidal and fungicidal action according to EN 1500. Recommended application of VAH: at least 3 ml	+				
Surgical hand disinfection according to EN 12791: at least 3 ml / 90 seconds. Keep hands moist throughout the application.				+	
Skin disinfection: on low-sebaceous skin	+				
Skin disinfection: on sebum-rich skin (keep moist)			+		
tuberculocidal (including M. terrae)		+			
MRSA	+				
Limited virucidal according to RKI / DVV	+				
BVDV / vaccinia (including HBV, HCV, H IV)	+				
Virucidal according to EN 14476					+
Adeno			+		
Polio					+
Noro (MNV) (suspension test)	+				
Rota (suspension test)	+				
Vaccinia Ankara		+			
Virucidal according to RKI / DVV SV 40 (with and without load)			+		

Medizid® AF

Our universal, aldehyde-free surface disinfectant concentrate

Product properties and application areas

Liquid surface wipe disinfectant concentrate for the disinfection and cleaning of non-invasive medical devices and medical inventory as well as all wipe-resistant (upper) surfaces of medical devices. Ideal for areas with high levels of dirt and odour.

Application notes

Medizid® AF is a concentrate! Pay attention to application concentrations and exposure time. Completely wet surfaces. Use dosing aids. Note dosing instructions. After contact with food contact surfaces, clean product residues with clean water after exposure time. Use biocides safely. Always read the label and product information before use. Biocide of product types PT 2 and PT 4. Biocidal reg no.: N-71965. Limit the use of biocides to a minimum.

Application example

10 l of a 1% solution = 100 ml of disinfectant concentrate + 9.9 l of water. The required amount of water is obtained by subtracting from the amount of ready-to-use solution the proportion of the disinfectant concentrate used.

Shelf life

Medizid® AF can be used for up to 12 months after opening the original container for the first time. The minimum shelf life is 3 years (see label printing). The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 - 25°C.

Composition

100g contains:

5 g didecyldimethylammonium chloride.

Further components:

Non-ionic surfactants, synergistic additives, auxiliaries, complexing agents, fragrances.

Application concentrations / Exposure time

Efficacies	15min.	60min.	240min.
According to VAH / DGHM:			
Bactericidal, yeasticidal	3,0 % with load	1,0 % with load	0,5 % with load
Limited virucidal acc. to RKI:			
BVDV / vaccinia (including HBV, HIV, HCV)	0,5 %	2,0 % with load	
Virus efficacies:			
Noroviruses		2,0 %	
Adeno		3,0 %	
Polyoma / SV 40		2,0 %	
Tuberculocide (M.-terrae)		2,0 %	
MRSA	0,5 % with load		



Preparation properties

- bactericidal, yeasticidal, tuberculocidal effects (incl. M.-terrae)
- limited virucidal against all enveloped viruses: BVDV / Vaccinia (incl. HBV / HIV)
- virucide effective against the non-enveloped viruses: Adeno, polyoma / SV40 and noro viruses (FCV)
- effective against MRSA
- excellent cleaning properties
- no exposures to be monitored
- very good material compatibility
- tested according to the standard methods of the DGHM and the RKI recommendations on effectiveness against viruses
- VAH listed
- IHO listed

REF	Box contents	Shipping unit	PZN
00972D	2 l	6 bottles	00068819
00971D	5 l	2 canisters	00069090

Accessories (see page 159)

Fillable dosing bottle

Canister pump for dosing

Canister tap

Canister key

Medizid® Rapid

Our fast acting, ready-to-use spray wipe disinfectant for cleaning and disinfecting alcohol resistant surfaces of medical devices

Product type

Medizid® Rapid is a ready-to-use, very fast acting preparation for cleaning and disinfecting alcohol resistant surfaces and objects in medical facilities. Ideal for rapid disinfection between treatments and procedures when very short exposure times are required.

Applications

For all alcohol resistant surfaces and objects of medical products in close to patient surroundings in clinics, surgeries, ambulances and care areas.

Application

Ready to use. Moisten objects and surfaces completely with Medizid® Rapid, allow to dry and, if necessary, wipe with a disposable cloth. Dries quickly.

Shelf life

Medizid® Rapid can be used for up to 12 months after opening the original container for the first time. The minimum shelf life is 3 years (see label printing). The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 - 25°C.

Exposure time:

Efficacies	completely wetted
According to VAH / DGHM:	
Bactericidal, yeasticidal	15 sec.
Tuberculocide (including M. terrae)	30 sec.
MRSA / mycobactericidal	1 min.
Limited virucidal (according to DVV):	
BVDV / vaccinia (including HBV, HCV, HIV)	30 sec.
Virus efficacies:	
Rota	30 sec.
Noroviruses	3 min.
Polio / Adeno	15 min.
Polyoma / SV 40	1 min.

REF	Box contents	Shipping unit	PZN
01002D	250 ml	12 bottles	11169624
00975	500 ml	18 bottles	10014658
00976D	1 l	12 bottles	00068794
00977D	5 l	2 canisters	11169647

Accessories (see page 159)

Canister tap

Canister key

Spray head for 1 litre bottle



Composition

100g contains:

30 g of 1-propanol, 15 g of 2-propanol
10 g of ethanol, 0.05 g of didecylmethyl ammonium propionate, 0.02 g of alkylaminoguanidine diacetate

Further components:

Fragrances, non-ionic surfactants.

Preparation properties

- ready to use
- aldehyde-free
- fast and effective, already in 15 sec.
- bactericidal, yeasticidal, tuberculocidal effects (incl. M.-terrae), mycobactericidal (including M. avium)
- effective against all enveloped viruses: BVDV / Vaccinia (incl. HCV, HBV / HIV)
- effective against all non-enveloped viruses: Rota, adeno, polyoma / SV40, polio and noro viruses
- MRSA effective
- residue-free
- pleasant smell
- good material compatibility, even with sensitive materials
- be careful with alcohol sensitive surfaces such as acrylic glass!
- tested according to the standard methods of the DGHM and the RKI recommendations for effectiveness against viruses
- VAH listed
- IHO listed

Medizid® Rapid Disinfection wipes

Our ready-to-use, alcohol based hygienic wipes for disinfecting and cleaning medical device surfaces and medical device items

Medizid® Rapid disinfectant wipes are impregnated with the fast acting liquid of the alcohol based preparation Medizid® Rapid. The hygienic wipes are therefore very well suited for disinfecting and cleaning medical devices such as device surfaces and objects in the medical field and in the patient friendly environment between treatments and interventions.

Applications

For all alcohol resistant medical devices, surfaces and objects, in particular for device surfaces and objects near the patient in clinics, surgeries, ambulances and care areas.

Application and dosage

Medizid® Rapid disinfectant wipes are ready to use. The surfaces and objects to be disinfected are to be carefully wiped off with the Medizid® Rapid wipe and completely wetted. Exposure time: 60 sec. according to EN 16615 (4 field test).

Please note

Be careful with all alcohol sensitive materials!

Shelf life

The minimum shelf life of Medizid® Rapid disinfectant wipes is 3 years (see label). To prevent premature drying of the wipes, reseal the packaging / can / container immediately after use. Shelf life / expiry after opening: 28 days. The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 - 25° C.

Composition

In 100 g of impregnation solution are included:

30 g of 1-propanol, 15 g of 2-propanol, 10 g of ethanol, 0.05 g of didecylmethylammonium propionate, 0.02 g of alkylamine guanidine diacetate.

Further components: Fragrances, non-ionic surfactants.



Product properties

- effectiveness of the impregnating solution Medizid® Rapid
- ready to use
- aldehyde-free
- very fast and comprehensively effective in 15 sec.
- bactericidal, yeasticidal, tuberculocidal effect (including *M. terrae*), mycobactericidal, virus-inactivating (including HBV / HIV)
- effective against enveloped viruses: BVDV (HCV) and vaccinia according to the current recommendation of the RKI
- MRSA effective and tested (1 min.)
- residue-free
- pleasant smell
- good material compatibility
- virucide effective according to 14476 (Polio / Adeno)
- noro viruses and SV40 tested
- tested according to the 4-field test / EN 16615
- impregnating solution VAH listed
- impregnating solution IHO listed

REF	Item	Cloth size	Box contents	Shipping unit	PZN
00941D	Cloths in PE dispenser	approx. 15 x 18 cm	150 cloths	12 cans	11169593
00942D	Refill pack can	approx. 15 x 18 cm	150 cloths	20 bags	00068802
00946	Cloths in the container	approx. 26 x 29 cm	70 cloths	1 container	11169601
00947	Refill pack container	approx. 26 x 29 cm	70 cloths	6 rolls	11169618

Medizid® Rapid QF Flowpack

Our ready-to-use, alcohol based hygienic wipes for disinfecting and cleaning medical device surfaces and medical device items

Product type

Medizid® Rapid QF Disinfectant Wipes (Flowpack) are saturated with the very fast acting liquid of the alcohol based preparation Medizid® Rapid. The hygienic wipes are therefore very well suited for the disinfection and cleaning of medical devices such as device surfaces and objects in the medical field and in the patient friendly environment between treatments and interventions.

Applications

For all alcohol resistant medical devices, surfaces and objects, in particular for device surfaces and objects near the patient in clinics, practices, ambulances and care areas.

Application and dosage

Medizid® Rapid QF disinfectant wipes (Flowpack) are ready to use. The surfaces and objects to be disinfected are to be carefully wiped off with the Medizid® Rapid wipe and completely wetted. Exposure time: 60 seconds according to EN 16615 (4 field test).

Please note

Be careful with all alcohol sensitive materials!

Shelf life

The minimum shelf life of Medizid® Rapid disinfectant wipes is 3 years (see label). To prevent the premature drying out of the wipes, reseal the packaging immediately after use. Shelf life / expiry after opening: 28 days. The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 - 25° C.

Composition

In Medizid Rapid QF in 100 g, the active ingredient is: 50 g ethanol.



Product properties

- effectiveness of the impregnation solution Medizid® Rapid
- ready to use
- aldehyde-free
- very fast and comprehensively effective in 15 sec.
- bactericidal, yeasticidal, tuberculocidal effect (including M. terrae), mycobactericidal, virus-inactivating (including HBV / HIV)
- effective against enveloped viruses: BVDV (HCV) and vaccinia according to the current recommendation of the RKI
- MRSA effective and tested (1 min)
- residue-free
- pleasant smell
- good material compatibility
- virucide effective according to 14476 (Polio / Adeno)
- noro viruses and SV40 tested

REF	Cloth size	Box contents	Shipping unit	PZN
00976FP	approx. 22 x 20 cm	80 pcs	12 boxes	-

Medizid® Alpha

Our ready-to-use spray wipe disinfectant for cleaning and disinfecting all surfaces of medical devices, especially for alcohol sensitive surfaces

Product type

Medizid® Alpha is a ready-to-use, alcohol-free and aldehyde-free disinfectant for cleaning and disinfecting all wipe resistant surfaces, especially sensitive surfaces and objects. Can also be applied as foam with appropriate spray head adjustment. Also suitable for polyacrylic glass!

Applications

For medical facilities, sanitary and bathing areas and nursing homes.

Application

Apply undiluted, completely wet surfaces and objects, leave for 1 minute. Exposure time: CONC. / 1 MINUTE.

Shelf life

Medizid® Alpha can be used for up to 12 months after opening the original container for the first time. The minimum shelf life is 3 years (see label printing). The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 - 25° C.

Exposure time:

Efficacies	completely wetted
According to VAH / DGHM:	
Bakterizid, levurozid	1 min.
MRSA	1 min.
Begrenzt viruzid BVDV / vaccinia	1 min.

REF	Box contents	Shipping unit	PZN
00964	1 l	12 bottles	PZN 11169653

Accessories (see page 159)

Spray head



Composition

100g contains:
0.2 g Polyhexamethylene biguanide hydrochloride, 0.64 g Didecyldimethylammoniumchloride.

Preparation properties

- fast and comprehensive bactericidal, yeasticidal effect, already after 60 sec.
- effective against all enveloped viruses: BVDV / vaccinia (incl. HCV / HBV / HIV)
- effective against the non-enveloped viruses: rota, polyoma / SV40
- good cleaning performance
- good material compatibility
- suitable for nursing tubs
- suitable for sensitive surfaces, such as acrylic glass
- tested according to the standard methods of the DGHM and the RKI recommendations for effectiveness against viruses

Medizid® Alpha Disinfection wipes

Our ready-to-use, alcohol-free hygienic wipes for disinfecting and cleaning surfaces and objects

Product type

Medizid® Alpha disinfectant wipes are ready-to-use, alcohol-free hygiene wipes for cleaning and disinfecting small surfaces and objects. Ideal for use between treatments and interventions in healthcare facilities. Bactericidal, yeasticidal, virucidal effect against all enveloped viruses including HBV / HIV.

Applications

For all surfaces and objects, in particular for device surfaces and objects near the patient in clinics, surgeries, ambulances and care areas.

Application

Ready to use. Completely wet objects, surfaces, allow to dry, if necessary wipe. Only use damp cloths for disinfection.

Shelf life

The minimum shelf life of Medizid® Alpha disinfectant wipes is 3 years (see label printing). To prevent premature drying of the wipes, reseal the packaging / can / container immediately after use. Shelf life / expiry after opening: 28 days. The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 - 25° C.

Composition

In 100 g of impregnation solution are included: 0.2 g Polyhexamethylene biguanide hydrochloride, 0.64 g didecyldimethylammonium chloride.



Preparation properties

- ready to use
- fast and comprehensive, effective already after 60 sec.
- comprehensive bactericidal and yeasticidal effect
- effective against all enveloped viruses: BVDV / Vaccinia (incl. HCV / HBV / HIV)
- effective against non-enveloped viruses: Rota, Polyoma / SV40
- pleasant smell
- very good material compatibility even with sensitive materials, such as acrylic glass
- 100% synthetic material with optimum active ingredient intake
- favourable cloth format with approx. 15 x 18 cm
- impregnation liquid: VAH listed
- impregnation liquid: IHO listed

REF	Item	Cloth size	Box contents	Shipping unit	PZN
00969	Cloths in the container	26 x 29 cm	70 cloths	1 container	11169676
00970D	Refill pack	26 x 29 cm	70 cloths	6 rolls	11169682

Medizid[®] Alpha Flowpack

Our ready-to-use, alcohol-free hygienic wipes for disinfecting and cleaning surfaces and objects

Product type

Medizid[®] Alpha disinfectant wipes (Flowpack) are ready-to-use, aldehyde and alcohol-free hygiene wipes for cleaning and disinfecting small surfaces and objects. Ideal for use between treatments and interventions in healthcare facilities. Bactericidal, yeasticidal, virucidal effect against all enveloped viruses including HBV / HIV. Medizid[®] Alpha disinfectant wipes (Flowpack) are impregnated with the very fast acting liquid of non-alcohol based Medizid[®] Alpha.

Applications

For all surfaces and objects, in particular for device surfaces and objects near the patient in clinics, surgeries, ambulances and care areas.

Application notes

Ready to use. Completely wet objects, surfaces, allow to dry, if necessary wipe. Only use damp cloths for disinfection.

Shelf life

The minimum shelf life of Medizid[®] Alpha disinfectant wipes (Flowpack) is 3 years (see label). To prevent the premature drying of wipes, close the packaging immediately after use. Shelf life / expiry after opening: 28 days. The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 - 25°C.

Composition

In 100 g of impregnation solution are included: 0.2 g Polyhexamethylene biguanide hydrochloride, 0.64 g Didecyldimethylammoniumchloride.



Product properties

- ready to use
- fast and comprehensive, effective already after 60 sec.
- comprehensive bactericidal and yeasticidal effect
- effective against all enveloped viruses: BVDV / Vaccinia (incl. HCV / HBV / HIV)
- effective against non-enveloped viruses: Rota, Polyoma / SV40
- pleasant smell
- very good material compatibility even with sensitive materials, such as acrylic glass
- 100% cellulose

REF	Cloth size	Box contents	Shipping unit	PZN
00964FP	approx. 22 x 20 cm	80 pcs	12 boxes	-

Instrusol® AF +

Our fully effective, aldehyde and phenol-free instrument disinfectant for manual reprocessing in the insertion process and in the ultrasonic cleaning device

Product type

Instrusol® AF + is suitable for the treatment of all thermostable and thermolabile instruments such as anaesthetic accessories and rigid endoscopes. An innovative formulation ensures fast, comprehensive and material saving instrument reprocessing with low application concentrations.

Applications

For all applications in which medical instruments are prepared in the insertion process and in the ultrasonic cleaning device.

Application notes

The instructions for use and reprocessing guidelines of the instrument manufacturers must be observed; in particular the instructions for the preparatory measures, for cleaning and disinfection, rinsing, drying and storage.



Application concentrations / Exposure time

Efficacies	5 min.	15 min.	30 min.	60 min.
According to VAH / DGHM:				
Bactericidal, yeasticidal, tuberculocidal (including <i>M. terrae</i>)		4,0 %	3,0 %	2,0 %
Limited virucidal (according to DVV): BVDV / vaccinia (including HBV, HCV, HIV)		2,0 %		
Virus efficacy:				
Adeno / Noro				4,0 %
SV 40			3,0 %	
Rota	2,0 %			
MRSA		3,0 %	2,0 %	1,0 %
Ultrasound bath	3,0 %			

REF	Box contents	Shipping unit	PZN
00973	2 l	6 bottles	11169699
00974	5 l	2 canisters	11169707

Accessories (see page 159)

Canister pump for dosing
Canister tap
Canister key

Composition

100g contains:

5 g of alkyldimethylethyl ammonium ethylsulfate
1 g Polyhexamethylene biguanide hydrochloride
1.8 g Coco propylene guanidinium diacetate

Further components: Surfactants, corrosion inhibitors, and suspending agents adjuvants and pH regulators.

Preparation properties

- aldehyde and phenol-free
- comprehensively effective
 - bactericidal, yeasticidal, tuberculocidal effect (including *M. terrae*)
 - virucide effective against all enveloped viruses: BVDV / Vaccinia (incl. HCV / HBV / HIV)
 - virucide effective against the non-enveloped viruses: rota, adeno, polyoma / SV40, noro
- economic
- good material compatibility
- pleasant odour, no emissions to be monitored (workplace values)
- good biological degradability
- tested according to the standard methods of the DGHM and the RKI recommendations
- VAH listed
- IHO listed

Processing instructions

Insertion process

Insert instruments immediately after use. All surfaces and cavities must be covered or filled.

Please observe the application concentrations and exposure times. Remove instruments, rinse thoroughly with water and dry.

Ultrasonic bath

The prescribed frequency range and the specified duration of disinfection must be observed.

Water temperatures above 45° C should be avoided.

Preparation change

Instrusol® AF + is not compatible with aldehyde containing disinfectants. During the first use and change of the preparation, the instruments and the accessories, the tubs and preparation devices as well as the application aids must be thoroughly cleaned and rinsed.

Instrusol® AF + use solution is protein resistant and can be used several times if necessary.

The working solution can be used for up to 7 days when inserting well pre-cleaned instruments.

In the event of cloudiness or visible contamination, it is recommended to replace the working solution sooner.

Stability / water quality

In the production of working solutions, drinking water quality should be used. For very hard water we recommend distilled water. For the last rinse before steam sterilisation we recommend demineralised water or distilled water. For semi-critical medical devices that are not sterilised, germ-free, preferably sterile, water should be used for final rinsing.

Instrument cleaning

When adding cleaners, the working solution should be renewed daily.

Shelf life

Instrusol® AF + can be used for at least 12 months after opening the original container for the first time.

The minimum shelf life is 3 years (see label printing). The products must only be stored in the original sealed container, contamination must be avoided while removing product, and the products must be kept and stored in accordance with standard storage conditions (protected from damp, direct sunlight, heat and frost) at a temperature of between 5 - 25°C.

Note

As a working solution Instrusol® AF + is not a dangerous preparation according to the Ordinance on Hazardous Substances. The working solution is not corrosive. However, we recommend that you follow the usual safety precautions when handling chemicals and wear the usual protective clothing.

Instrusol® Enzyme

Our enzymatic instrument cleaning for endoscopy and surgery

Product type

Instrusol® Enzyme is an enzyme cleaner for the careful, residue-free cleaning of medical instruments, especially flexible and rigid endoscopes, endoscope accessories and surgical and orthopaedic instruments. Highly active enzymes, supported by an innovative surfactant system, eliminate deposits of blood, secretions, detritus, tissue residues and fats. Instrusol® enzyme prevents biofilms that harbour resistant bacteria and interfere with hygienic treatment. With effective, residue-free cleaning in 2-3 minutes, Instrusol® Enzyme creates the condition for a safe disinfection and sterilisation. Instrusol® Enzyme prevents coagulation, thus reducing device failures and ensuring the functionality of moving instrument parts. The excellent properties facilitate the residue-free cleaning, especially with Englumigem instruments, endoscopes and endoscopic accessories. Instrusol® Enzyme is foam-free and ideal for ultrasonic preparation. Instrusol® provides a microbial protection of the use solution by bactericidal and virus-inactivating properties and thus ensures maximum user safety.

Applications

For all areas of application in which medical instruments can be carefully, thoroughly and residue-free cleaned in the insertion process, circulation process and in the ultrasonic bath.

Application notes

The instructions for use and reprocessing guidelines of the instrument, endoscope, endoscope accessories and device manufacturers must be observed, in particular the preparatory measures for cleaning.

Composition

Proteases and specifically amylases, glycols, non-ionic surfactants, preservatives, fragrances, purified water.

REF	Box contents	Shipping unit	PZN
02001B	500 ml	6 bottles	13816269
02002B	1 l	6 bottles	13816275
02003B	2 l	6 bottles	13816281
02004B	5 l	2 canisters	13816298

Accessories (see page 159)

Canister pump for dosing
Canister tap
Canister key



Preparation properties

- effective even at cold water temperatures; warm water improves the cleaning effect (30° - 45° C)
- prevents biofilms that harbour bacteria resistant to disinfectants and can cause infections
- a bactericidal and virus-inactivating component protects the handling solution contaminated with dirt
- yeasticidal
- Instrusol® Enzyme does not replace any disinfectant

Application and dosage:

0.5 - 5% depending on the degree of soiling and the type of application. If the instrument is cleaned after each use with Instrusol® Enzyme, 1 - 3% as application concentration is sufficient. If necessary, the application concentration can be increased to 5%. The working solution can be used several times, however it must be renewed regularly. Under normal conditions, the working solution should be renewed once or twice a day. In difficult conditions and low application concentrations, renewal is recommended after each treatment.

Processing instructions

Insertion process

Manual preparation in the insertion process and ultrasonic bath: Insert instruments after use. All surfaces and cavities must be covered or filled. Thoroughly clean the instruments, observing the instructions of the instrument manufacturer. For encrustations, let Instrusol® Enzyme act and preferably use warm water. In the ultrasonic bath, observe the specified frequency and duration of cleaning. Avoid water temperatures over 45° C. Then remove the instruments and disinfect or rinse thoroughly with water, dry well and sterilise. For disinfection we recommend Instrusol® AF + or InstruPlus® forte.

Circulation procedures

Insert instruments after use. Follow the instructions of the device manufacturer.

Compatibility / stability / water quality:

Instrusol® Enzyme is compatible with the aldehyde-free instrument disinfectants Instrusol® AF +. Mixing with aldehyde-containing disinfectants generally is not useful. Mixing with anionic surfactants and amphoteric surfactants is not recommended. Instrusol® Enzyme should preferably be kept cool. Storage temperatures above room temperature and direct sunlight should be avoided. Instrusol® Enzyme is effective at all water temperatures, including cold water. Warm water accelerates the cleaning performance. The optimum temperature range is between 30 and 45° C. For the production of working solutions, drinking water quality water is required. For very hard water we recommend distilled water. For the final rinse we recommend preferably germ-free or sterile water. For the last rinse before steam sterilisation we recommend DI water or distilled water.

Chemical-thermal endoscope reprocessing

If endoscopes are processed both manually and automatically, InstruPlus® forte can be used for manual reprocessing.

Biodegradability

Instrusol® Enzyme is very biodegradable. The working solution can be disposed of with general waste water. A disturbance of the biological treatment plants is not produced. Instrusol® Enzyme is not classified as environmentally hazardous according to EC Directives. Shelf life after opening: 6 months

InstruPlus® forte

Our highly effective instrument disinfectant with virucidal effectiveness for thermolabile and thermostable instruments.

Product type

InstruPlus® forte is an innovative aldehydic preparation based on acetals and aldehydes. InstruPlus® forte with low application concentrations is excellently compatible with materials and is particularly suitable for the treatment of very sensitive instruments, such as flexible endoscopes, endoscopic additional instruments, sensitive optics and anaesthesia accessories. The excellent material compatibility is well proven. InstruPlus® forte is suitable for manual reprocessing in the insertion process and for circulation processes in disinfection equipment. InstruPlus® forte is comprehensively assessed in accordance with the DIN EN standards for instrument disinfection, the standard methods of the DGHM and the RKI / DVV recommendations for the "Testing and declaration of disinfectants against viruses".

Applications

InstruPlus® forte is suitable for all areas in which medical instruments are processed: In the inserting and circulation process at normal water temperatures (25° C).

Application notes

The instructions for use and reprocessing instructions of the instrument, endoscope, endoscope accessory and device manufacturers must be observed, in particular the instructions for the preparatory measures, for cleaning and disinfection, for rinsing and drying as well as for storage.

Semi-critical medical devices that are not disinfected or sterilised thermally (92° C - 95° C) must be treated as virucidal.

Composition

100g contains: 5.75 g formacetal, 8.0 g glutaric dialdehyde (pentanedial)

Further components: Surfactants, corrosion inhibitors, preservatives, suspending agents and adjuvants as well as pH regulators.



Preparation properties

- economical due to very low application concentrations (1-3%)
- comprehensively effective
 - bactericidal, yeasticidal / fungicidal (incl. A. niger) tuberculocidal (incl. M. terrae)
 - effective against all enveloped viruses (limited virucidal), including HBV, HIV, HCV
 - effective against all non-enveloped viruses (virucidal): adeno, polyoma / SV40, polio
 - MRSA effective
- excellent material compatibility
- low exposures
- good biological degradability
- use solution can be used multiple times (up to 7 days)
- VAH listed
- IHO listed

REF	Box contents	Shipping unit	PZN
02005B	100 ml	6 bottles	13816909
02006B	2 l	6 bottles	13816915
02007B	5 l	2 canisters	13816921

Accessories (see page 159)

Canister pump for dosing

Canister tap

Canister key

Application concentrations / Exposure time:

Efficacies	15 min.	30 min.	60 min.
According to VAH / DGHM:			
Bactericidal, yeasticidal, fungicidal, tuberculocidal tuberculocidal (including <i>M. terrae</i>)	1,5 %	1,0 %	1,0 %
Limited virucidal (according to DVV): BVDV / vaccinia (including HBV, HCV, HIV)	1,5 %	1,0 %	
Virucidal (according to DVV):	3,0 %		
Adeno	1,0 %		
Polyoma / SV 40	2,0 %	1,0 %	
Polio	3,0 %		
MRSA	1,5 %	1,0 %	

InstruPlus® forte use solution is protein resistant and can be used several times if necessary. The working solution can be used for up to 7 days when inserting well pre-cleaned instruments. In the event of cloudiness or visible contamination, it is recommended to replace the working solution sooner. If surfactant cleaners are used, they must only contain non-ionic surfactants and should be neutral pH formulated.

Stability / water quality

In the production of working solutions, drinking water quality should be used. For very hard drinking water we recommend distilled water. For the last rinse before steam sterilisation we recommend DI water or distilled water. For semi-critical medical devices that are not sterilised, germ-free, preferably sterile, water should be used for final rinsing.

Instrument cleaning

For the pre-cleaning of endoscopes, accessories, surgical instruments, especially narrow bore instruments. For heavy soiling we recommend the enzyme cleaner Instrusol® Enzyme. If cleaning agents are added, the working solution should be renewed daily.

Shelf life

InstruPlus® forte can be used 12 months after onset.

Biodegradability

InstruPlus® forte is readily biodegradable as an aldehyde formulated product. The working solution can be disposed of with general waste water. If properly initiated, it is unlikely that the biological treatment system will be disturbed. In accordance with OECD guidelines, InstruPlus® forte is readily biodegradable.

Processing instructions**Insertion process**

Insert instruments immediately after use. All surfaces and cavities must be covered or filled. Observe the application concentrations and reaction times. Remove instruments, rinse thoroughly with water and dry.

Circulation procedure with endoscopes (disinfection devices)

The instructions of the device manufacturer must be observed, in particular the notes on dosing settings, program sequences and product changes. After disinfection, attention must be paid to sufficient final rinsing and adequate drying.

Preparation change

InstruPlus® forte is not compatible with aldehyde-free disinfectants. During initial use and when changing the specimen, thoroughly clean and rinse the instruments, accessories, baths and preparation equipment as well as the application aids. For support in the pre-cleaning we recommend Instrusol® Enzym.

Gentle Med[®] wash lotion

pH skin neutral and nourishing washing lotion

Characteristics

A washing lotion dermatologically tested by "excellent" for daily gentle full body application especially for sensitive skin, with fresh scent of chamomile extract. Prevents the skin from drying out. Mild surfactant combination. Paraben-free Alkali-free, free of silicone and paraffin oil. Neutral skin. Economical in use due to optional dosing pump.

Not tested on animals.

Application

For all areas with increased hygiene requirements in hospitals, surgeries, nursing homes or private households.

REF	Box contents	Shipping unit	PZN
00980D	500 ml	12 bottles	11170314
00981D	1 l	10 bottles	11170320
00982D	5 l	2 canisters	11170337



Gentle Med[®] MRSA effective wash lotion

pH skin neutral and nourishing wash lotion MRSA

Characteristics

Dermatologically tested as "excellent", colourless washing lotion for daily gentle full body application in the case of MRSA. The fast and comprehensive effect occurs after just 30 seconds. Economical in use due to optional dosing pump. Not tested on animals.

Application

For all areas with increased hygiene requirements in clinics, surgeries, ambulances or nursing homes in the case of MRSA. The washing lotion is suitable for the antiseptic treatment of skin, hands and whole body including head hair in MRSA / ORSA colonisation.

REF	Box contents	Shipping unit	PZN
00961D	500 ml	12 bottles	11170277
00962D	1 l	10 bottles	11170283



Gentle Med[®] Skin care lotion (O / W)

The skin care lotion for protection and care stressed skin, with urea

Characteristics

Skin care lotion for sensitive and dehydrated skin, pH skin neutral, enriched with urea. Absorbs quickly and produces a non-greasy protective film with a simultaneous vitalisation of the skin. A new and particularly effective formulation that nourishes, moisturises and smooths stressed skin with precious natural oils (almond and jojoba oils), urea and hyaluronic acid.

With a light and pleasant scent. Dermatologically tested as "excellent". PEG-free. Free of silicone and paraffin oil. Not tested on animals.



REF	Box contents	Shipping unit	PZN
01011	100 ml	30 tubes	11170308
00990D	500 ml	12 bottles	11170231

Gentle Med[®] Cream lotion (W / O)

Barrier cream and skin protection balm for very dry and heavily stressed skin

Characteristics

Cream lotion with a particularly nourishing formula to protect and support the regeneration of extremely stressed skin. The active ingredients used support the regeneration of brittle and chapped skin with a moisture repellent function. Gentle Med[®] cream lotion is dermatologically tested as "excellent" and is suitable for whole body care, pH skin neutral enriched with panthenol and vitamin E, and is thereby soothing. Free of silicone oil. Not tested on animals.



REF	Box contents	Shipping unit	PZN
01013	100 ml	30 tubes	11169794
00988D	500 ml	12 bottles	11170053

Gentle Med[®] Active gel

Vitalising active gel with essential oils
(menthol, camphor, rosemary)

Characteristics

Active gel for circulation stimulating, invigorating and skin friendly rubbing in. Very suitable for muscle relaxation. Deodorising and cooling. Very suitable for massage in case of tension. The sustainably cooling gel provides a soothing sensation of freshness and vitality. Has a refreshing effect on the respiratory tract. Dermatologically tested as "excellent". Alkali-free. Paraben-free PEG-free. Free of silicone and paraffin oil. Only natural fragrances. Not tested on animals.

REF	Box contents	Shipping unit	PZN
00991D	500 ml	12 bottles	11169788



Gentle Med[®] Oil bath

Oil bath with natural oils

Characteristics

Ready-to-use, medicinal oil bath enriched with high quality soya bean oil. Promotes the physiological structure of the skin. Very beneficial for sensitive and itchy skin. Suitable for shower and full baths. Not tested on animals.

REF	Box contents	Shipping unit	PZN
00992D	500 ml	12 bottles	11170260



Gentle Med[®] Skin cleansing and care foam

Gentle cleaning and care with reliable skin protection

Characteristics

Fast and gentle cleansing as well as skin care are of main importance for patients with incontinence, anus praeter and those that are bedridden. The skin cleansing care foam contains moisturising and odour binding substances. The fine pored foam offers a very good odour coverage, an excellent cleaning without irritation and leaves a fresh and nourished skin feeling. Paraben-free PEG-free. Free of silicone and paraffin oil. Without allergenic fragrances. Not tested on animals.

Application

Gentle, thorough, odour binding cleansing, especially for incontinence.

For foaming, hold the can (valve) diagonally down to ensure the residue is emptied.

Allow to act for a short time and absorb with disposable cloths. A nourishing and deodorising protective film remains. Shake the can vigorously before use.



REF	Box contents	Shipping unit	PZN
00994D	500 ml	12 bottles	11170248

Gentle Med[®] Moist care wipes

Our high quality wipes for normal and sensitive skin

Characteristics

- moist care wipes with witch hazel extract (natural plant extract for additional care)
- natural Betaine protects against dehydration of the skin and prevents skin irritation
- texture based on mild coconut oil-based surfactants
- PEG-free
- paraben-free
- paraffin and silicone oil free
- resealable plastic closure for better and more hygienic removal; protects the wipes from drying out



Application

For effective, fast and particularly gentle care and skin cleansing of soiled areas of the skin (e.g. incontinence). Protection against germ transmission by single use. Use without water. Warning: Please close the packaging carefully after removing the wipe.

REF	Size	Box contents	Shipping unit	PZN
2356	approx. 18 x 20 cm	80 pcs	12 boxes	12357658



Application aids

Our application aids for the removal and dosing of disinfectants and care products

Dosing pumps / attachments

Allow convenient and hygienic removal of products from original packaging.

REF	Item		PZN
00929	Dosing pump	500 ml	-
00930	Dosing pump	1 l	-
00931	Canister dosing pump	2 l, 5 l and 10 l	-

(Foam) spray heads

allow a uniform wetting of even hard to reach areas, e.g. in surface spraying of surface disinfectants. Also suitable for overhead application.

REF	Item	PZN
00934	Foam spray head 1 l	-

Canister tap

Enables the clean and drip-free dispensing of products from canisters.

REF	Item	PZN
00996	Canister tap 2 l, 5 l and 10 l	-

Canister key

For effortless opening of canisters.

REF	Item	PZN
00932	Canister key 5 l and 10 l	-

Dosing bottle / measuring beaker

For user friendly dispensing during transfer or for the preparation of working solutions or concentrates.


REF	Item		PZN
00943	Dosing bottle	1 l	-
00945	Measuring beaker	1000 ml	-
00927	Measuring beaker	2000 ml	-



Disinfectant dispenser

For the dosed dispensing of liquids

REF	Item		PZN
00936KST	Plastic dispenser	500 ml	-
00937KST	Plastic dispenser	1000 ml	-
00938AL	Aluminium dispenser	500 ml	-
00939AL	Aluminium dispenser	1000 ml	-



Material

Area A: ASA (acrylonitrile-styrene-acrylic ester)

Area B: ABS (acrylonitrile-butadiene-styrene)

Area C: ABS (acrylonitrile-butadiene-styrene)

Area D: SAN (styrene-acrylonitrile)



Our hygiene and cleaning products support you in all areas

One of the primary aims in the healthcare sector and other business facilities is achieving professional cleaning standards in all areas. We will help you achieve hygiene standards with our wide range of cleaning and hygiene products.

Note:

The recommendations for our products are based on scientific tests and are given to the best of our knowledge. Further recommendations, e.g. in terms of material compatibility, are only possible separately. Our recommendations do not exempt our customers from testing our products independently to ensure they are suitable for the intended purposes and procedures. In that regard, we can not assume any liability.

beopur WC fresh

Acidic, ready to use and viscous toilet cleaner

Characteristics

WC fresh cleaner for quick and thorough cleaning of toilets, urinals, tubs, shower cubicles and all acid resistant surfaces in sanitary areas.



Application

For stubborn limescale deposits, for limescale loosening and for intensive cleaning in heavy soiling, also able to be used undiluted. Rinse thoroughly with water afterwards.

REF	Box contents	Shipping unit	PZN
00926	750 ml	15 bottles	-

BioSan® Fresh scent

Highly concentrated for the removal of unpleasant and bad odours

Characteristics

This biologically fully degradable, high quality high concentrate easily removes any unpleasant odour. Safe for the environment




Application


Area of application In areas with strong odour nuisance. For unpleasant smelling rooms, spray as needed.

REF	Box contents	Shipping unit	PZN
6854	500 ml	6 spray bottles	-

Reliable hygiene in patient rooms



Compliance with strict hygiene requirements in patient rooms is justifiably a high priority. KLEEN PURGATIS offers a complete range of cleaning and disinfecting agents which make it easy to comply with legal regulations and pass hygiene tests. The application areas of the products are clearly defined. Cleaning staff and patients can be sure of safe hygiene.



Reliable hygiene in the sanitary sector

The sanitary sector is considered a hygiene challenge in all areas. This is especially true in health and geriatric care. KLEEN PURGATIS supplies hygiene safety in the sanitary sector with the high quality and user friendliness of its products.

**Shelves, washbasins, taps
PREMIUM No 1**

**Chemo-thermal laundry disinfection at 60°C (RKI / VAH)
LAVO DES 60 KOMPAKT**

**Radiators
PREMIUM No 1**




























**Basic disinfection
Medizid® AF**

**Rapid disinfection
Medizid® Rapid**

**WC, floor, wall tiles, toilet brush
PREMIUM No 1**

**Mirror
ULTRA GLASREIN**

Cleaning and hygiene plan

what	with what	how	hazardous substances
 Disinfect hands	 Alcoman® + Hand disinfection, VAH listed biocide	For hygienic hand disinfection rub 3-5 ml Alcoman® / Alcoman® + into the hands, make sure to wet completely	
 Handles, switches	 Medizid® Rapid Disinfectant cleaner, VAH listed	If necessary, according to specification / instructions	
 Basic disinfection	 Medizid® AF Surface disinfectant cleaner, VAH listed	If necessary, according to specification / instructions	
 Floor cleaning	 BLUE STAR For all waterproof surfaces and floors	clean, if necessary wipe dry, or with a machine	
 Furniture surfaces	 BLUE STAR For all waterproof surfaces and floors	clean, if necessary wipe dry	
 Glass and mirror	 ULTRA GLASREIN For all glass and mirror surfaces	spray on, wipe dry with a cloth	
 Textile surfaces	 NOVATEP For all water and colour resistant textile surfaces	e.g. spot cleaning, spray extraction, stain removal	
 Laundry disinfection	 LAVO DES 60 KOMPAKT Chemo-thermal laundry disinfection at 60°C, RKI / VAH listed	Introduce into the rinsing tray of the main washing compartment by means of the dosing system or according to the size of the washing machine and the water hardness according to the dosing table.	
 Taps	 PREMIUM No.1 Sanitary maintenance cleaner	apply, allow to act, rinse well with water	
 Toilets, urinals	 PREMIUM No.1 Sanitary maintenance cleaner	apply, allow to act, rinse well with water	
 Sanitary floor	 PREMIUM No.1 Sanitary maintenance cleaner	apply, allow to act, rinse well with water	

Cleaning in commercial kitchens



only vertical stainless steel surfaces
PRESTAN STAINLESS STEEL CARE




Commercial dishwashers
PRESTAN INTENSIV
PRESTAN KLAR FORTE



Washing up liquid
VANULAN



Ground
SPONTAN SUPER, PURINA



Grill
GRANIT S



Cutting appliances
Food preparation work surfaces
Medizid® Rapid

Are you a professional in the food service sector? Using professional cleaning products from KLEEN PURGATIS offers you the opportunity to comply with all statutory hygiene regulations with minimum expense. You can print out all relevant data for the products listed at www.kleen-purgatis.de. With these you lay the groundwork to be able to face every hygiene test.

Stoves and surroundings
PURINA

Water heaters
PARA LIMESCALE SOLVENT

General surface disinfection
Medizid® AF

Extractor hoods
PURINA

Convection ovens
GRANIT S

Hand disinfection
Alcoman®

Refrigerators
Medizid® Rapid

Cleaning in other areas

The taproom and the sanitary areas are the calling card of any establishment. Also here, KLEEN PURGATIS offers a large number of optimum products for use in commercial operations. Show your guests your business in a perfect, hygienically clean condition. The guests will thank you and feel good.

Glass and plastics
ULTRA GLASREIN

All waterproof surfaces and floors
BLUE STAR

Laminate and parquet cleaner
Premium N° 3

Taproom

Disinfection dispensing equipment
Medizid® AF

Glass dishwasher
Glass detergent

Ideal hand detergent
Vanulan

Mirror, glass
ULTRA GLASREIN

Washing cream for Euro soap dispensers
COPELIA GASTRO

Floors, tiles, ceramic surfaces
PUSH AMI

Sanitary areas

Shelves, supports, partitions (e.g. Resopal)
BLUE STAR








































Disinfection
Medizid® Rapid

WC and urinals
TRO WC-TABS

PUSH WC-CLEANER

172

Cleaning and hygiene plan

what	with what	how	hazardous substances
 Disinfect hands	 BUEDESIN® LOTIO HD VAH listed	hygienic: Rub in 3 ml / 30 sec undiluted, lather with a little water, clean hands, rinse with water.	
 Dishes, cutlery	 PRESTAN® INTENSIV Dishwashing detergent	machine dishwashing	
 Dishes, cutlery	 PRESTAN®KLAR FORTE Rinse aid	machine dishwashing	
 Dishes, pots	 VANULAN Washing up liquid - ideal for glasses -	manual dishwashing	
 Stove, extractors, surroundings	 PURINA Degreaser alkaline	Apply cleansing solution evenly, allow to act, work in, rinse off and, if necessary, wipe dry	
 Water heaters	 PARA LIMESCALE SOLVENT Descaler and decrusting agent	Add approx. 100 ml descaler per litre of water. Fill the water tank of the appliance with the solution 2/3. Heat up the appliance. Important! After each descaling, rinse the appliance sufficiently.	
 Surfaces, floors	 Medizid® AF Surface disinfection cleaner, VAH listed	Pre-clean surfaces after each use, apply the disinfectant solution, clean, reaction time min. 15 min	  
 Floor cleaning	 Change cleaning PURINA (5x / wk), SPONTAN SUPER (2x / wk)	Apply product solution, work on it with a scrubber, rinse and wipe dry	
 Stainless steel surfaces	 PRESTAN® STAINLESS STEEL CARE only vertical surfaces, not for food processing surfaces	Apply with clean, dry cloth; distribute the product evenly on the surface, do not polish again; do not use on surfaces for food processing	
 Convection oven, grill	 GRANIT S Grill cleaner with metal protection  	If necessary, pre-heat the surface, pre-rinse, pre-steam, spray on the product, allow to act	
 Surfaces, appliances	 Medizid® Rapid Surface disinfection VAH listed	After each use, pre-clean and disinfect equipment and surfaces, spray on clean and apply evenly, reaction time 60 sec	  



Everyday products, which is why they are essential

Healthcare workflows should primarily focus on the patient. That is why we have an extensive and well thought out range of ward and patient supplies in our portfolio to meet this need on a daily basis. Whether it be using our care products and wash gloves for the daily hygiene of immobile patients, or the essentials like cotton buds, wooden tongue depressors or protective towels.

Our innovative disposable BagBath® body wash system which has proven itself millions of times to patients and nursing staff is particularly noteworthy.

BagBath® Haircare

Hair washing without water

Characteristics

Innovative hair washing system without rinsing the hair. For easy handling and higher cost efficiency. Releases any tangled hair and makes the hair easier to comb.

Can be used both heated and cooled. Contains no fragrance and is therefore particularly skin friendly. Excellent for use in combination with the BagBath body wash system.



Application

Skin friendly hair washing for bedridden and people in need of care.

REF	Box contents	Shipping unit	PZN
APO021012	100 ml	30 bottles	-

BagBath®

Innovative disposable body wash system without soap and water

Characteristics

The wipes contained in BagBath® have a pleasant and very subtle fragrance. They consist of a special fabric blend that guarantees high tear resistance. Moisturising ingredients and clinically tested formulas for the ingredients ensure the skin is protected against irritation and promote its regeneration. Additional creaming is no longer necessary after use.



Application

Easy handling for the nursing staff and the patient. The box is easy to open and close, even with gloves on.

No more need for wash basins, soap, water, wash cloths, towels and care lotions! Eye and mucous membrane compatible.


REF	Size	Box contents	Shipping unit	PZN
2359	25 x 20 cm	8 cloths	30 boxes	06149252

BagBath® enables hospitals and care centres to improve the quality of the skin care of their patients while saving time and money.

By applying it without soap and water, the natural acid mantle of the patient's skin is preserved and no bacterial nutrient remains on the skin due to soap residue.

Cover all body areas:

BagBath® ensures complete body hygiene.
One box contains 8 cloths - one for each body region.



Benefits for the patient:

- Avoidance of smear infection by conventional washing utensils such as cloths and towel.
- Prevents cross-contamination by eliminating reusable laundry materials.
- Washing with BagBath® poses less risk to operated, intubated, traumatised, and immobile patients, as there is no need for additional transfer due to wet bedding.

Special, heatable packaging:

→ Microwave + BagBath®

- Thermal diffusion ensures a deeper cleansing of the skin.
- Most patients prefer the feeling of a warm bath.
- Warmth prevents the patient from cooling down.



Special, coolable packaging:

→ Refrigerator + BagBath®

- Refreshing and cooling effect on patients in summer temperatures.
- Intensive care units can thus gently wash patients with fevers.



Easy disposal:

→ Refuse bin

As a disposable product, BagBath® can be easily and simply disposed of along with the normal household waste. When washing quarantined patients, the in-house disposal regulations must be adhered to.



The BagBath® cleansing wipes contain a pleasant and vitalising fragrance. They consist of a special fabric blend that guarantees high tear resistance. Moisturising ingredients and clinically tested formulas of ingredients ensure the protection of the skin from irritation and promote its regeneration. Additional creaming is no longer necessary after use.

BeeSana® Wash glove

Hygienic disposable wash glove

Characteristics

Disposable wash glove in various designs and qualities according to the individual requirements of users and patients as well as those in need of care. Maintain their dimensional stability even when wet and have impressive abrasion resistance.



Application

For skin cleansing and body care. Applicable in all nursing areas, especially nursing homes and departments at risk of infection.

REF	Description	Execution	Box contents	Shipping unit	PZN
3709	Foil	Ultrasonically sealed	50 pcs	20 boxes	05461042
3723	Molton / soft	Ultrasonically sealed	50 pcs	20 boxes	11703466
3724	Molton / soft	stitched	25 pcs	40 boxes	05529001
3725	Molton 80 g	stitched	50 pcs	20 boxes	09336931

Protective towel

Characteristics

Protective towel made from a layer of cellulose and PE film. The PE inner coating prevents the penetration of moisture on clothing. With punched straps for fixation around the neck and a foldable collection bag in the lower area. Due to the soft cellulose outer layer, these can also be used as napkins after eating. Free of latex, PVC and DEHP, colour: Apricot.

Application

For single use.

REF	Size	Box contents	Shipping unit	PZN
2997	38 x 68 cm	100 pcs	6 boxes	11727314



Protective bibs

Characteristics

- protective bib made from medical gauze
- 6 x
- white

REF	Size	Box contents	Shipping unit	PZN
4239	20 x 30 cm	10 pcs	10 boxes	12724401



BeeSana® Non-woven disposable panty

Our pants for secure intimate areas

Characteristics

Disposable pants made of soft non-woven fabric

- with rubber spring
- white, with coloured cuffs for size differentiation



Application

For patients to preserve their privacy during examinations or after operations.

REF	Size	Waistband colour	Waist size	Box contents	Shipping unit	PZN
2251	M	blue	51 - 55 cm	20 pcs	30 boxes	11857318
2252	L	brown	54 - 58 cm	20 pcs	30 boxes	11857324
2253	XL	green	58 - 62 cm	20 pcs	30 boxes	11857330
2254	XXL	red	62 - 64 cm	20 pcs	30 boxes	11857301

Anti-slip sock universal

Our sock for fall prevention

Characteristics

- ABS disposable sock (laid)
- with nubs as anti-slip coating
- made of cotton
- for pleasantly warm and dry feet
- adapts perfectly to every foot
- black
- universal size
- washable at 30°C



Application

For single, versatile use in hospitals and nursing homes, especially in all areas, which require high hygiene standards.

REF	Size	Box contents	Shipping unit	PZN
9605	40 x 9 cm	20 pcs	8 boxes	12671811

Surgical sock

Our cotton sock for versatile use

Characteristics

- single sock (laid) without heel
- cuffs with rubber band
- for pleasantly warm and dry feet
- adapts perfectly to every foot
- white
- washable at 30° C



Application

For versatile use, especially in all areas requiring high hygiene standards. For single use in surgical areas, laboratories or clean rooms to protect staff and patients from infection and to prevent the spread of bacteria and diseases.

REF	Size	Box contents	Shipping unit	PZN
9837	39 x 8 cm	40 pcs	20 boxes	12671805
9836	44 x 8 cm	40 pcs	20 boxes	12671797

PE shoe covers

Shoe covers made of CPE with rubber band

Characteristics

Blue, with rubber band. For single use.

Dimensions (approx.): Length: 41 cm
 Height: 15 cm
 Thickness: 0.03 mm



REF	Box contents	Shipping unit	PZN
70-151	1000 pcs	1 box	11363007

Cellulose swabs

Characteristics

- bleached according to DAB
- absorbent and soft

Application

Cleaning and hygiene cloth for medical areas.



REF	Size	Box contents	Shipping unit	PZN
12431	20 x 20 cm	300 g	26 boxes	13655576

Layered cellulose

Crepe cuts in loose layers for versatile use

Characteristics

Layered cellulose, unbleached from 100% recycled cellulose as an environmentally friendly versatile consumable material.

Application

Applicable as underlay for non-sterile activities, as absorbent material or for cleaning and wiping.



REF	Size	Box contents	Shipping unit	PZN
102200D	20 x 20 cm	5 kg	3 bags	06803349
102300D	20 x 30 cm	5 kg	3 bags	06803355
104600D	40 x 60 cm	5 kg	3 bags	06803361

Rollicel® / Medical crepe

Bed covers on a roll



Characteristics

Bed cover made of cellulose, white, very absorbent and micro-perforated. In different sizes and designs as well as with and without coating.

Application

For all types of treatment beds and partially undressed examinations and treatments.

REF	Size	Layers	Coating	Shipping unit	PZN
3739-50	39 cm x 50 m	2 ply	No	9 rolls	11350476
3739-50B	39 cm x 50 m	2 ply	Yes	6 rolls	11350571
3739-100	39 cm x 100 m	2 ply	No	6 rolls	11350536
3745-100B	45 cm x 100 m	2 ply	Yes	6 rolls	11350884
3750-50	50 cm x 50 m	2 ply	No	9 rolls	11350482
3750-50B	50 cm x 50 m	2 ply	Yes	6 rolls	11350855
3750-100	50 cm x 100 m	2 ply	No	6 rolls	11350542
3755-50	55 cm x 50 m	2 ply	No	9 rolls	11350499
3755-50B	55 cm x 50 m	2 ply	Yes	6 rolls	11350861
3755-100	55 cm x 100 m	2 ply	No	6 rolls	11350559
3759-50	59 cm x 50 m	2 ply	No	9 rolls	11350513
3759-50B	59 cm x 50 m	2 ply	Yes	6 rolls	11350878
3759-100	59 cm x 100 m	2 ply	No	6 rolls	11350565
3759-100B	59 cm x 100 m	2 ply	Yes	6 rolls	11350890

Protective sheets

Disposable product as sitting and lying surface

Characteristics

Disposable product for use on intact skin. 2 ply cellulose with transparent PE film, which prevents the leakage of liquids. Free of latex, DEHP and PVC. Longitudinal reinforcement by resistant polyester seam as well as lateral blocking edges. 45 to 48 ply.



Application

Disposable protective sheet for covering of carrying, examination and transport beds, work tables and mattresses.

REF	Size (cm)	Shipping unit	PZN
1988	80 x 175	100 pcs	04353787
1991	80 x 100	200 pcs	04353741
1992	80 x 200	100 pcs	04353793
1995	80 x 210	100 pcs	04353801

Protective covering

Disposable product as sitting and lying surface, the ideal alternative to the protective sheet

Characteristics

Disposable product for use on intact skin. 1 ply cellulose with transparent PE film, which prevents the leakage of liquids. Free of latex, DEHP and PVC.



Application

Disposable protective sheet for covering of carrying, examination and transport beds, work tables and mattresses.

REF	Size (cm)	Lateral blocking edges	Shipping unit	PZN
1993	38 x 40	No	1000 pcs	04353681
1990	80 x 90	Yes	200 pcs	04353735
2982	80 x 140	Yes	100 pcs	01082431
2984	80 x 210	Yes	100 pcs	01082454

Protective covering "Eco-Line"

2983	80 x 180	Yes	100 pcs	01082448
------	----------	-----	---------	----------

Multi-purpose tissues

Characteristics

Multi-purpose cloth made of densely needled non-woven material. As a result, it is tear resistant, absorbent and does not fuzz.

Application

For all cleaning work in the kitchen, laboratory, hospital and hotel rooms.

REF	Cloth format	Box contents	Shipping unit	PZN
3539	32 x 38 cm	50 pcs	20 boxes	-



Bed covers

HDPE transparent

Characteristics

Protective film made of HDPE (High density polyethylene)

Application

For hospital and nursing beds, to protect the freshly made bed.



REF	Cloth format	Box contents	Shipping unit	PZN
1985	320 x 95 cm	200 pcs	1 roll	-

Universal tissues

Skin friendly universal tissues

Characteristics

Disposable product for cleaning and hygiene on intact skin. With optional version with higher viscose content for a more skin friendly sensation.

Application

For single use in all care areas.



Universal tissues

REF	Size	Box contents	Shipping unit	PZN
6051	approx. 30 x 30 cm (unfolded) approx. 15 x 15 cm (folded)	50 pcs	10 boxes	04355415

Universal soft tissues (thicker quality)

REF	Size	Box contents	Shipping unit	PZN
6052	approx. 30 x 30 cm (unfolded) approx. 15 x 15 cm (folded)	50 pcs	10 boxes	11722759

Micro-fibre tissues

Micro-fibre tissues in different colours, absolutely lint-free, with high absorption capacity and washable

Characteristics

Micro-fibre tissues for the highest demands, in different colours for better differentiation in the application. Fine microfibres remove even the finest dirt particles on virtually any surface, even without the addition of detergent additives.

Application

For lint-free applications with maximum absorption capacity.



REF	Size	Colour	Box contents	Shipping unit	PZN
6942	40 x 40 cm	blue	50 pcs	4 boxes	-
6943	40 x 40 cm	yellow	50 pcs	4 boxes	-
6945	40 x 40 cm	red	50 pcs	4 boxes	-
6944	40 x 40 cm	green	50 pcs	4 boxes	-

Wipes

Recyclable cellulose wipes

Characteristics

Multi-ply and semi-bleached cellulose wipes, free of latex, PVC and DEHP.

Application

Disposable product for use as a cleaning material.



REF	Size	Number of layers	Box contents	Shipping unit	PZN
6020	39 x 30 cm	3 ply	50 pcs	20 boxes	11857270
6021	39 x 30 cm	4 ply	50 pcs	20 boxes	11857264

Cosmetic tissues

Characteristics

Made of 2 ply cellulose, bright white.

REF	Size	Box contents	Shipping unit	PZN
6701	20.5 x 21 cm	100 pcs	35 boxes	-



V-fold paper hand towels

Characteristics

- 25 x 23 cm

1 ply, natural

REF	Box contents	Shipping unit	PZN
5722	20 x 250 pcs	1 pack	-

2 ply, white

REF	Box contents	Shipping unit	PZN
5721	20 x 160 pcs	1 pack	-



BaSick Bag

Sick bag

Characteristics

The safe and hygienic solution for nausea.
Semi-transparent, tear resistant polyethylene bag.
With scaling and practical locking ring.

Application

For single use in case of nausea.



REF	Size	Box contents	Shipping unit	PZN
9401	1.500 ml	25 pcs	20 bags	13818825

Kidney dish

non-sterile, cardboard

REF	Size	Colour	Shipping unit	PZN
9400	255 x 150 x 48 mm	grey	300 pcs	01082164



double sterile, plastic, pack of 1

REF	Size	Box contents	Shipping unit	PZN
4249	700 ml	12 x 1 pc	2 boxes	-



Disposable comb

Hair comb made of plastic, individually cellophane wrapped

REF	Colour	Box contents	Shipping unit	PZN
9679	brown	100 x 1 pc	20 boxes	11857672



Disposable razor

Characteristics

Stainless and non-sterile disposable razor for a gentle, safe and skin friendly shave for any hair length.



REF	Blade	Box contents	Shipping unit	PZN
9610	single edged	100 pcs	10 boxes	02501189

Urine bottle

Urine bottle for men 1000 ml

Characteristics

Opaque surface with scaling and clip closure. Therefore spill proof and odourless.

REF	Shipping unit	PZN
3006	25 pcs	-



Urine bottle holder

Urine bottle holder with chain and lid

Characteristics

Urine bottle holder for secure holding of the urine bottle.

Application

The urine bottle holder should be mounted (hung) on the bed frame of the bed within reach of the patient. Cleaning as required. Can be disinfected with a commercial disinfectant suitable for plastics. **WARNING:** When using disinfectants the instructions for use of the manufacturer must be observed.

REF	Unit	Shipping unit	PZN
3001	1 pc	20 pcs	-



Cotton buds

Characteristics

In different versions and variants for different applications. In sterile as well as non-sterile presentation and different head sizes as well as wood or flexible plastic version.



Application

For wide ranging use in all medical applications, such as for cleaning and application of liquid substances or ointments and for work in the laboratory.

non-sterile, in ZIP bag

REF	Length	Head	Handle	Box contents	Shipping unit	PZN
4046	7.5 cm	small head	plastic	10 x 200 pcs	10 boxes	11857620
1260	15 cm	small head	wood	10 x 100 pcs	20 boxes	06177225
1261	15 cm	medium head	wood	10 x 100 pcs	10 boxes	06177231
1263	15 cm	large head	wood	10 x 50 pcs	10 boxes	01363762
1268	15 cm	large head	plastic	10 x 50 pcs	10 boxes	11857494

sterile, pack of 2

1262	15 cm	small head	wood	100 x 2 pcs	6 boxes	06177254
1264	15 cm	large head	wood	100 x 2 pcs	3 boxes	01363779
1270	15 cm	large head	plastic	100 x 2 pcs	3 boxes	13711168

Lemon glycerine stick

Characteristics

- sugar free
- soaked with a solution of glycerine USP, citric acid USP, lemon aroma, sodium benzoate

Application

For professional use in everyday clinical practice.

non-sterile

REF	Box contents	Shipping unit	PZN
6804	25 x 3 pcs	20 boxes	03916886



Wooden tongue depressor

Our splinter-free spatula for examinations and the application of ointments



non-sterile

REF	Size	Box contents	Shipping unit	PZN
9802	150 x 18 mm	100 pcs	20 boxes	00822334

sterile, pack of 1

9665	150 x 18 mm	100 x 1 pc	5 boxes	01333821
------	-------------	------------	---------	----------

Bandage clips

For practical fixation of bandage ends

Characteristics

- elastic rubber band
- firmly grooved with sturdy, bend-proof metal hooks
- white



Application

For fast and secure fixation of non-cohesive bandages.

REF	Box contents	Shipping unit	PZN
4055	100 pcs	10 boxes	-

Nurse's stethoscope

Nurse's stethoscope for auscultation in adults and children

Characteristics

Acoustic flat head stethoscope, mechanical, consisting of ear hooks with ear tips, tube and chest piece.

Tube length: approx. 60 cm, latex-free.

Application

The use of stethoscopes may only be carried out by authorised specialist personnel. Depending on the examination, the stethoscope should be placed on the desired body region and used under low pressure. Please observe the enclosed operating instructions!



REF	Box contents	Shipping unit	PZN
9805	1 pc	100 pcs	11857459

Tourniquet

High quality tourniquets with innovative operation

Characteristics

Single hand operation without grasping. No pinching, clamping or squeezing of the skin. With safety lock. Latex and metal free.



REF	Item	Colour	PZN
5020	Meditrade tourniquets	Designer strap colourful	06177892
N20200	Jetpull	Designer strap blue	-

Dosicair®

Our infusion regulator with injection adapter

Characteristics

- Luer-lock connection system
- connectable to all infusion systems
- increased accuracy and uniformity

Application

For infusion therapy.



REF	Box contents	Shipping unit	PZN
DF100	50 pcs	8 boxes	-

Urine beaker

REF	Size	Box contents	Shipping unit	PZN
3009	125 ml	100 pcs	5 boxes	04354427
3010	100 ml	500 pcs	1 pack	-



Ultrasound scan gel

Characteristics

Ultrasound scan gel for sonography with ultrasound probes. Outstanding contact creation for an optimum representation of even the smallest parts and structures. Highest penetration rate, sonic head protection, micro-bubble-free and skin-friendly.

Application

For all forms of sonography.

REF	Box contents	Shipping unit	PZN
6831	250 g	30 bottles	01904778




Repose®

An advanced concept in the fight against decubitus and pain relief



Repose® : Prevent and treat pressure ulcers

Repose® consists of air-filled, interconnected polyurethane membranes. The polyurethane film with a thickness of 50 µ is extremely elastic: It can be stretched up to 6 times without tearing or losing its elastic properties. The combination of air with polyurethane guarantees the best possible pressure distribution by means of sinking and covering. The material is breathable, waterproof and temperature neutral. Repose® is very easy to use and can be implemented very quickly both preventively and therapeutically in all situations. With Repose®, pressure, shear forces, frictional forces and the micro-climate can be effectively controlled.



- R** evolutionary
- E** fficient
- P** roprious
- O** ptimal pressure distribution
- S** wiftly and easy handling
- E** xtraordinarily

Advantages:

- Declining prevalence and incidence of pressure ulcers
- Optimum protection in bed and on a chair
- Excellent care of all patients
- Cost saving
- Time saving (care and logistics)
- Space saving
- Easy to use

REF	Item	Dimension	PZN
7951	Mattress	190 x 77 x 5 cm	4353646
7953	Seat cushions	46 x 43 x 7 cm	4353652
7955	Heel wedge	70 x 40 cm	4353617
7954	Heel protector (1 pair)	38 x 20 cm	4353623
7964	Heel protector (1 piece)	38 x 20 cm	6640551

Cold-warm compress

Compress for cold and heat application

Characteristics

Permanent compress with gel filling, easy to model. For freezer compartment (flexible up to -30° C), water bath and microwave suitable. Follow instructions for use or information on the label.

Application

For sprains, strains, swelling and burns. Also for relief of fever, headaches and bruises. Follow instructions for use or information on the label. There is the right size and shape for every application.



	REF	Size	Box contents	Shipping unit	PZN
Gel colour: transparent	9668	13 x 14 cm	50 pcs	1 pack	04352500
	9870	12 x 29 cm	25 pcs	1 pack	04352463
	9666	16 x 26 cm	25 pcs	1 pack	04352517
	9661	30 x 40 cm	6 pcs	1 pack	04352546
	9663	7.5 x 35 cm	40 pcs	1 pack	04352552
Gel colour: blue	9465	16 x 25 cm	20 pcs	1 pack	-

Non-woven cover for compresses

Non-woven cover for cold-warm compresses

Characteristics

Our non-woven covers for cold-warm compresses protect the skin, since direct contact between the compress and the skin is avoided. For reasons of hygiene only intended for single use.

REF	Size	Box contents	Shipping unit	PZN
9682	13 x 14 cm	500 pcs	1 pack	-
9683	12 x 29 cm	500 pcs	1 pack	-
9684	16 x 26 cm	500 pcs	1 pack	-
9687	21 x 38 cm	500 pcs	1 pack	-
9693	30 x 40 cm	500 pcs	1 pack	-

Rescue blanket

Characteristics

- colour: gold / silver
- reflects body heat or sun rays
- danger of hypothermia is reduced or heat kept away from the body
- packed individually with instructions for use

Application

Cold and heat protection for injured persons. Suitable for accidents, wind, rain and heat.



REF	Size	Box contents	Shipping unit	PZN
4070	160 x 210 cm	50 pcs	5 boxes	13510895

Medicine beaker

Medicine beaker, 30 ml, graduated



Characteristics

Made of highly transparent polypropylene, content optimally visible, unbreakable with lateral graduation.



REF	Colour	Box contents	Shipping unit	PZN
Medicine beaker				
9729	transparent	75 pcs	50 boxes	-
9727	blue	75 pcs	50 boxes	-
9725	red	75 pcs	50 boxes	-
9726	yellow	75 pcs	50 boxes	-
9728	green	75 pcs	50 boxes	-

REF	Colour	Box contents	Shipping unit	PZN
Medicine beaker lid				
9700	white	750 pcs	5 boxes	-

Medication dispenser

Characteristics

Tablet dispenser made of PP / PVC with transparent sliding lid. Supports regular and timely medication in the morning, at noon, in the evening and at night.



REF	Colour	Box contents	PZN
9707	white	200 pcs	-
9708	blue	200 pcs	-

Medi-Box

Container for safe cannula and syringe disposal

Characteristics

Round disposal container, with lid and round opening, flat closure, pull-off and calibration aid.

Application

For safe disposal of cannulas and syringes.

REF	Size	Box contents	PZN
9694	0.5 l	1 pc	-
9702	2.5 l	1 pc	-



Wall holder

Wall holder for glove boxes of all sizes, made of acrylic or polystyrene

Characteristics

- environmentally friendly materials
- very easy assembly handling
- glued stable



REF	Item	Description	Size (in mm)	Shipping unit	PZN
4040	Dracula	PS wall holder without imprint (white)	universal	1 pc	-
4050	Ghost	Acrylic wall holder with imprint	250 x 145 x 100	1 pc	-
4059	Ghost	Acrylic wall holder without imprint	250 x 140 x 70	1 pc	-
4060	Ghost	Acrylic wall holder without imprint	250 x 145 x 100	1 pc	-

Refuse bags HD

PE refuse bags

Characteristics

HDPE (High density polyethylene) refuse bags for collection of waste products. For single use.



REF	Size	Colour	Box contents	Shipping unit	PZN
9747	49 x 60 cm (approx. 24 Litres)	transparent	50 pcs	40 rolls	-
9748	63 x 74 cm (approx. 60 Litres)	transparent	50 pcs	40 rolls	-

Refuse bags LD

PE refuse bags

Characteristics

Refuse bags made of LDPE (Low density polyethylene) for the collection of waste products. For single use.



REF	Size	Colour	Box contents	Shipping unit	PZN
9646	70 x 110 cm (approx. 120 Litres)	black	15 pcs	10 rolls	-
9651	70 x 110 cm (approx. 120 Litres)	blue	25 pcs	10 rolls	-
9746	50 x 60 cm (approx. 30 Litres)	Grey	50 pcs	20 rolls	-
9749	63 x 72 cm (approx. 60 Litres)	transparent	50 pcs	20 rolls	-
9750	63 x 85 cm (approx. 65 Litres)	white	40 pcs	9 rolls	-
9760	57.5 x 100 cm (approx. 70 Litres)	blue	25 pcs	10 rolls	-
9761	70 x 110 cm (approx. 120 Litres)	blue	25 pcs	10 rolls	-
9762	70 x 110 cm (approx. 120 Litres)	transparent	25 pcs	10 rolls	-
9764	60 x 80 cm (approx. 60 Litres)	Grey	25 pcs	10 rolls	-
9765	70 x 110 cm (approx. 120 Litres)	red	25 pcs	10 rolls	-
9766	57.5 x 100 cm (approx. 70 Litres)	yellow	25 pcs	10 rolls	-
9768	70 x 110 cm (approx. 120 Litres)	blue	15 pcs	10 rolls	-

Kitchen towel

Characteristics

- bright white
- made of cellulose
- 50 sheets per roll
- roll length: 11.2 m
- sheet width: 25.7 cm
- sheet length: 22.4 cm



3 ply

REF	Box contents	Shipping unit	PZN
57301	4 rolls	8 boxes	-

Toilet paper

Characteristics

- bright white
- made of cellulose
- 250 sheets per roll
- roll length: 28 m
- sheet width: 9.5 cm
- sheet length: 11.2 cm



2 ply

REF	Box contents	Shipping unit	PZN
9200	8 rolls	8 boxes	-

3 ply

REF	Box contents	Shipping unit	PZN
9201	8 rolls	9 boxes	-

Toilet paper

Eco-friendly premium toilet paper


Characteristics

- white
- made from recycled cellulose
- Eco-label
- 200 sheets per roll
- roll length: 27.6 m
- sheet width: 9.6 cm
- sheet length: 13.8 cm




3 ply

REF	Box contents	Shipping unit	PZN
1210	8 rolls	9 boxes	-



Notes:



Notes:

Index

A

Abdominal cloths 64 / 65
 ABE® Absorbent adhesive dressing 87
 ABE® Absorbent dressing pad 84
 ABE® Body bandage 110
 ABE® Cannula retention dressing I.V.
 89
 ABE® compression bandage roll 107
 ABE® derm (plus) 88
 ABE® drainal 83
 ABE® extra / ABE® extra adhesion
 100 / 101
 ABE® fix adhesion c 99
 ABE® fix color 96
 ABE® Fixation non-woven fabric 89
 ABE® Gauze bandage 93
 ABE® Gilchrist bandage 109
 ABE® Head bandage 91
 ABE® Heart catheter pressure bandage
 111
 ABE® Ideal 100 102
 ABE® Injection plaster 86
 ABE® Last 94
 ABE® Last adhesive 97
 ABE® Last crepe 95
 ABE® lastic 103
 ABE® light 86
 ABE® Limb holder 111
 ABE® Nose bandage 92
 ABE® Padding bandage 106
 ABE® Ribbon gauze 73
 ABE® silk 90
 ABE® tracheal 83
 ABE® Tubular bandage 108
 ABE® Tubular net bandage 91
 ABE® universal 104
 ABE® Ven 105
 Absorbent dressing pads 84
 Absorbent pads 134
 Active gel 157
 Alcohol swabs 72
 Alcoman® + 138
 Alcoman® Gel 139

All purpose bandages 96 - 100
 All purpose cleaner 164 / 165
 Anti-decubitus heel protector 192
 Anti-decubitus heel wedge 192
 Anti-decubitus mattresses 192
 Anti-decubitus seat cushion 192
 Anti-slip sock 180
 Application aids 159
 Arm bandages 109
 Arm cuffs 111
 Astronaut hoods 55

B

Baby nappies 135
 BagBath® (Hair care) 176
 Bandage change kit 120
 Bandage clips 189
 Bandages 91 - 111
 BaSick bag 186
 Bed covers 183
 BeeSana® Absorbent dressing pad 84
 BeeSana® Absorbent pad 127 / 128
 BeeSana® Bandage change kit 120
 BeeSana® Cellulose swabs 72
 BeeSana® Comflock Super 133
 BeeSana® Eye pad 85
 BeeSana® Gauze ball 66 / 67
 BeeSana® Gauze cutting 74
 BeeSana® Gauze (slit) swab 75 - 80
 BeeSana® Lap sponge 64 / 65
 BeeSana® Male incontinence pad 129
 BeeSana® Medical gauze 74
 BeeSana® Non-woven (slit) swabs
 80 / 82
 BeeSana® Non-woven
 disposable panty 179
 BeeSana® PEG bandage kit 119
 BeeSana® Pointed swab 70 / 71
 BeeSana® PP gown 57
 BeeSana® PP / PE gown 58

BeeSana® PP protective gown 57
 BeeSana® Preparation balls 69
 BeeSana® Rectangular pad 128
 BeeSana® SafeCut 98
 BeeSana® SMS gown 58 / 59
 BeeSana® SSS gown 59
 BeeSana® Suture removal kit 120
 BeeSana® Thread swabs 68
 BeeSana® Ward catheter set 118
 BeeSana® Ward dialysis kit 121
 BeeSana® Ward protection set 121
 BeeSana® Wash glove 178
 BeeSana® Women's incontinence pad
 129
 beopur WC fresh 164
 BioSan® Fresh scent 165
 Body bandages 110
 Body washing systems 176
 Box holder 196

C

Canister dosing pump 159
 Canister / drain taps 159
 Canister key 159
 Cannula collecting container 195
 Cannula retention dressing / plasters
 89
 Cardiac catheter pressure bandages
 110
 Care foam 158
 Care wipes 158
 Cellulose swabs 181
 Cleaning agents 166 - 173
 Cold-warm compresses 193
 Combs 186
 Comflock Super 133
 Compress cover 193
 Compresses 75 - 82, 193
 Compression pads 100 - 107
 Copolymed™ 30

Copolymer gloves 30
 Cosmetic tissues 185
 Cotton buds 188
 Cotton gloves 26
 Cream lotion 156

D

Dental masks 50 - 54
 Disinfectant dispenser 160
 Disinfection cleaner 166 - 173
 Disinfection (products) 136 - 154
 Disinfection wipes
 144 / 145, 147 / 148
 Disposable apron 60
 Disposable body wash system 176
 Disposable combs 186
 Disposable gloves 9 - 47
 Disposable razor 187
 Disposal container 195
 Dolicair® 191
 Dosing aids 159 / 160
 Dosing bottle 159
 Dosing pumps 159
 Dressing pads 75 - 85, 128, 193
 Dressing plasters 86 - 88
 Dressings 112 - 115 / 86 - 88

E

Élastify® 27
 Ethasept® 140
 Examination gloves (non-)sterile
 9 - 25, 31 - 33
 Eye pads 85

F

Face mask 48 - 54
 Film wound dressings 88
 Finessis® Aegis 44 / 45
 Finessis® Corium 46
 Finessis® Zero 43
 Finger bandages 98
 Finger cots 28
 Fixation bandages 91 - 99
 Fixation non-woven fabric 89
 Fixation pants 130 / 131
 Fixation plasters 89 / 90
 Flexylon gloves 42 - 46
 Flock (medical) undersheets 132 / 133
 Foam spray heads 159
 Fresh scent 165
 Froggy™ 27

G

Gauze ball 66 / 67
 Gauze bandages 93
 Gauze (slit) 74
 Gauze swabs 73 - 78
 Gentle Med® Active gel 157
 Gentle Med® Care foam 158
 Gentle Med® Cream lotion (W/O) 156
 Gentle Med® Moist care wipes 158
 Gentle Med® MRSA wash lotion 155
 Gentle Med® Oil bath 157
 Gentle Med® Skin care lotion (O/W)
 156
 Gentle Med® Skin cleansing and care
 158
 Gentle Med® wash lotion 155
 Gentle Skin® Aloecare™ 14
 Gentle Skin® black 15
 Gentle Skin® Classic (x-long) 10
 Gentle Skin® Compact + 11
 Gentle Skin® Grip 11
 Gentle Skin® HiRisk 13
 Gentle Skin® Isopretex® OP 40

Gentle Skin® Micro OP® 38
 Gentle Skin® Premium OPT™ 36
 Gentle Skin® Securitex OP® 39
 Gentle Skin® Sensitive 9
 Gentle Skin® Smooth OP 37
 Gentle Skin® steril 31
 Gentle Tec® 12
 Gilchrist bandages 109
 Glove box holder 196
 Gloves 8 - 47
 Gowns 57 - 60

H

Hair combs 186
 Hair washing system 176
 Hand disinfection 136 - 140
 Hand spray head 159
 Head bandages 91
 Heel wedge 192
 Heel protector 192
 Hoods 55
 Household gloves 29
 Hygiene gloves 26 / 27

I

Incontinence pads 127 - 129
 Incontinence products 122 - 133
 Incontinence slips 124 / 125
 Incontinence trousers 130 / 131
 Incontinence undersheets 132 / 133
 Ideal bandages 102
 Infusion regulator 191
 Injection plasters 86
 Instrument disinfection 149 - 154
 Instruplus® Forte 153 / 154
 Instrusol® AF + 149 / 150
 Instrusol® Enzym 151 / 152

Index

K

Kidney dishes 186
 Kitchen towel 198
 Kleen Purgatis cleaning and
 disinfectant agents 166 - 173
 Knitted tube bandages 108

L

Latex finger cots 28
 Latex gloves
 9 - 15, 28, 29, 31, 35 - 39, 47
 Layered underpad 132
 Layered cellulose 181
 Leg holder 111
 Lemon glycerine sticks 188
 Liberty Light 126
 Light fixation pants 130
 Liner gloves 26
 Long stretch bandages 105 / 106

M

Male incontinence pad 129
 Masks 48 - 54
 Maternity inner pads 128 / 129
 Measuring beaker 159
 Medi-Box 195
 Medical crepe 182
 Medical gauze 74
 Medication dispenser 195
 Medicine beaker 194
 Medi-Pants 125
 Medizid® AF 141
 Medizid® Alpha 146
 Medizid® Alpha (disinfection wipes)
 147
 Medizid® Alpha Flowpack 148
 Medizid® Rapid 143
 Medizid® Rapid (disinfection wipes)
 144

Medizid® Rapid Optima 142
 Medizid® Rapid QF Flowpack 145
 Membrane bandages 112 - 115
 Mesh tube bandages 91 / 92
 Micro-fibre tissues 184
 Middle stretch bandages 104
 Moist (care) wipes 158
 Mouth spatula 189
 MRSA-effective wash lotion 155
 Multi-purpose tissues 183

N

Nappies / nappy pants 124 - 126, 135
 Neoprene gloves 41
 Neopretex® OP 41
 New Quik 135
 Nitril® 3000 16
 Nitril® 3000 X-Long Blue 19
 Nitril® BestGen® 18
 Nitril® black 20
 Nitril® LIOX 22, 23
 Nitril® magenta 21
 Nitril® NextGen® 16
 Nitril® Sensory® 17
 Nitril® steril 32
 Nitrile gloves 16 - 23, 32
 Nose bandages 92
 Nose sling 92
 Non-woven cover for compresses 193
 Non-woven disposable panty 179
 Non-woven swabs 81 / 82
 Nurse caps 55
 Nurse's stethoscope 190
 Nursing pads 134

O

Oil bath 157

P

Padding bandages 106
 Panty fixation pants 131
 Paper towels 185
 PE aprons 58
 PE disposable aprons 60
 PEG-Bandage kit 119
 PE-Gloves 27, 30
 PE shoe covers 181
 Plasters 86 - 88
 Pointed swab 70 / 71
 Polychloroprene gloves 41
 Polyisoprene gloves 40
 PolyMem 112 - 115
 PP / (PE) Gowns / Jackets / Aprons
 57 - 60
 Preparation balls 69
 Pressure bandages 107
 Protective bibs 179
 Protective covering 183
 Protective Jackets / Gowns / Aprons
 57 - 60
 Protective sheets 182
 Protective towel 178

R

Radiation protection gloves 47
 Razor 187
 Rectangle pads 128
 Recycling toilet paper 199
 Reference™ 9
 Reference™ OP 35
 Refuse bags 197
 Repose® 192
 Rescue blankets 194
 Ribbon gauze 73
 Rollicel® / Medical crepe 182
 Roll plaster 90

S

SafeCut 98
 Sets 118 - 121
 Shoe covers 181
 Shoulder bandage 109
 Skin care lotion 156
 Skin cleansing foam 158
 Slit (hole) swabs 81 - 83
 SMS gown 58 / 59
 Soap dispenser 160
 Socks 180
 Special bandages 108 - 111
 Spit bag (BaSick bag) 186
 SSS gown 59
 Stethoscopes 190
 Stretch bandages 100 - 106
 Suavel® Antifluid 53
 Suavel® Antifog 53
 Suavel® Apollo 55
 Suavel® Astrid 55
 Suavel® Comfort 52
 Suavel® Comfort Plus 52
 Suavel® Protec 51
 Suavel® Protec Plus 51
 Suavel® Sensima 54
 Support bandages 107 - 111
 Superform® 126
 Superform® Air 127
 Superslip® 124
 Superslip® Air Comfort 125
 Surface disinfection 141 - 148
 Surgical - Dental masks 48 - 54
 Surgical - Face mask 48 - 54
 Surgical - Gloves 34 - 47
 Surgical - Jackets / Gowns / Aprons
 56 - 60
 Surgical - Hoods 55
 Surgical - Masks 48 - 54
 Surgical - Sets 118 - 121
 Surgical - Socks 180
 Suspensory 92
 Suture removal kits 120
 Swabs 66 - 73

T

Tablet dispenser 195
 Tamponade 73
 Testicular roll 107
 Thread swab 68
 Toilet paper 198 / 199
 Tourniquet 190
 Tubular bandages 91, 108

U

Ultrasound scan gel 191
 Underpads 132 / 133
 Universal tissues 184
 Urine beaker 191
 Urine bottle (holder) 187

V

V-fold paper hand towels 185
 Vinyl 2000 24
 Vinyl 2000 Stretch PF 25
 Vinyl gloves 24 / 25
 Visitor gown 57 - 60

W

Wall holder 196
 Ward dialysis kits 121
 Ward catheter sets 118
 Ward gowns / aprons 56 - 60
 Ward protection set 121
 Wash glove 178
 Washing lotions 155
 Wipes 144 / 145, 147 / 148, 185
 Women's incontinence pad 129
 Women's sanitary towels 129
 Wooden tongue depressor 189
 WC cleaner 164

X

X-Guard 47
 X-ray protective gloves 47



Senator Hans Hermann Rösner Mautby

(* 27 June 1917 in Beuthen, † 24 June 1999 in Oberaudorf)

For his services to the development of the clinic in Oberaudorf he was appointed as an honorary senator of the Ludwig Maximilian University of Munich. In addition, he was awarded the Federal Cross of Merit and the Bavarian Order of Merit because of his entrepreneurial spirit in health care.



Touch. Help. Create opportunities.

With a comprehensive framework of action and targeted public relations work, the Senator Rösner Foundation collects funds for people who have been in financial difficulties due to cancer. In addition, the founders Peter Rösner and his wife Eveline Rösner use their foundation to support cancer research and the further education of medical professionals.

As a successor to the former "Relief organization of the Bad Trissl clinic", the Senator Rösner Foundation can now draw on more than 25 years of experience in the field of foundations.

These projects are currently being sustainably supported.

HEAT study

The Senator Rösner Foundation supports the clinical research study of Prof. Dr. med. Dipl. Biochem. Rolf D. Issels at the University of Munich Hospital with innovative approaches in the field of "pancreatic carcinoma" with several hundred thousand Euro.

The power of fast help

Cancer leaves some people in unexpected financial and social emergencies. The Senator Rösner Foundation "Living with Cancer" helps them quickly and without the red tape.

Competence promotion and further education

The Senator Rösner Foundation uses its proximity to all major oncology centres in Munich and Bavaria. It maintains permanent contacts with scientists and researchers and provides funding for research and development.

We run - we help

With this slogan and under the auspices of FIS World Cup winner Felix Neureuther, the fund raising run of the Senator Rösner Foundation is held, which has now become established nationwide. Alone in 2013, the magnificent distance of 6756 km was achieved by almost 800 runners and rewarded with a record breaking donation amount of € 105,000.

WSG Swarovski Wattens

SOS Children's villages

A trend setting aid project beyond the state borders. Thanks to an amicable connection between the Rösner family and the club president Diana Langes-Swarovski, children and adolescents are also supported by the 'SOS Kinderdörfer' beyond the development of sports to be able to cope with difficult life situations.

www.wsg-wattens-fussball.at

Contact: senator-roesner-stiftung.de

Corporate development

1982	Founding of Hermes Trans Atlantic Trading as an offshoot of our Canadian companies	1995	Asset deal of the company MCM Landhut
1988	Founding of Meditrade GmbH in its current orientation with the company motto: "Protection across the board" - hand protection, face protection, caps, surgical gowns, disinfection, incontinence products, dressings, compresses, swabs, abdominal cloths, etc.	1994-1995	Reaching the break even point
1989	Participation with 50% in the company Meditrade Austria (separation 1992) participation in the company Meditrade Switzerland (separation 1992)	1996	Participation in Eurocare GmbH and Mediclean Bayern GmbH
1990	Sales acquisition of the companies Cosmotec and Sauter Trading	1996	Foundation of Meditrade International Thailand Ltd. with subsequent share purchase in the joint venture Siamount Plastic Ltd. Bangkok with Mitsubishi Group (vinyl factory)
1990-1992	Purchase of the first glove factory in Thailand, Hycare Ltd. Hatyai	until 2002	Continuous range expansion with exclusive contracts and minority interests in Italy
1992	Founding of Meditrade Dresden (today Praximed Meditrade DD)	2003	Construction of another factory in Thailand Hatyai, Rösner-Mautby Thailand Ltd (DVD)
1994	In addition, acquisition of AiA Leinberger GmbH Grossräschen	until 2006	Sale of all Thai factories
1992	Relocation to the new company premises in Thierseestr. 196-198, Kiefersfelden	2008	Equity investment in Adventa Health Group AG, Malaysia (world market leader in surgical glove production)
1993	Acquisition of the company Hese GmbH in Nuremberg-Feucht, cellulose pulp and underpads, production VP-Schickedanz / Temca	2009-2011	Acquisition of Beese GmbH & Co. KG in Hamburg with production - Subsidiary in China
1992-1994	Re-establishment of the companies Hycare Ltd. UK / Barrier Healthcare Ltd. UK in Cambridge	2012	Moving into new company premises Medipark 1, Kiefersfelden
		2016	Majority share in the MaiMed Group
		2018	Expansion of the logistics centre in Neuenkirchen, 30 years of Meditrade GmbH





Meditrade GmbH

Medipark 1
83088 Kiefersfelden
Germany
Phone +49(0)8033/9760-0
Fax +49(0)8033/9760-70

**For questions about our product range
and further information please do not
hesitate to contact our customer service:**

Phone +49(0)8033/9760-0
Email: info@meditrade.de

**You can also find us on the internet
at www.meditrade.de**

On the website you can find
the current product range
and more up to date information.



Förderung nachhaltiger
Waldwirtschaft
www.pefc.de