

ABE[®] derm

Transparenter Folienwundverband zur Fixierung von Kathetern, steril / Transparent film wound bandage for fixation of catheter, sterile



Varianten / variants

REF	Größe / Sizes	Stück pro Packung/ Piece per pack	Sterilisation/ Sterilization	PZN
9571	6 x 7 cm	100	EO	06790204
9572	10 x 12 cm	50	EO	06790210
9573	15 x 20 cm	50	EO	06790227

Allgemeine Informationen / General information

Materialzusammensetzung / Material composition

Trägermaterial: Polyurethanfolie, Kleber: hypoallergenes Polyacrylat, beidseitiger Papierträger: weißes Pergamentpapier, einseitig silikonisiert. / holder material: polyurethane film, glue: hypoallergenic polyacrylate, paper backing on both sides: white parchment paper, siliconized on one side.

Code nach UMDNS / UMDNS code

Heftpflaster 10-026 / adhesive plaster 10-026

Qualitätszertifikate / Quality certificates

EN ISO 13485:2016; 93/42/EWG / EEC: Medizinprodukt der Klasse Is gem. Anhang IX (Regel 4) / Medical product of class I acc. to appendix IX (rule 4)

Anwendungshinweise / Application notes

Zur Fixierung von Intravenösen Kathetern. / For fixing of intravenous catheters in place.

Haltbarkeit / Shelf life

5 Jahre bei produktgerechter Lagerung: kühl und trocken lagern. / 5 years if stored properly: store in a cool and dry place.

Wichtiger Hinweis: Nichtbeachtung der Informationen, insbesondere zur Verwendung, Verträglichkeit und Beständigkeit kann zu Personen- und Sachschäden führen. Die MaiMed GmbH übernimmt keine Haftung für falschen Einsatz der Produkte. Es handelt sich um ein Produkt zum Einmalgebrauch – Materialveränderungen / Infektionsrisiken können bei Mehrfachgebrauch nicht ausgeschlossen werden. Die hier enthaltenen Informationen entsprechen dem Zeitpunkt des letzten Änderungsstands. Irrtümer, Druckfehler sowie Änderungen behalten wir uns vor. Holen Sie im Zweifelsfall fachkundigen Rat ein.

Important Note: Non-observance of the information, especially in view of use, compatibility and durability may lead to damage to persons and / or material. MaiMed does not take any liability in case of improper use. This is a single use product – changes in material / infection risk cannot be excluded in case of multiple use. All information included correspond with the last document amendment date. Errors, literal mistakes and changes expected. In case of doubt please seek expert advice.